**ENDOSCOPIC DRAINAGE STENTS
ENDOPROTHÈSE DE DRAINAGE
ENDOSCOPISCHE DRAINAGE-STENTS
STENT ENDOSCÓPICOS DE DRENAJE**

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EN - INSTRUCTION FOR USE

Intended Use

Endoscopic Drainage Stents are intended to drain biliary, cysto and/or pancreatic obstructions via endoscopes.

Product Description

The device is made of a radiopaque stents with or without flaps and may or not come pre-mounted with it's application device. If the application device is not present with your stent, please consider having one at hand. (Refer to the current catalogue for more information)

Advices before first application

Please read carefully and follow all safety operating instructions and warnings before first application of the device.

A previous knowledge regarding handling and operation is required and essential. Unpack the instrument carefully and examine the devices for any possibility of damage. In case of any damage or missing items contact your distributor immediately.

If the package present any damaged, the sterility of the device is not guarantee.

Handling and Operation

The device is placed in conjunction with an endoscope, a guide wire, a pusher and, when necessary, a guiding catheter

Procedure

To facilitate the application, an endoscopic sphincterotomy prior to guide wire placement is recommended.

For the application please proceed as follows:

- When implanting stents with a diameter of more than 7 Fr., it is strongly recommended to apply a guiding catheter. Introduce this catheter over the guide wire and position it on top of the stenosis also.
- Choose type and length of stent to be placed.
- Use the included positioning sleeve for easier introduction by sliding it onto the proximal end of the stent (the rear flap should fit back into it).
- Introduce the device (with positioning sleeve) and the pusher over the proximal end of the guide wire (guiding catheter) into the operating channel. The positioning sleeve does not fit into the channel and will stay outside.
- Place the device by using the pusher and carefully place the device in the desired position.

- Pull the guide wire (and guiding catheter) back into the operating channel while maintaining the position of the stent with the help of the pusher. After removing guide wire, guiding catheter and pusher the stent will return to its original shape.

Storage

This device is delivered sterilized and in order to maintain this state until its first use it must be kept in its original sterilized package.

Do not put any objects on the instrument or its package!
Do not store the instruments near aggressive chemical products!

Do not expose the instruments to direct or indirect sunlight or other ultra-violet rays!

Keep in dry area

Complaints will not be taken under consideration if the instruments have been stored improperly.

Cleaning and Sterilization

This device is designed and warranted for single use!

Risk in case of reuse

This device is for single use only. G-Flex denied all responsibility in case of reuse. The reuse of single use device might present risk to the safety of the user or patient, due to possible uncontrolled contamination and / or unreliability of the device performance.

User

The users of G-FLEX instruments must be specialists in their fields. An appropriate and specific training for preparation, care and maintenance of the flexible instruments is required.

Repairing and complains

In health protection of our employees only disinfected or sterilized instruments will be accepted for analysis or repairing. In any case disinfection or sterilization dates and validity must be labelled outside the package. If this requirement has not been fulfilled the instruments will be returned to the sender without analysis or repairing.

Legal foundation

The Law of the European Union applied.

Product support

In case of questions or difficulties concerning our instruments please contact your local distributor or G-FLEX directly during regular working hours.

Times available: Monday to Friday; 9am to 5pm (MET)

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G-Flex USA, Inc., 74 Southwood Pkwy Suite 600