**DISPOSABLE GUIDING CATHETER, PUSHER
AND STENT APPLICATION SYSTEM
SYSTÈME JETABLE DE CATHÉTER DE
GUIDAGE, POUSSOIR ET APPLICATION
D'ENDOPROTHÈSE
EENMALIG TE GEBRUIKEN BEGELEIDENDE
KATHETER, DRUKKER EN
STENTPLAATSINGSSYSTEEM
INTRODUCTORES**

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Current Name: IFU-
GuidingCatheterPusherStenApplicatio
nSystem-Disposable

Creation Date: 23/02/2017 Version: 6

This document is only controlled in electronic format by the QM.

EN - INSTRUCTION FOR USE

Intended Use

Disposable Guiding Catheter, Pusher and Stent Application System are intended to place transpapillary biliary stents in patients with biliary obstructions via endoscopes.

Product Description

The device is made of an inner guiding catheter with stainless steel marker at the distal end and a Y connector at the proximal end. They also have an outer pusher with stainless steel marker at the distal end and a Luer-lock adapter at the proximal end.

Advices before first application

Please read carefully and follow all safety operating instructions and warnings before first application of the device.

A previous knowledge regarding handling and operation is required and essential. Unpack the instrument carefully and examine the devices for any possibility of damage. In case of any damage or missing items contact your distributor immediately.

If the package present any damaged, the sterility of the device is not guarantee.

Handling and Operation

The device can only be used with a duodenoscope and a guide wire. The device is not suitable for placing Pigtail- or Double-Pigtail-Stents

Endoscopic sphincterotomy and selective duct cannulation prior to guide wire placement is recommended. (Please observe HF-generator and sphincterotome manufacturer's instructions).

Procedure

- Steer the tip of the duodenoscope close to the papilla.
- Insert the guide wire through the operating channel into the biliary duct and pass the stenosis.
- Use the included positioning sleeve for easier introducing by sliding it to the proximal end of the stent. So the rear flap will fit back into it.
- Ensure that the Luer-Lock connection between pusher and guiding catheter is firmly closed. Loosen the sealing cap at the Y-adapter.
- Place the biliary stent on to the device by sliding it over the distal end of guiding catheter up to the pusher.
- Now introduce the loaded the device over the proximal end of the guide wire into the operating channel. Positioning sleeve does not fit to the channel and will stay outside.

- Advance the device under constant fluoroscopic control carefully until the radiopaque metal tip of the guiding catheter passes the stenosis.
- Closing the sealing cap at the Y-adapter prevents movement of the guide wire inside of the guiding catheter. Contrast injection is possible through the Luer-Lock-Port at the Y-adapter.
- Open up the Luer-Lock connection between guiding catheter and pusher and advance the biliary stent to its final position by moving forward the pusher carefully. **ATTENTION:** Keep guiding catheters position while placing the stent !
- Pull back the guiding catheter and guide wire into the operating channel. Prevent movement of the stent by keeping the pusher in its last position. After that the stent will return to its original shape.
- Finally, remove pusher, guiding catheter and guide wire out of the duodenoscope.

Storage

This device is delivered sterilized and in order to maintain this state until its first use it must be kept in its original sterilized package.

Do not put any objects on the instrument or its package!
Do not store the instruments near aggressive chemical products!

Do not expose the instruments to direct or indirect sunlight or other ultra-violet rays!

Keep in dry area

Complaints will not be taken under consideration if the instruments have been stored improperly.

Cleaning and Sterilization

This device is designed and warranted for single use!

Risk in case of reuse

This device is for single use only. G-Flex denied all responsibility in case of reuse. The reuse of single use device might present risk to the safety of the user or patient, due to possible uncontrolled contamination and / or unreliability of the device performance.

User

The users of G-FLEX instruments must be specialists in their fields. An appropriate and specific training for preparation, care and maintenance of the flexible instruments is required.

Repairing and complains

In health protection of our employees only disinfected or sterilized instruments will be accepted for analysis or repairing. In any case disinfection or sterilization dates

and validity must be labelled outside the package. If this requirement has not been fulfilled the instruments will be returned to the sender without analysis or repairing.

Legal foundation

The Law of the European Union applied.

Product support

In case of questions or difficulties concerning our instruments please contact your local distributor or G-FLEX directly during regular working hours.

Times available: Monday to Friday; 9am to 5pm (MET)

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G-Flex USA, Inc., 74 Southwood Pkwy Suite 600