PUNCTION NEEDLE DISPOSABLE AGUJA DE PUNCION DESECHABLES

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EN - INSTRUCTION FOR USE

Intended Use

Punction Needle are intended to perforate intragastric balloon (either air or liquid) and if needed to aspirate the content inside the balloon.

Product Description

The device is made of a flexible catheter with an inner tube connected to a whole needle at the distal end and a connector at the proximal end

Advices before first application

Please read carefully and follow all safety operating instructions and warnings before first application of the device.

A previous knowledge regarding handling and operation is required and essential.

Unpack the instrument carefully and examine the devices for any possibility of damage. In case of any damage or missing items contact your distributor immediately.

If the package present any damaged, the sterility of the device is not guarantee.

Handling and Operation

Check to ensure that the needle is operating properly by pushing the needle piston in and out while holding the catheter.

Procedure

- Introduce the catheter inside the endoscope with the needle retracted
- Push the needle piston forward and screw it in order to lock the needle
- Perforate the balloon with the needle
- If a drainage of the balloon is necessary, place the whole needle and part of its outer catheter inside the balloon, unscrew the needle and remove it
- Apply a suction at the proximal end, when finish remove the catheter from the endoscope

Storage

This device is delivered sterilized and in order to maintain this state until its first use it must be kept in its original sterilized package.

Do not put any objects on the instrument or its package!
Do not store the instruments near aggressive chemical products!

Do not expose the instruments to direct or indirect sunlight or other ultra-violet rays!

Keep in dry area

Complaints will not be taken under consideration if the instruments have been stored improperly.

Cleaning and Sterilization

This device is designed and warranted for single use!

Risk in case of reuse

This device is for single use only. G-Flex denied all responsibility in case of reuse. The reuse of single use device might present risk to the safety of the user or patient, due to possible uncontrolled contamination and / or unreliability of the device performance.

User

The users of G-FLEX instruments must be specialists in their fields. An appropriate and specific training for preparation, care and maintenance of the flexible instruments is required.

Repairing and complains

In health protection of our employees only disinfected or sterilized instruments will be accepted for analysis or repairing. In any case disinfection or sterilization dates and validity must be labelled outside the package. If this requirement has not been fulfilled the instruments will be returned to the sender without analysis or repairing.

Legal foundation

The Law of the European Union applied.

Product support

In case of questions or difficulties concerning our instruments please contact your local distributor or G-FLEX directly during regular working hours.

Times available: Monday to Friday; 9am to 5pm (MET)

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G-Flex USA, Inc., 74 Southwood Pkwy Suite 600