EN - BIOPSY FORCEPS
FR - PINÇA À BIOPSIE
NL - BIOPSIETANG
ES - FORCEPS DE BIOPSIA
DE - BIOPSIEZANGE
PT - FÓRCEPS DE BIOPSIA
IT - PINZA PER BIOPSIA
EL - ΛΑΒΉΔΑ ΒΙΟΨΊΑΣ
AR - ملقط الخزعة
TR - BIYOPSI FORSEPSI

G-Flex Europe S.P.R.L. - 20, Rue de l'industrie 1400
Nivelles - Belgium | **Phone:** +32 67.88.36.65
Email: info@g-flex.com | **Website:** www.g-flex.com

Current Name: IFU-BiopsyForceps

Creation Date: 01/10/2018 **Version** 6

This document is only controlled in electronic format by the QM.

EN - INSTRUCTION FOR USE	3
FR - MODE D'EMPLOI	6
NL - GEBRUIKSAANWIJZING	9
ES - INSTRUCCIONES DE USO	12
DE - GEBRAUCHSANWEISUNG	15
PT - INSTRUÇÕES DE USO	19
IT - ISTRUZIONI PER L'USO	22
EL - ΟΔΗΓΊΕΣ ΧΡΉΣΗΣ	25
AR - العربية - تعليمات الاستخدام	28
TR - KULLANIM TALIMATLARI	32

EN - INSTRUCTION FOR USE

Intended Use

Biopsy Forceps are used for taking out biopsy specimens (tissues samples) from the human body, through endoscopes.

Product Description

The device is composed by a jaw (in different shapes and sizes), a spiral spring hose (in different lengths and diameters) and a operating handle.

Advices before first application

Please read carefully and follow all safety operating instructions and warnings before first application of the device.

A previous knowledge regarding handling and operation is required and essential.

Unpack the instrument carefully and examine the devices for any possibility of damage. In case of any damage or missing items contact your distributor immediately.

Considering that device are not delivered sterilized the user must do it before the first application (See Cleaning and Sterilization)

If the package present any damaged, the sterility of the device is not guarantee.

Handling and Operation

The device must be passed through the endoscope operating channel with its cups closed. Please observe compatibility between the biopsy forceps and the endoscope operating channel diameter.

Unwind the biopsy forceps, manipulate backward - forward smoothly the plastic spool of the handle to verify the opening / closing of the cups.

Biopsy Forceps are composed by a spiral and should never be manipulated with its catheter winded as this may damage the device and make its usage impossible!

Procedure

- Approach the biopsy forceps near to the tissue (lesion) to be taken from the patient's body.
- Open the cups by pushing the plastic spool at the handle.
- Place the forceps on the tissue and take the sample gently by pulling back the plastic spool (this manipulation will close the forceps cups).
- Keep the forceps cups closed and remove biopsy forceps out of the endoscope carefully.
- This tissue specimen taken will be ready for further histological examination.

Storage

Do not put any objects on the instrument or its package!
Do not store the instruments near aggressive chemical products!

Do not expose the instruments to direct or indirect sunlight or other ultra-violet rays!

Keep in dry area

Complaints will not be taken under consideration if the instruments have been stored improperly.

Cleaning and Sterilization

G-Flex recommend to follow one of the following sequences when sterilizing the device:

- Manual Cleaning, Ultrasonic Cleaning, Automatic/Thermal Cleaning and Autoclave; or
- Manual Cleaning, Ultrasonic Cleaning, Automatic/Thermal Cleaning and Ethylene Oxide Gas Sterilization

Either sterilize by Autoclave or by ETO. There is no need to do both.

Manual Cleaning

The cleaning of the device should be done immediately after each use:

- Immerse the device in a suitable liquid detergent or disinfectant (Always observe the minimum diameter of 20 cm when winding the instrument and applicability of the cleaning product!).
- Wipe the instrument with a soft towel.

Ultrasonic Cleaning

The ultrasonic cleaning facilitates the dislodging of residual material and must be done before sterilization. The ultrasonic cleaning must last for at least 30 minutes.

- Clean the device right after its use.
- Put the instrument in an ultrasonic cleaning tank always observing the winding minimum diameter 20 cm (Please observe the ultrasonic equipment manufacturer's instructions!).
- Use only tap water and the detergents which dissolve albumen, always according to ultrasonic equipment manufacturer's instructions.
- Rinse the instrument with tap water and wipe dry by using gauze pads.

Automatic / Thermal Cleaning

Alternatively or additionally to the manual cleaning, this instrument may also be cleaned in automatic / thermal cleaning machines by using approved alkaline detergents. Please refer to relevant instructions of the cleaning machine supplier.

Sterilization - Autoclave

Before autoclaving sterilization the instrument should be cleaned carefully as described in Manual Cleaning, Ultrasonic Cleaning and Automatic/Thermal Cleaning.

This symbol means that the product can be sterilized by autoclave according to the below parameters

Parameters:

Temperature: 134 °C / 273°F

Pressure: 3 bar

Exposure: ≥ 18 min

Also in here the minimum 20 cm diameter for winding the instrument must be noted.

The original packages of G-FLEX are not autoclavable!

Sterilization - Ethylene Oxide Gas Sterilization

Please follow manufacturer's warnings regarding to the biological indicators

Parameters:

Temperature: 57 °C

Pressure: 1,7 bar

Exposure: max. 4 h

Relative Humidity: 50 %

Gas Concentration: 12 %

Aeration Time: 7 days in ambient temperature or 12 hours by 50-57 °C

User

The users of G-FLEX instruments must be specialists in their fields. An appropriate and specific training for preparation, care and maintenance of the flexible instruments is required.

Repairing and complaints

In health protection of our employees only disinfected or sterilized instruments will be accepted for analysis or repairing. In any case disinfection or sterilization dates and validity must be labelled outside the package. If this requirement has not been fulfilled the instruments will be returned to the sender without analysis or repairing.

Legal foundation

The Law of the European Union applied.

Product support

In case of questions or difficulties concerning our instruments please contact your local distributor or G-FLEX directly during regular working hours.

Times available: Monday to Friday; 9am to 5pm (MET)

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