**EN - CLEANING BRUSH**  
**FR - BROSSE DE NETTOYAGE**  
**NL - REINIGINGSBORSTEL**  
**ES - CEPILLO DE LIMPIEZA**  
**DE - REINIGUNGSBÜRSTE**  
**PT - ESCOVA DE LIMPEZA**  
**IT - SPAZZOLA PER LA PULIZIA**  
**EL - ΒΟΥΡΤΣΑ ΚΑΘΑΡΙΣΜΟΥ**  
**AR - فرشاة التنظيف**  
**TR - TEMIZLIK FIRÇASI**

G-Flex Europe S.P.R.L. - 20, Rue de l'industrie 1400  
Nivelles - Belgium | Phone: +32 67.88.36.65  
Email: [info@g-flex.com](mailto:info@g-flex.com) | Website: [www.g-flex.com](http://www.g-flex.com)

**Current Name:** IFU-CleaningBrush

**Creation Date:** 22/10/2019 **Version** 10

*This document is only controlled in electronic format by the QM.*

EN - INSTRUCTION FOR USE .....	3
FR - MODE D'EMPLOI .....	6
NL - GEBRUIKSAANWIJZING .....	9
ES - INSTRUCCIONES DE USO .....	13
DE - GEBRAUCHSANWEISUNG .....	16
PT - INSTRUÇÕES DE USO .....	19
IT - ISTRUZIONI PER L'USO .....	23
EL - ΟΔΗΓΊΕΣ ΧΡΉΣΗΣ .....	26
AR - العربية - تعليمات الاستخدام .....	30
TR - KULLANIM TALIMATLARI .....	33

# EN - INSTRUCTION FOR USE

## Intended Use

The Cleaning Brush is meant to clean the biopsy channel of endoscopes.

## Product Description

The device is made of a long catheter (or spiral) with a brush on each side.

## Advices before first application

Please read carefully and follow all safety operating instructions and warnings before first application of the device.

A previous knowledge regarding handling and operation is required and essential.

Unpack the instrument carefully and examine the devices for any possibility of damage. In case of any damage or missing items contact your distributor immediately.

## Handling and Operation

The Cleaning Brush has to be passed through the full length of the endoscope biopsy channel several times with the solution described in your endoscope manufacturer notice.

## Procedure

1. Introduce the brush through one side of the biopsy channel and pull it out from the other side of the same.
2. Repeat this process until the working channel is free from dirt. (Recommended a minimum of 3 times)
3. Ensure to clean the brush before introducing it in the channel each time.

Please note that some models have the distal and proximal ends different from each other. In such models, make sure to introduce the brush with the thinnest diameter tip into the channel.

## Storage

Do not put any objects on the instrument or its package!  
Do not store the instruments near aggressive chemical products!

Do not expose the instruments to direct or indirect sunlight or other ultra-violet rays!

Keep in dry area

Complaints will not be taken under consideration if the instruments have been stored improperly.

## Cleaning and Sterilization

The disposable version of this device is designed and warranted for single use!

### **For the reusable version of this device, please consult the following**

G-Flex recommend to follow one of the following sequences when sterilizing the device:

- Manual Cleaning, Ultrasonic Cleaning, Automatic/Thermal Cleaning and Autoclave; or
- Manual Cleaning, Ultrasonic Cleaning, Automatic/Thermal Cleaning and Ethylene Oxide Gas Sterilization

Either sterilize by Autoclave or by ETO. There is no need to do both.

### **Manual Cleaning**

The cleaning of the device should be done immediately after each use:

- Immerse the device in a suitable liquid detergent or disinfectant (Always observe the minimum diameter of 20 cm when winding the instrument and applicability of the cleaning product!).
- Wipe the instrument with a soft towel.

### **Ultrasonic Cleaning**

The ultrasonic cleaning facilitates the dislodging of residual material and must be done before sterilization. The ultrasonic cleaning must last for at least 30 minutes.

- Clean the device right after its use.
- Put the instrument in an ultrasonic cleaning tank always observing the winding minimum diameter 20 cm (Please observe the ultrasonic equipment manufacturer's instructions!).
- Use only tap water and the detergents which dissolve albumen, always according to ultrasonic equipment manufacturer's instructions.
- Rinse the instrument with tap water and wipe dry by using gauze pads.

### **Automatic / Thermal Cleaning**

Alternatively or additionally to the manual cleaning, this instrument may also be cleaned in automatic / thermal cleaning machines by using approved alkaline detergents. Please refer to relevant instructions of the cleaning machine supplier.

### **Sterilization - Autoclave**

Before autoclaving sterilization the instrument should be cleaned carefully as described in Manual Cleaning, Ultrasonic Cleaning and Automatic/Thermal Cleaning.

This symbol means that the product can be sterilized by autoclave according to the below parameters

**Parameters:**

Temperature:	134 °C / 273°F
Pressure:	3 bar
Exposure:	≥ 18 min

Also in here the minimum 20 cm diameter for winding the instrument must be noted.

The original packages of G-FLEX are not autoclavable!

**Sterilization - Ethylene Oxide Gas Sterilization**

Please follow manufacturer's warnings regarding to the biological indicators

**Parameters:**

Temperature:	57 °C
Pressure:	1,7 bar
Exposure:	max. 4 h
Relative Humidity:	50 %
Gas Concentration:	12 %
Aeration Time:	7 days in ambient temperature or 12 hours by 50-57 °C

**Risk in case of reuse**

The disposable version of the device is for single use only. G-Flex denied all responsibility in case of reuse. The reuse of single use device might present risk to the safety of the user or patient, due to possible uncontrolled contamination and / or unreliability of the device performance.

**User**

The users of G-FLEX instruments must be specialists in their fields. An appropriate and specific training for preparation, care and maintenance of the flexible instruments is required.

**Repairing and complains**

In health protection of our employees only disinfected or sterilized instruments will be accepted for analysis or repairing. In any case disinfection or sterilization dates and validity must be labelled outside the package. If this requirement has not been fulfilled the instruments will be returned to the sender without analysis or repairing.

**Product disposal**

After use, these products may be a potential biohazard. Handle and dispose of in accordance with accepted medical practice and applicable local, state and federal laws and regulations.

**Legal foundation**

The Law of the European Union applied.

## Product support

In case of questions or difficulties concerning our instruments please contact your local distributor or G-FLEX directly during regular working hours.

**Times available:** Monday to Friday; 9am to 5pm (MET)

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