CYSTO GASTRO SET

SET DE GASTRO-QUISTE

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EN - INSTRUCTION FOR USE

Intended Use

Cysto-Gastro-Sets are intended to be used to puncture pancreatic cysts or pseudo-cysts endoscopically as an alternative to surgical or percutaneous treatments. During the procedure the cyst is punctured by a electrical high frequency needle and dilated to a diameter of 6, 8.5 or 10 Fr. (depending on which set is used).

Product Description

The device is made of an outer flexible teflon tube equipped with a metal tip at the distal end and an HF connector and the proximal end (all models). There is an inner catheter with a puncture needle at the distal end and both a Luer-lock and an HF connector and the proximal end (except on the 6Fr set).

Advices before first application

Please read carefully and follow all safety operating instructions and warnings before first application of the device.

A previous knowledge regarding handling and operation is required and essential. Unpack the instrument carefully and examine the devices for any possibility of damage. In case of any damage or missing items contact your distributor immediately.

Considering that device are not delivered sterilized the user must do it before the first application (See Cleaning and Sterilization)

If the package present any damaged, the sterility of the device is not guarantee.

Handling and Operation

Ensure that both components are assembled and fixed at the Luer-Lock connection. Please note: The inner catheter is approx. 20 cm longer than the outer Teflon tube. Both parts of the device are supplied with HF connectors that need to be fixed to and appropriate active cords.

Cysto-Gastro-Sets require a endoscope's working channel of at least 3.7 mm Ø for 10Fr Set, 3.2 mm Ø for the 8.5 Fr Set or 2.8 mm for the 6 Fr.

Use cut frequency only since the use or excess of coagulation frequency can damage the product. Recommended frequency is between 80 and 120watts.

Procedure

After the size and position of the cyst or pseudo-cyst is determined by endosonography, the device is inserted into the endoscope's operating channel.

- Connect the active cord to the HF connector at the proximal end of the inner catheter and into the HF generator.

- Carefully move out the puncture needle of the inner catheter (on 8.5 or 10 Fr sets only).
- Puncture the stomach or duodenal wall and the cyst with the puncture needle by using electric current (For the 6 Fr set use directly the outer tube as there is no inner catheter).
- Disconnect the electric current and remove the active cord from the inner catheter.
- If necessary, contrast can be injected through the lateral Luer-Lock adapter at the proximal end of the inner catheter (on 8.5 or 10 Fr sets only).
- Now connect the active cord to the lateral connector at the outer tube and to the HF generator.
- Loosen the Luer-Lock fixation between inner catheter and outer tube. Push forward the outer tube by using the inner catheter as a guiding catheter until the distal metal tip reaches the stomach wall.
- Switch on electricity and dilate the stomach's and cyst's opening in the wall by pushing forward the outer tube's metal tip carefully.
- Remove the active cord and leave the outer tube at its position.
- Terminate the HF -connection and withdraw the inner catheter completely (out of the outer tube and the endoscope).
- Place a guide wire (.025" or .035") through the outer tube into the cyst.
- Finally, remove the outer tube out of the working channel, while ensuring that the guide wire stays in position.
- Use the guide wire to place a Pigtail or Double Pigtail stent in order to drain the cyst, after the cystostomy.

Storage

This device is delivered sterilized and in order to maintain this state until its first use it must be kept in its original sterilized package.

Do not put any objects on the instrument or its package!
Do not store the instruments near aggressive chemical products!

Do not expose the instruments to direct or indirect sunlight or other ultra-violet rays!

Keep in dry area

Complaints will not be taken under consideration if the instruments have been stored improperly.

Cleaning and Sterilization

G-Flex recommend to follow one of the following sequences when sterilizing the device:

- Manual Cleaning, Ultrasonic Cleaning, Automatic/Thermal Cleaning and Autoclave; or
- Manual Cleaning, Ultrasonic Cleaning, Automatic/Thermal Cleaning and Ethylene Oxide Gas Sterilization

Either sterilize by Autoclave or by ETO. There is no need to do both.

Manual Cleaning

The cleaning of the device should be done immediately after each use:

- Immerse the device in a suitable liquid detergent or disinfectant (Always observe the minimum diameter of 20 cm when winding the instrument and applicability of the cleaning product!).
- Wipe the instrument with a soft towel.

Ultrasonic Cleaning

The ultrasonic cleaning facilitates the dislodging of residual material and must be done before sterilization. The ultrasonic cleaning must last for at least 30 minutes.

- Clean the device right after its use.
- Put the instrument in an ultrasonic cleaning tank always observing the winding minimum diameter 20 cm (Please observe the ultrasonic equipment manufacturer's instructions!).
- Use only tap water and the detergents which dissolve albumen, always according to ultrasonic equipment manufacturer's instructions.
- Rinse the instrument with tap water and wipe dry by using gauze pads.

Automatic / Thermal Cleaning

Alternatively or additionally to the manual cleaning, this instrument may also be cleaned in automatic / thermal cleaning machines by using approved alkaline detergents. Please refer to relevant instructions of the cleaning machine supplier.

Sterilization - Autoclave

Before autoclaving sterilization the instrument should be cleaned carefully as described in Manual Cleaning, Ultrasonic Cleaning and Automatic/Thermal Cleaning.

This symbol means that the product can be sterilized by autoclave according to the below parameters

Parameters:

Temperature: 134 °C / 273°F

Pressure: 3 bar

Exposure: ≥ 18 min

Also in here the minimum 20 cm diameter for winding the instrument must be noted.

The original packages of G-FLEX are not autoclavable!

Sterilization - Ethylene Oxide Gas Sterilization

Please follow manufacturer's warnings regarding to the biological indicators

Parameters:

Temperature:	57 °C
Pressure:	1,7 bar
Exposure:	max. 4 h
Relative Humidity:	50 %
Gas Concentration:	12 %
Aeration Time:	7 days in ambient temperature or 12 hours by 50-57 °C

User

The users of G-FLEX instruments must be specialists in their fields. An appropriate and specific training for preparation, care and maintenance of the flexible instruments is required.

Repairing and complains

In health protection of our employees only disinfected or sterilized instruments will be accepted for analysis or repairing. In any case disinfection or sterilization dates and validity must be labelled outside the package. If this requirement has not been fulfilled the instruments will be returned to the sender without analysis or repairing.

Legal foundation

The Law of the European Union applied.

Product support

In case of questions or difficulties concerning our instruments please contact your local distributor or G-FLEX directly during regular working hours.

Times available: Monday to Friday; 9am to 5pm (MET)

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G-Flex USA, Inc., 74 Southwood Pkwy Suite 600