**EN - DILATION BALLOON**  
**FR - BALLONNET DE DILATATION**  
**NL - DILATATIEBALLON**  
**ES - BALÓN DE DILATACIÓN**  
**DE - DILATIONSBALLON**  
**PT - BALÃO DE DILATAÇÃO**  
**IT - PALLONCINO DI DILATAZIONE**  
**EL - ΜΠΑΛΟΝΙ ΔΙΑΣΤΟΛΗΣ**  
**AR - بالون التوسيع**  
**TR - DİLATASYON BALONU**

G-Flex Europe S.P.R.L. - 20, Rue de l'industrie 1400  
Nivelles - Belgium | Phone: +32 67.88.36.65  
Email: [info@g-flex.com](mailto:info@g-flex.com) | Website: [www.g-flex.com](http://www.g-flex.com)

Current Name: IFU-DilationBalloon

Creation Date: 30/11/2018 Version 9

*This document is only controlled in electronic format by the QM.*

EN - INSTRUCTION FOR USE .....	3
FR - MODE D'EMPLOI .....	6
NL - GEBRUIKSAANWIJZING .....	9
ES - INSTRUCCIONES DE USO .....	13
DE - GEBRAUCHSANWEISUNG .....	16
PT - INSTRUÇÕES DE USO .....	20
IT - ISTRUZIONI PER L'USO .....	23
EL - ΟΔΗΓΪΕΣ ΧΡΪΗΣΗΣ .....	26
AR العربية - تعليمات الاستخدام .....	30
TR - KULLANIM TALIMATLARI .....	33

# EN - INSTRUCTION FOR USE

## Intended Use

The Dilation Balloon is used to dilate stenosis in gastrointestinal and/or bronchial tract areas via endoscopes. Models with balloons in 20 or 30 mm long are for biliary dilation, sphincteroplasty or broncho treatments, models with balloons in 55 mm long are for oesophageal, colonic or broncho treatments.

## Product Description

The device is composed of a double lumen teflon tube with a plastic connector at the proximal end in “Y” shape. This connector has two entries, one to inflate/deflate the balloon (Luer-lock type) and the other to insert a 0.035” guide wire. The dilation balloon is on the distal end.

## Advices before first application

Please read carefully and follow all safety operating instructions and warnings before first application of the device.

A previous knowledge regarding handling and operation is required and essential. Unpack the instrument carefully and examine the devices for any possibility of damage. In case of any damage or missing items contact your distributor immediately.

If the package present any damaged, the sterility of the device is not guarantee.

## Handling and Operation

- If damage or an open seal is detected, sterilization may be compromised and the instrument should be discarded and replaced.
- Keep an appropriate inflation device with pressure gauge available.
- Check the inflation device for appropriate pressure and volume.
- Keep enough fill liquid (e. g. sterile low viscosity contrast media; sterile, isotonic saline solution; or a mixture of both) available.
- To support application with a guide wire, keep an appropriate guide wire available.
- Open packaging and carefully remove the pouch.
- Inspect catheter to any damages. If you note any problems, exchange the balloon for a new one.
- Connect a syringe to the Luer-Lock and apply a moderate vacuum and remove protective tubing.
- To improve the slide properties inside the endoscope, spray the balloon with silicone. Only use silicone, other lubricants may damage the product.

## Procedure

### Filling

- To remove the air from the balloon, please follow the instructions below:
- Connect an inflation device to the Luer-Lock and hold the catheter tip down.

- Apply vacuum and keep for at least 30 seconds.
- Lock the stop-lock lever to keep both balloon and catheter free from air.
- Disconnect the inflation device from the balloon and fill it completely with liquid.
- Recommended liquids are sterile, isotonic saline solution; Sterile low viscosity contrast media or a mixture of both.
- Keep vacuum during insertion in the endoscope

 Essential volume of air in the system may result in inconsistent expansion of the balloon. If this is recognized, repeat filling steps.

 Allowing air or liquid to fill the balloon before entering the working channel of the endoscope may prevent the device to be inserted and therefore the procedure can not be completed.

## **Application**

- Slowly introduce the balloon catheter into the endoscope's working channel. It can be used over a 0,035"guide wire.
- Due to different endoscope specifications and types, there may be higher resistance when initially inserting the catheter into the endoscope and approx. 2 to 3 cm from the working channel's distal end.
- Advance the balloon under endoscopic control to its final position.
- Fill the balloon using an appropriate inflating device.

 Never fill the balloon with any gas!

 Never exceed the given maximum inflation pressure. In case of a balloon burst or a liquid leakage, empty the balloon and remove it carefully together with the endoscope. Don't try to pull burst balloon back into the endoscope. Start with a new balloon again.

- After reaching the indicated balloon diameter adjust the balloon pressure as necessary (do not exceed maximum inflation pressure). After the inflator device has been set to working pressure the pressure gauge will indicate slightly decrease of pressure because of pressure equalization in the catheter system. After a few seconds the indication should stabilize.
- Keep pressure till the requested dilation is reached. The pressure indication may fluctuate during dilation.
- After dilation deflate the balloon slowly. This may take approx. 20 to 30 seconds depending on the size and liquid.
- To deflate completely the proximal end of the balloon must be visible in the endoscope during vacuum application. Remove the catheter only if completely empty.
- Keep the endoscope as straight as possible. Every curve increases wall friction and difficulty to remove the catheter.
- Pull the catheter slowly out of the endoscope. If an unusual resistance is recognized, remove endoscope and balloon together to avoid patient injury or damages on the endoscope.
- After use, dispose the entire instrument according to legal requirements and your facility's infection control protocols.

## Storage

This device is delivered sterilized and in order to maintain this state until its first use it must be kept in its original sterilized package.

Do not put any objects on the instrument or its package!  
Do not store the instruments near aggressive chemical products!

Do not expose the instruments to direct or indirect sunlight or other ultra-violet rays!

Keep in dry area

Complaints will not be taken under consideration if the instruments have been stored improperly.

## Possible complication

Perforation, hemorrhage, hematoma, septicemia/infection, allergic reaction to contrast medium

## Cleaning and Sterilization

This device is designed and warranted for single use!

## Risk in case of reuse

This device is for single use only. G-Flex denied all responsibility in case of reuse. The reuse of single use device might present risk to the safety of the user or patient, due to possible uncontrolled contamination and / or unreliability of the device performance.

## User

The users of G-FLEX instruments must be specialists in their fields. An appropriate and specific training for preparation, care and maintenance of the flexible instruments is required.

## Repairing and complaints

In health protection of our employees only disinfected or sterilized instruments will be accepted for analysis or repairing. In any case disinfection or sterilization dates and validity must be labelled outside the package. If this requirement has not been fulfilled the instruments will be returned to the sender without analysis or repairing.

## Legal foundation

The Law of the European Union applied.

## Product support

In case of questions or difficulties concerning our instruments please contact your local distributor or G-FLEX directly during regular working hours.

**Times available:** Monday to Friday; 9am to 5pm (MET)

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