EN - DISPOSABLE BIOPSY FORCEPS

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NL - WEGWERP-BIOPSIETANG

ES - FORCEPS DE BIOPSIA DESECHABLE

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Current Name: IFU-DisposableBiopsyForceps

Creation Date: 01/10/2018 **Version** 9

This document is only controlled in electronic format by the QM.

EN - INSTRUCTION FOR USE	3
FR - MODE D'EMPLOI	5
NL - GEBRUIKSAANWIJZING	7
ES - INSTRUCCIONES DE USO	9
DE - GEBRAUCHSANWEISUNG	11
PT - INSTRUÇÕES DE USO	13
IT - ISTRUZIONI PER L'USO	15
EL - ΟΔΗΓΪΕΣ ΧΡΪΗΣΗΣ	17
AR - العربية - تعليمات الاستخدام	19
TR - KULLANIM TALIMATLARI	21

EN - INSTRUCTION FOR USE

Intended Use

Disposable Biopsy Forceps are used for taking out biopsy specimens (tissues samples) from the human body, through endoscopes.

Product Description

The device is composed by a jaw (in different shapes and sizes), a spiral spring hose (in different lengths and diameters) and an operating handle.

Advices before first application

Please read carefully and follow all safety operating instructions and warnings before first application of the device.

A previous knowledge regarding handling and operation is required and essential.

Unpack the instrument carefully and examine the devices for any possibility of damage. In case of any damage or missing items contact your distributor immediately.

If the package present any damaged, the sterility of the device is not guarantee.

Handling and Operation

The biopsy forceps must be passed through the endoscope operating channel with its cups closed. Please observe compatibility between the biopsy forceps and the endoscope operating channel diameter.

Unwind the disposable biopsy forceps, manipulate backward - forward the plastic spool of the handle to verify the opening / closing of the cups.

Procedure

- Approach the biopsy forceps near to the tissue (lesion) to be taken from the patient's body.
- Open the cups by pushing the plastic spool at the handle.
- Place the forceps on the tissue and take the sample gently by pulling back the plastic spool (this manipulation will close the forceps cups).
- Keep the forceps cups closed and remove biopsy forceps out of the endoscope carefully.
- This tissue specimen taken will be ready for further histological examination.

Storage

This device is delivered sterilized and in order to maintain this state until its first use it must be kept in its original sterilized package.

Do not put any objects on the instrument or its package!
Do not store the instruments near aggressive chemical products!

Do not expose the instruments to direct or indirect sunlight or other ultra-violet rays!

Keep in dry area

Complaints will not be taken under consideration if the instruments have been stored improperly.

Cleaning and Sterilization

This device is designed and warranted for single use!

Risk in case of reuse

This device is for single use only. G-Flex denied all responsibility in case of reuse. The reuse of single use device might present risk to the safety of the user or patient, due to possible uncontrolled contamination and / or unreliability of the device performance.

User

The users of G-FLEX instruments must be specialists in their fields. An appropriate and specific training for preparation, care and maintenance of the flexible instruments is required.

Repairing and complaints

In health protection of our employees only disinfected or sterilized instruments will be accepted for analysis or repairing. In any case disinfection or sterilization dates and validity must be labelled outside the package. If this requirement has not been fulfilled the instruments will be returned to the sender without analysis or repairing.

Legal foundation

The Law of the European Union applied.

Product support

In case of questions or difficulties concerning our instruments please contact your local distributor or G-FLEX directly during regular working hours.

Times available: Monday to Friday; 9am to 5pm (MET)

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