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# EN - INSTRUCTION FOR USE

## Intended Use

Disposable Hemoclip are intended to place clip within the Gastrointestinal (GI) tract for the following purposes:

- Endoscopic marking;
- Hemostasis:
  - Mucosal/sub-mucosal defects (< 3 cm)
  - Bleeding ulcers
  - Arteries (< 2 mm)
  - Polyps (< 1.5 cm in diameter)
  - Diverticula in the colon
- Anchor to attached jejunal feeding tubes to the wall of the small bowel;
- As a supplementary method to close luminal perforations (< 20 mm) in the GI tract that can be treated conservatively

## Product Description

The device is made of a pre-loaded, rotatable, radiopaque, single-use clip in a flexible delivery system. There is a handle in the proximal part to manipulate the rotation and the deployment, a long flexible tube and the clip at the distal part. In the U2525 series the handle is also used for opening and closing of the clip. It is designed to work in gastroscope and colonoscopes (see model working length) with a working channel equal to or greater than 2.8 mm diameter (See picture on page 3).

## Advices before first application



Please read carefully and follow all safety operating instructions and warnings before first application of the device.

A previous knowledge regarding handling and operation is required and essential.

Unpack the instrument carefully and examine the devices for any possibility of damage. In case of any damage or missing items contact your distributor immediately.

If the package present any damaged, the sterility of the device is not guarantee.

## Handling and Operation

Ensure the device and its packaging are intact before use. Remove the device from its packaging before use.

Please observe compatibility between the operating channel and device diameter.

U2025 Series - Before inserting the device into the biopsy channel ensure that the clip is inside the over-sheath tube.

U2525 Series - Before inserting the device into the biopsy channel ensure that the clip is closed.

The U2025 Series of clip can NOT be opened and closed prior to final deployment.


The U2525 Series of clip can be opened and closed **no more** than 5 times prior to final deployment.

## Procedure

### Insertion into the Biopsy Channel



- Insert the device into the biopsy channel until it comes out of the distal part of the endoscope.
- **U2025 Series** - Once under direct visualization, remove the stopper located between the handle and the over-sheath grip (see picture below).


- **U2525 Series** - Once under direct visualization, you are ready for the next step

## Deployment of the clip

- **U2025 Series** - To deploy the clip hold the over-sheath grip with one hand and push the handle toward the over-sheath grip with the other hand (see picture below). The clip is then deployed open.


- **U2525 Series** - The clip is deployed closed. To open the clip, move the slider distally (away from the thumb ring). The opening and closing can be done no more than 5 times prior to final deployment (see picture below).


## Rotation of the clip

- **U2025 Series** - To rotate the clip, in order to adjust the angle between the clip and the lesion, move the slider distally, hold the rotating seat with one hand and turn the complete handle in any direction with the other hand (see picture below).

- **U2525 Series** - To rotate the clip, in order to adjust the angle between the clip and the lesion, just rotate the complete handle in any direction (See picture below)


## Closing the clip

### • U2025 Series


- To close the clip jaws on the desired location move the slider proximally (close to the thumb ring) until the first tactile resistance is felt in the handle (see picture below).

**After this resistance point, the clip will be permanently deployed.**


- To permanently deploy the clip, continue moving the slider proximally (toward the thumb ring) beyond the second tactile resistance point. At this point, a first click may be heard or felt. Continue moving the slider toward the thumb ring until a second click and/or tactile resistance point is heard or felt. Continue moving the slider until it reaches the thumb ring (see picture below)


## • U2525 Series

- To close the clip jaws on the desired location move the slider proximally (close to the thumb ring) until the first tactile resistance is felt in the handle (see picture below).

- To permanently deploy the clip, continue moving the slider proximally (toward the thumb ring) beyond the tactile resistance point until the Hole in the slider reaches the Locking System. At this point, click may be heard or felt. The Locking System ensures that the tracking wire stays inside the spring tube (spiral) and does not damage the endoscope biopsy channel. (see picture below).

## Removal of the Delivery System

- Remove the whole delivery system from the biopsy channel and check the result with the endoscope, if need repeat the whole procedure with other clips.

## Storage

This device is delivered sterilized and in order to maintain this state until its first use it must be kept in its original sterilized package.

Do not put any objects on the instrument or its package!  
Do not store the instruments near aggressive chemical products!

Do not expose the instruments to direct or indirect sunlight or other ultra-violet rays!

Keep in dry area

Complaints will not be taken under consideration if the instruments have been stored improperly.

## Cleaning and Sterilization

This device is designed and warranted for single use!

## Risk in case of reuse

This device is for single use only. G-Flex denied all responsibility in case of reuse. The reuse of single use device might present risk to the safety of the user or patient, due to possible uncontrolled contamination and / or unreliability of the device performance.

## Contraindications

- Do not use this device when hemostasis cannot be verified visually with an endoscopic field of view.
- Do not perform MRI procedures on patients who have had clips placed within their gastrointestinal tract, as this could be harmful to the patient.
- Arteries greater than 2 mm.
- Polyps greater than 1.5 cm in diameter.
- Mucosal / Sub- Mucosal defects greater than 3 cm.

## User

The users of G-FLEX instruments must be specialists in their fields. An appropriate and specific training for preparation, care and maintenance of the flexible instruments is required.

## Repairing and complains

In health protection of our employees only disinfected or sterilized instruments will be accepted for analysis or repairing. In any case disinfection or sterilization dates and validity must be labelled outside the package. If this requirement has not been fulfilled the instruments will be returned to the sender without analysis or repairing.

## **Legal foundation**

The Law of the European Union applied.

## **Product support**

In case of questions or difficulties concerning our instruments please contact your local distributor or G-FLEX directly during regular working hours.

**Times available:** Monday to Friday; 9am to 5pm (MET)

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G-Flex USA, Inc., 74 Southwood Pkwy Suite 600