EN - DISPOSABLE POLYPECTOMY SNARES

FR - ANSE À POLYPECTOMIE JETABLE

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DESECHABLES

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PT - ALÇA PARA POLIPECTOMIA

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IT - ANSE DA POLIPECTOMIA MONOUSO

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Nivelles - Belgium | Phone: +32 67.88.36.65

Email: info@g-flex.com | Website: www.g-flex.com

Current Name: IFU-Disposable Polypectomy Snares **Creation Date:** 01/10/2018 **Version** 9

This document is only controlled in electronic format by the QM.

EN - INSTRUCTION FOR USE	3
FR - MODE D'EMPLOI	5
NL - GEBRUIKSAANWIJZING	7
ES - INSTRUCCIONES DE USO	10
DE - GEBRAUCHSANWEISUNG	12
PT - INSTRUÇÕES DE USO	15
IT - ISTRUZIONI PER L'USO	17
EL - ΟΔΗΓΊΕΣ ΧΡΉΣΗΣ	19
AR - العربية - تعليمات الاستخدام	22
TR - KULLANIM TALIMATLARI	25

EN - INSTRUCTION FOR USE

Intended Use

Disposable Polypectomy Snares are intended to be used to remove and or cauterize any kind of polyps in the gastrointestinal tract.

Product Description

Disposable Polypectomy Snares are made of a flexible tube with a loop in different sizes and shapes at the distal end and a handle to expose and retract the loop into the catheter. The handle has a HF connector to be attached to a high frequency (HF) generator and an active cord.

Advices before first application

Please read carefully and follow all safety operating instructions and warnings before first application of the device.

A previous knowledge regarding handling and operation is required and essential. Unpack the instrument carefully and examine the devices for any possibility of damage. In case of any damage or missing items contact your distributor immediately.

If the package present any damaged, the sterility of the device is not guarantee.

Handling and Operation

Insertion into the Endoscope

The Disposable Polypectomy Snare must be passed through the endoscope operating channel with its snare withdrawn into the tube. Please observe compatibility between the Disposable Polypectomy Snare and the operating channel diameters.

HF active cord Connection

Connect the active cord to electricity connector at the handle and to the HF-generator.

HF Generators

Disposable Polypectomy Snares can be connected to all licensed HF monopolar generators by using the active cord.

Always check compatibility between HF generator and active cord. In doubt contact the manufacturer of generator or your local distributor.

Polypectomy snares are composed by a catheter and should never be manipulated with its catheter winded as this may damage the device and make its usage impossible!

Procedure

- Place the distal tip of the tube right in front of the polyp to be snared-cut.
- By pushing the handle spool expose the snare fully out of the tube and capture the polyp into the open snare.
- Pull back the spool until the snare closes on the polyp stalk and stays tighten around the polyp pedicle.
- Apply HF-electric current on the snare (Please observe HF generator manufacturer's instructions!).
- Cut the polyp by pulling the snare completely into the tube.
- Stop the HF-electric current and remove carefully the Disposable Polypectomy Snare out of the endoscope.

Do not touch other tissue areas with the snare while using HF-electric current (high hurt risk for the patient)!

Storage

This device is delivered sterilized and in order to maintain this state until its first use it must be kept in its original sterilized package.

Do not put any objects on the instrument or its package!
Do not store the instruments near aggressive chemical products!

Do not expose the instruments to direct or indirect sunlight or other ultra-violet rays!

Keep in dry area

Complaints will not be taken under consideration if the instruments have been stored improperly.

Cleaning and Sterilization

This device is designed and warranted for single use!

Risk in case of reuse

This devices is for single use only. G-Flex denied all responsibility in case of reuse. The reuse of single use device might present risk to the safety of the user or patient, due to possible uncontrolled contamination and / or unreliability of the device performance.

User

The users of G-FLEX instruments must be specialists in their fields. An appropriate and specific training for preparation, care and maintenance of the flexible instruments is required.

Repairing and complaints

In health protection of our employees only disinfected or sterilized instruments will be accepted for analysis or repairing. In any case disinfection or sterilization dates and validity must be labelled outside the package. If this requirement has not been fulfilled the instruments will be returned to the sender without analysis or repairing.

Legal foundation

The Law of the European Union applied.

Product support

In case of questions or difficulties concerning our instruments please contact your local distributor or G-FLEX directly during regular working hours.

Times available: Monday to Friday; 9am to 5pm (MET)

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