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Current Name: IFU-Endoscopic Drainage Stent and Application System **Creation Date:** 30/11/2018 **Version:** 10

This document is only controlled in electronic format by the QM.

EN - INSTRUCTION FOR USE	3
FR - MODE D'EMPLOI	7
NL - GEBRUIKSAANWIJZING.....	11
ES - INSTRUCCIONES DE USO	16
DE - GEBRAUCHSANWEISUNG	20
PT - INSTRUÇÕES DE USO	25
IT - ISTRUZIONI PER L'USO	29
EL - ΟΔΗΓΊΕΣ ΧΡΉΣΗΣ	34
AR العربية - تعليمات الاستخدام.....	38
TR - KULLANIM TALIMATLARI	43

EN - INSTRUCTION FOR USE

Intended Use

Endoscopic Drainage Stents and Application System are intended to place stent via endoscopes and drain biliary, cysto or pancreatic obstructions.

Product Description

The device is made of a radiopaque stent with or without flaps and may or not come pre-mounted with its application device, which is made of an inner guiding catheter with stainless steel marker at the distal end and a Y connector at the proximal end. They also have an outer pusher with stainless steel marker at the distal end and a Luer-lock adapter at the proximal end. If the application device is not present with your stent, please consider having one at hand. (Refer to the current catalogue for more information)

Advices before first application

Please read carefully and follow all safety operating instructions and warnings before first application of the device.

A previous knowledge regarding handling and operation is required and essential.

Unpack the instrument carefully and examine the devices for any possibility of damage. In case of any damage or missing items contact your distributor immediately.

If the package present any damaged, the sterility of the device is not guarantee.

Handling and Operation

The Endoscopic Drainage Stents are placed in conjunction with an endoscope, a guide wire, a pusher and, when necessary, a guiding catheter

To facilitate the application, an endoscopic sphincterotomy prior to guide wire placement is recommended. (Please observe HF-generator and sphincterotome manufacturer's instructions).

Procedure

For the Endoscopic Drainage Stent, please proceed as follows:

- When implanting stents with a diameter of more than 7 Fr., it is strongly recommended to apply a guiding catheter. Introduce this catheter over the guide wire and position it on top of the obstructed area.
- Choose type and length of stent to be placed.
- Use the included positioning sleeve for easier introduction by sliding it onto the proximal end of the stent (the rear flap should fit back into it).
- Introduce the device (with positioning sleeve) and the pusher over the proximal end of the guide wire (guiding catheter) into the operating channel. The positioning sleeve does not fit into the channel and will stay outside.

- Place the device by using the pusher and carefully pass the obstructed area (under constant fluoroscopic control).
- Pull the guide wire (and guiding catheter) back into the operating channel while maintaining the position of the stent with the help of the pusher. After removing guide wire, guiding catheter and pusher the Stent will return to its original shape.

For the Application System, please proceed as follows:

- Use the included positioning sleeve for easier introducing by sliding it to the proximal end of the stent. So the rear flap will fit back into it.
- Ensure that the Luer-Lock connection between pusher and guiding catheter is firmly closed. Loosen the sealing cap at the Y-adapter.
- Place the biliary stent on to the device by sliding it over the distal end of guiding catheter up to the pusher.
- Now introduce the loaded device over the proximal end of the guide wire into the operating channel. Positioning sleeve does not fit to the channel and will stay outside.
- Advance the device under constant fluoroscopic control carefully until the radiopaque metal tip of the guiding catheter passes the obstructed area.
- Closing the sealing cap at the Y-adapter prevents movement of the guide wire inside of the guiding catheter. Contrast injection is possible through the Luer-Lock-Port at the Y-adapter.
- Open up the Luer-Lock connection between guiding catheter and pusher and advance the biliary stent to its final position by moving forward the pusher carefully.

Keep guiding catheters position while placing the stent

- Pull back the guiding catheter and guide wire into the operating channel. Prevent movement of the stent by keeping the pusher in its last position. After that the stent will return to its original shape.
- Finally, remove pusher, guiding catheter and guide wire out of the duodenoscope.

Injection of contrast medium into the biliary duct is possible via the guiding catheter by using the included Luer-Lock adapter.

Storage

This device is delivered sterilized and in order to maintain this state until its first use it must be kept in its original sterilized package.

Do not put any objects on the instrument or its package!
Do not store the instruments near aggressive chemical products!

Do not expose the instruments to direct or indirect sunlight or other ultra-violet rays!

Keep in dry area

Complaints will not be taken under consideration if the instruments have been stored improperly.

Cleaning and Sterilization

Most of these products are designed and warranted for single use!

Some pushers are reusable, for those, please consult the following

G-Flex recommend to follow one of the following sequences when sterilizing the device:

- Manual Cleaning, Ultrasonic Cleaning, Automatic/Thermal Cleaning and Autoclave; or
- Manual Cleaning, Ultrasonic Cleaning, Automatic/Thermal Cleaning and Ethylene Oxide Gas Sterilization

Either sterilize by Autoclave or by ETO. There is no need to do both.

Manual Cleaning

The cleaning of the device should be done immediately after each use:

- Immerse the device in a suitable liquid detergent or disinfectant (Always observe the minimum diameter of 20 cm when winding the instrument and applicability of the cleaning product!).
- Wipe the instrument with a soft towel.

Ultrasonic Cleaning

The ultrasonic cleaning facilitates the dislodging of residual material and must be done before sterilization. The ultrasonic cleaning must last for at least 30 minutes.

- Clean the device right after its use.
- Put the instrument in an ultrasonic cleaning tank always observing the winding minimum diameter 20 cm (Please observe the ultrasonic equipment manufacturer's instructions!).
- Use only tap water and the detergents which dissolve albumen, always according to ultrasonic equipment manufacturer's instructions.
- Rinse the instrument with tap water and wipe dry by using gauze pads.

Automatic / Thermal Cleaning

Alternatively or additionally to the manual cleaning, this instrument may also be cleaned in automatic / thermal cleaning machines by using approved alkaline detergents. Please refer to relevant instructions of the cleaning machine supplier.

Sterilization - Autoclave

Before autoclaving sterilization the instrument should be cleaned carefully as described in Manual Cleaning, Ultrasonic Cleaning and Automatic/Thermal Cleaning.

This symbol means that the product can be sterilized by autoclave according to the below parameters

Parameters:

Temperature: 134 °C / 273°F

Pressure: 3 bar

Exposure: ≥ 18 min

Also in here the minimum 20 cm diameter for winding the instrument must be noted.

The original packages of G-FLEX are not autoclavable!

Sterilization - Ethylene Oxide Gas Sterilization

Please follow manufacturer's warnings regarding to the biological indicators

Parameters:

Temperature: 57 °C

Pressure: 1,7 bar

Exposure: max. 4 h

Relative Humidity: 50 %

Gas Concentration: 12 %

Aeration Time: 7 days in ambient temperature or 12 hours by 50-57 °C

Risk in case of reuse

Most of these devices are for single use only. G-Flex denied all responsibility in case of reuse. The reuse of single use device might present risk to the safety of the user or patient, due to possible uncontrolled contamination and / or unreliability of the device performance.

User

The users of G-FLEX instruments must be specialists in their fields. An appropriate and specific training for preparation, care and maintenance of the flexible instruments is required.

Repairing and complains

In health protection of our employees only disinfected or sterilized instruments will be accepted for analysis or repairing. In any case disinfection or sterilization dates and validity must be labelled outside the package. If this requirement has not been fulfilled the instruments will be returned to the sender without analysis or repairing.

Legal foundation

The Law of the European Union applied.

Product support

In case of questions or difficulties concerning our instruments please contact your local distributor or G-FLEX directly during regular working hours.

Times available: Monday to Friday; 9am to 5pm (MET)

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