**MUCOSECTOMY SNARES**  
**ANSE À MUCOSECTOMIE**  
**MUCOSECTOMIE-SNAAR**  
**ASAS DE MUCOSECTOMÍA**

EN - INSTRUCTION FOR USE .....	2
FR - MODE D'EMPLOI .....	4
NL - GEBRUIKSAANWIJZING .....	6
ES - INSTRUCCIONES DE USO .....	8

**G-Flex Europe S.P.R.L.** - 20, Rue de l'industrie 1400  
Nivelles - Belgium | **Phone:** +32 67.88.36.65  
**Email:** [info@g-flex.com](mailto:info@g-flex.com) | **Website:** [www.g-flex.com](http://www.g-flex.com)  
**US Commercial Office Only:**  
G-Flex USA, Inc., 74 Southwood Pkwy Suite 600  
Atlanta, GA 30354 – USA

 **0120**

# EN - INSTRUCTION FOR USE

## Intended Use

Mucosectomy Snares are intended to be used with high frequency electric current to cut polyps in the gastrointestinal tract.

## Product Description

Mucosectomy Snares are made of a flexible tube with a loop in different size and shape at the distal end and a handle to expose and retract the loop into the catheter. The handle has a HF connector to be attached to a high frequency (HF) generator and an active cord.

## Advices before first application

Please read carefully and follow all safety operating instructions and warnings before first application of the device.

A previous knowledge regarding handling and operation is required and essential. Unpack the instrument carefully and examine the devices for any possibility of damage. In case of any damage or missing items contact your distributor immediately.

If the package present any damaged, the sterility of the device is not guarantee.

## Handling and Operation

### Insertion into the Endoscope

The Mucosectomy Snares must be passed through the endoscope operating channel with its snare withdrawn into the tube. Please observe compatibility between the Mucosectomy Snares and the operating channel diameters.

### HF active cord Connection

Connect the active cord to electricity connector at the handle and to the HF-generator.

### HF Generators

Mucosectomy Snares can be connected to all licensed HF generators by using the active cord.

Always check compatibility between HF generator and active cord. In doubt contact the manufacturer of generator or your local distributor.

## Procedure

- Place the distal tip of the tube right in front of the polyp to be snared.
- By pushing the handle spool expose the snare fully out of the tube and capture the polyp into the open snare.
- Pull back the spool until the snare closes on the polyp stalk and stays tighten around the polyp pedicle.

- Apply HF-electric current on the snare (Please observe HF generator manufacturer's instructions!).
- Cut the polyp by pulling the snare completely into the tube.
- Stop the HF-electric current and remove carefully the Mucosectomy Snares out of the endoscope.

Do not touch other tissue areas with the snare while using HF-electric current (high hurt risk for the patient)!

## Storage

This device is delivered sterilized and in order to maintain this state until its first use it must be kept in its original sterilized package.

Do not put any objects on the instrument or its package!  
Do not store the instruments near aggressive chemical products!

Do not expose the instruments to direct or indirect sunlight or other ultra-violet rays!

Keep in dry area

Complaints will not be taken under consideration if the instruments have been stored improperly.

## Cleaning and Sterilization

This device is designed and warranted for single use!

## Risk in case of reuse

This device is for single use only. G-Flex denied all responsibility in case of reuse. The reuse of single use device might present risk to the safety of the user or patient, due to possible uncontrolled contamination and / or unreliability of the device performance.

## User

The users of G-FLEX instruments must be specialists in their fields. An appropriate and specific training for preparation, care and maintenance of the flexible instruments is required.

## Repairing and complains

In health protection of our employees only disinfected or sterilized instruments will be accepted for analysis or repairing. In any case disinfection or sterilization dates and validity must be labelled outside the package. If this requirement has not been fulfilled the instruments will be returned to the sender without analysis or repairing.

## Legal foundation

The Law of the European Union applied.

## Product support

In case of questions or difficulties concerning our instruments please contact your local distributor or G-FLEX directly during regular working hours.

**Times available:** Monday to Friday; 9am to 5pm (MET)

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