NASAL BILIARY DRAINAGE DRAIN BILIAIRE NASAL DRENAJE BILIAR NASAL

| | |
|---------------------------------|---|
| EN - INSTRUCTION FOR USE | 2 |
| FR - MODE D'EMPLOI | 4 |
| ES - INSTRUCCIONES DE USO | 6 |

G-Flex Europe S.P.R.L. - 20, Rue de l'industrie 1400
Nivelles - Belgium | **Phone:** +32 67.88.36.65
Email: info@g-flex.com | **Website:** www.g-flex.com
US Commercial Office Only:
G-Flex USA, Inc., 74 Southwood Pkwy Suite 600
Atlanta, GA 30354 – USA

CE 0120

EN - INSTRUCTION FOR USE

Intended Use

Nasal Biliary Drainage are intended for the extracorporeal diversion of the biliary and pancreatic ducts. They provide an effective drainage and thereby reduce the risk of cholangitis. They are mostly used as a pre-operative preparation.

Product Description

The device is made of a radiopaque nasal biliary probe, a nasal catheter, a drainage connection tube and a Luer-lock adapter.

Advices before first application

Please read carefully and follow all safety operating instructions and warnings before first application of the device.

A previous knowledge regarding handling and operation is required and essential. Unpack the instrument carefully and examine the devices for any possibility of damage. In case of any damage or missing items contact your distributor immediately.

If the package present any damaged, the sterility of the device is not guarantee.

Handling and Operation

The insertion of the device is done via the operating channel of a duodenoscope. The usage of a guide wire (not included) during the insertion is necessary.

Procedure

- Place the distal end of the duodenoscope in front of the papilla.
- Pass the guide wire through the duodenoscope's operating channel into the biliary duct, up to the common hepatic duct.
- Slide the Nasal Biliary Probe over the guide wire with its distal tip straightened.
- Place the Nasal Biliary Probe into the common hepatic duct (under fluoroscopic control).
- Carefully remove the duodenoscope and ascertain that the Nasal Biliary Drainage Probe was left in place (under fluoroscopic control).
- Gently remove the guide wire from the probe.
- Insert the Nasal Catheter through the patient's nose and direct it into the pharynx.
- Grasp the Nasal Catheter's end with a forceps and remove it through the patient's mouth.
- The proximal end of the Biliary Drainage Probe, which is protruding from the patient's mouth, is led out of the nose via the Nasal Catheter.
- Remove the Nasal Catheter.
- Connect the enclosed Luer-Lock adapter to the proximal end of the Biliary Drainage Probe.

- Fix the end of the Drainage Probe on the patient's cheek and forehead applying adhesive bandage.
- The probe can now be connected to a receptacle (using the drainage connection tube).

Storage

This device is delivered sterilized and in order to maintain this state until its first use it must be kept in its original sterilized package.

Do not put any objects on the instrument or its package!
Do not store the instruments near aggressive chemical products!

Do not expose the instruments to direct or indirect sunlight or other ultra-violet rays!

Keep in dry area

Complaints will not be taken under consideration if the instruments have been stored improperly.

Cleaning and Sterilization

This device is designed and warranted for single use!

Risk in case of reuse

This device is for single use only. G-Flex denied all responsibility in case of reuse. The reuse of single use device might present risk to the safety of the user or patient, due to possible uncontrolled contamination and / or unreliability of the device performance.

User

The users of G-FLEX instruments must be specialists in their fields. An appropriate and specific training for preparation, care and maintenance of the flexible instruments is required.

Repairing and complains

In health protection of our employees only disinfected or sterilized instruments will be accepted for analysis or repairing. In any case disinfection or sterilization dates and validity must be labelled outside the package. If this requirement has not been fulfilled the instruments will be returned to the sender without analysis or repairing.

Legal foundation

The Law of the European Union applied.

Product support

In case of questions or difficulties concerning our instruments please contact your local distributor or G-FLEX directly during regular working hours.

Times available: Monday to Friday; 9am to 5pm (MET)

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