EN - NON-VASCULAR GUIDEWIRES
FR - FILS GUIDE NON VASCULAIRES
NL - NIET-VASCULAIRE VOERDRADEN
ES - ALAMBRES GUÍA NO VASCULAR
DE - NICHT-VASKULÄRE
FÜHRUNGSDRÄHTE
PT - FIOS-GUIA NÃO VASCULARES
IT - FILI GUIDA NON VASCOLARI
EL - ΟΔΗΓΉΟ ΣΎΡΜΑ ΜΗ ΑΓΓΕΙΑΚΉΣ
ΠΡΟΣΠΉΛΑΣΗΣ
AR - أسلاك التوجيه اللاوعائية
TR - NONVASKÜLER KILAVUZ
TELLERİ

G-Flex Europe S.P.R.L. - 20, Rue de l'industrie 1400
Nivelles - Belgium | Phone: +32 67.88.36.65
Email: info@g-flex.com | Website: www.g-flex.com

Current Name: IFU-NonVascularGuideWire

Creation Date: 01/10/2018

Version 8

This document is only controlled in electronic format by the QM.

EN - INSTRUCTION FOR USE	3
FR - MODE D'EMPLOI	6
NL - GEBRUIKSAANWIJZING	9
ES - INSTRUCCIONES DE USO	12
DE - GEBRAUCHSANWEISUNG	15
PT - INSTRUÇÕES DE USO	18
IT - ISTRUZIONI PER L'USO	21
EL - ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ	24
AR العربية - تعليمات الاستخدام	27
TR - KULLANIM TALIMATLARI	30

EN - INSTRUCTION FOR USE

Intended Use

Non-vascular Guidewires are intended to be placed in the human body through a natural body orifice via endoscopes; to be guided by controllable distal and proximal ends to a certain area of the human body and to be used as a guide for catheters or any other medical devices, which are intended to be used in the certain area of the body. Those areas do not include the blood vessels and/or the central circulatory system.

Product Description

The device can either be made of nitinol or stainless steel core and can be recovered of plastic components. It can have one or 2 flexible ends. The distal end is radiopaque for ROX vision.

Advices before first application

Please read carefully and follow all safety operating instructions and warnings before first application of the device.

A previous knowledge regarding handling and operation is required and essential.

Unpack the instrument carefully and examine the devices for any possibility of damage. In case of any damage or missing items contact your distributor immediately.

If the package present any damaged, the sterility of the device is not guarantee.

Handling and Operation

Prior to the extraction of the Non-vascular Guidewires from the dispenser, physiological saline solution has to be injected in the dispenser to wetting the complete hydrophilic surface of the Non-vascular Guidewires. After the injection of the saline solution, the Non-vascular Guidewires has to be extracted carefully out of the dispenser. If the Non-vascular Guidewires cannot be extracted easily out of the dispenser, additional saline solution has to be injected into the dispenser. After extraction, the Non-vascular Guidewires must not be introduced into the dispenser again.

Depending on the model, the Non-vascular Guidewires can be used from both ends. On the contrary, the distal end is always the flexible tip.

If not used correctly Non-vascular Guidewires can lead to perforation of tissue. The damages caused by perforation may be serious and can lead to death of the patient. Tissue perforation likely occurs due to kinking of the Non-vascular Guidewires. During the use of Non-vascular Guidewires it is very important to pay attention to the insertion technique, to avoid a kinking of the wire guide. An insertion guide has to be used.

 Non-vascular Guidewires contain magnetising materials such as high-grade steel and must therefore not be used under MR conditions, as heating and movement of the Non-vascular Guidewires can occur, which may lead to severe complications.

 According to clinical data, a passing through strictures of > 1cm length may lead to damage or breakage of the Non-vascular Guidewires and thus should be taken into thorough consideration during surgical planning. In rare cases, repeated leading back and forth of the Non-vascular Guidewires may lead to knot building at the distal tip of the wire. The user should have appropriate knowledge in such techniques in order to prevent serious consequences related to such events.

 If coated Non-vascular Guidewires are used special attention has to be paid to the combination product (e.g. catheter) to avoid damages to the integrity of the coating.

Procedure

- Insert the device in the endoscope operating channel with the help of a catheter or catheter made device.
- Advance the device to the desired location under fluoroscopic guidance.

 The stiffness of the device must be taken into consideration to avoid organ or tissue perforation. It is known that rigid devices should be avoided in some cases and flexible should then be preferred.

Contraindication

- Contraindications for guide wires > 0.014" (0.35 mm) can be Chronical Total Occlusion (CTOs).
- Guide wires containing magnetising materials such as high-grade steel and may therefore not be used under MR conditions, as heating and movement of the guide wire can occur, which may lead to catastrophic complications.
- Certain cases of entrapped or breakdown Non-vascular Guidewires are reported as well as possible techniques to retrieve the entrapped or broken parts. The user should have knowledge in such techniques in order to prevent serious consequences related to such events.
- If coated Non-vascular Guidewires are used special attention has to be paid to the combination product (e.g. catheter) to avoid damages to the integrity of the coating.
- Malposition of the catheterizing Non-vascular Guidewires may lead to serious complications.
- This device shall not be used in contact with blood vessels and/or the central circulatory system.

Storage

This device is delivered sterilized and in order to maintain this state until its first use it must be kept in its original sterilized package.

Do not put any objects on the instrument or its package!
Do not store the instruments near aggressive chemical products!

Do not expose the instruments to direct or indirect sunlight or other ultra-violet rays!

Keep in dry area

Complaints will not be taken under consideration if the instruments have been stored improperly.

Cleaning and Sterilization

This device is designed and warranted for single use!

Risk in case of reuse

This device is for single use only. G-Flex denied all responsibility in case of reuse. The reuse of single use device might present risk to the safety of the user or patient, due to possible uncontrolled contamination and / or unreliability of the device performance.

User

The users of G-FLEX instruments must be specialists in their fields. An appropriate and specific training for preparation, care and maintenance of the flexible instruments is required.

Repairing and complaints

In health protection of our employees only disinfected or sterilized instruments will be accepted for analysis or repairing. In any case disinfection or sterilization dates and validity must be labelled outside the package. If this requirement has not been fulfilled the instruments will be returned to the sender without analysis or repairing.

Legal foundation

The Law of the European Union applied.

Product support

In case of questions or difficulties concerning our instruments please contact your local distributor or G-FLEX directly during regular working hours.

Times available: Monday to Friday; 9am to 5pm (MET)

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