EN - SPHINCTEROTOME PAPILOTOME
FR - SPHINCTÉROTOME PAPILOTOME
NL - SFINCTEROTOOM PAPILOTOOM
ES - ESFINTERÓTOMO PAPILÓTOMO
DE - SPHINKTEROTOM/PAPILOTOM
PT - ESFINCTERÓTOMO/PAPILÓTOMO
IT - SFINTEROTOMA O PAPILOTOMA
EL - ΣΦΙΓΚΤΗΡΟΤΌΜΟΣ ΘΗΛΟΤΌΜΟΣ
AR - مبضع المصرّة/مبضع الحُلَيْمة
TR - SFINKTEROTOM PAPILOTOM

G-Flex Europe S.P.R.L. - 20, Rue de l'industrie 1400
Nivelles - Belgium | Phone: +32 67.88.36.65
Email: info@g-flex.com | Website: www.g-flex.com

Current Name: IFU-SphincterotomePapillotome

Creation Date: 01/10/2018

Version 7

This document is only controlled in electronic format by the QM.

EN - INSTRUCTION FOR USE	3
FR - MODE D'EMPLOI	6
NL - GEBRUIKSAANWIJZING	10
ES - INSTRUCCIONES DE USO	14
DE - GEBRAUCHSANWEISUNG	17
PT - INSTRUÇÕES DE USO	21
IT - ISTRUZIONI PER L'USO	25
EL - ΟΔΗΓΪΕΣ ΧΡΪΗΣΗΣ	28
AR - العربية - تعليمات الاستخدام	32
TR - KULLANIM TALIMATLARI	36

EN - INSTRUCTION FOR USE

Intended Use

Sphincterotome / Papillotome are intended to perform sphincterotomy (papillotomy) through the use of high frequency electric current.

Product Description

The device is made of a teflon flexible tube with a cutting wire at the distal end in different shapes and sizes. At the proximal end, there is a handle with an HF connector. Some models have a second and/or third lumen that can be used to inject contrast through a lateral port and/or others have a lumen that can be used with a guide wire. Please refer to the catalogue for more information

Advices before first application

Please read carefully and follow all safety operating instructions and warnings before first application of the device.

A previous knowledge regarding handling and operation is required and essential.

Unpack the instrument carefully and examine the devices for any possibility of damage. In case of any damage or missing items contact your distributor immediately.

If the package present any damaged, the sterility of the device is not guarantee.

Handling and Operation

The device must be passed through the endoscope operating channel with its cutting wire free from tension.

Please observe compatibility between the operating channel and device diameter (see catalogue).

The surgery room should contain non flammable substances.

Remove the stiffening wire before use.

The device when applied to a patient with a pacemaker implanted may cause malfunction or failure of the pacemaker, seriously affecting the patient. Before proceeding, always confirm with a cardiologist or the manufacturer of the pacemaker that it is safe.

After the device is connected with the equipment, the cutting wire cannot touch any other equipment. Spark discharge may cause clinician burns.

When using an electrocardiograph or other physiological monitoring equipment simultaneously with the device on a patient, any monitoring equipment should be placed as far away as possible from the electrodes used with the device. Needle monitoring electrodes should not be used, as they may cause patient burns.

-  Do not bend the tube (< 15mm). It may cause tube to be deformed permanently, thus, the device cannot be used.
-  When using the device in the vicinity of the heart, be sure to use it with the minimum necessary output. Spark discharge during operation may affect the heart.
-  During the preparation, do not quickly push and pull the thumb ring causing the distal end of the tube to be bent less than 90°. Do not do anything that might cause tube to kink. Rotating the handle while inserting the device into endoscope channel would cause tube to kink.
-  Sphincterotomes / Papillotomes are composed by a catheter and for that should always be straighten when used to avoid device damage.
-  **For the K18N, K18NDL and K18NTL models, the handle need to be operated forward. Operating it backward may damage the device.**

Procedure

- Connect the active cord to HF connector at the handle and to the HF-generator. (refer to the HF generator manufacturer instruction for proper connection)
- After connecting the active cord to the HF-generator the sphincterotomy (papillotomy) could be initiated.

-  Do not touch other tissue areas with the cutting wire when using HF electric current (high risk of hurt for the patient).
-  Do not inject contrast when using HF electric current.
-  For devices that accept guide wire, it is strongly recommended to pull the guide wire back into the tube while using HF current in order to avoid a possible induction of voltage on the guide wire

Procedure for Knife Snake®

Whenever a precut papillotome is needed after a regular (lateral cut) wire papillotome is already in place, instead of exchanging the current papillotome for a needle-knife papillotome, an alternative solution can be the Knife-Snake® precut papillotome. For its use, please follow the instructions below:

- Remove completely the guidewire from the guidewire lumen of your double/triple lumen papillotome.
- With the papillotome and duodenoscope elevator untensioned, through the guide wire lumen, insert the Knife Snake® until its end can be seen through endoscopic view.
- Attach the yellow connector of the device to the Luer-Lock connector from the guide wire lumen entry of the papillotome and ensure it is well fixed.
- Now bend the papillotome to the desired precut position.
- Adjust the length of the Knife Snake® to the desired length of cutting wire and then fix the length at the yellow connector by pushing it distally.
- Connect the active cord to HF connector at the handle and to the HF-generator. (refer to the HF generator manufacturer instruction for proper connection)

- After connecting the active cord to the HF-generator the sphincterotomy (papillotomy) could be initiated.
- Once the desired cut opening has been reached, disconnect the HF active cord from the device and remove it from the papillotome replacing it again by a guide wire.

Ensure that your papillotome guide wire lumen accepts a 0,035" guide wire

Do not bend the papillotome or endoscope while inserting the device as this may block the device.

Storage

This device is delivered sterilized and in order to maintain this state until its first use it must be kept in its original sterilized package.

Do not put any objects on the instrument or its package!
Do not store the instruments near aggressive chemical products!

Do not expose the instruments to direct or indirect sunlight or other ultra-violet rays!

Keep in dry area

Complaints will not be taken under consideration if the instruments have been stored improperly.

Cleaning and Sterilization

This device is designed and warranted for single use!

Risk in case of reuse

This device is for single use only. G-Flex denied all responsibility in case of reuse. The reuse of single use device might present risk to the safety of the user or patient, due to possible uncontrolled contamination and / or unreliability of the device performance.

Contraindications

- Patients with severe cardiopulmonary or renal dysfunction.
- Patients with acute pancreatitis or acute exacerbation of chronic pancreatitis.
- Patients with severe biliary tract infection.
- Patients who are allergic to iodine contrast agents.
- Patients who have upper gastrointestinal stenosis or obstruction.
- Contraindications for this device are those specific to endoscopic retrograde cholangiopancreatography(ERCP) and endoscopic sphincterotomy(ES).

User

The users of G-FLEX instruments must be specialists in their fields. An appropriate and specific training for preparation, care and maintenance of the flexible instruments is required.

Repairing and complaints

In health protection of our employees only disinfected or sterilized instruments will be accepted for analysis or repairing. In any case disinfection or sterilization dates and validity must be labelled outside the package. If this requirement has not been fulfilled the instruments will be returned to the sender without analysis or repairing.

Legal foundation

The Law of the European Union applied.

Product support

In case of questions or difficulties concerning our instruments please contact your local distributor or G-FLEX directly during regular working hours.

Times available: Monday to Friday; 9am to 5pm (MET)

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