EN - STENTS APPLICATION SYSTEM
FR - SYSTÈME D'APPLICATION D'ENDOPROTHÈSE
NL - STENTPLAATSINGS SYSTEEM
ES - INTRODUCTORES DE STENT
DE - STENTS-APPLIKATIONSSYSTEM
PT - SISTEMA DE APLICAÇÃO DE PRÓTESE
IT - SISTEMA DI APPLICAZIONE DE STENT
EL - ΣΥΣΤΗΜΑ ΕΦΑΡΜΟΓΉΣ ΤΩΝ
ΕΝΔΟΠΡΟΘΈΣΕΩΝ
AR - نظام تطبيق الدعامة
TR - STENTI UYGULAMA SISTEMI

G-Flex Europe S.P.R.L. - 20, Rue de l'industrie 1400
Nivelles - Belgium | **Phone:** +32 67.88.36.65
Email: info@g-flex.com | **Website:** www.g-flex.com

EN - INSTRUCTION FOR USE	3
FR - MODE D'EMPLOI	6
NL - GEBRUIKSAANWIJZING	10
ES - INSTRUCCIONES DE USO	14
DE - GEBRAUCHSANWEISUNG	18
PT - INSTRUÇÕES DE USO	22
IT - ISTRUZIONI PER L'USO	26
EL - ΟΔΗΓΊΕΣ ΧΡΉΣΗΣ	30
AR العربية - تعليمات الاستخدام	35
TR - KULLANIM TALIMATLARI	38

EN - INSTRUCTION FOR USE

Intended Use

Stents Application System are intended to place stent via endoscopes.

Product Description

The device is made of an inner guiding catheter with stainless steel marker at the distal end and a Y connector at the proximal end. They also have an outer pusher with stainless steel marker at the distal end and a Luer-lock adapter at the proximal end.

Advices before first application

Please read carefully and follow all safety operating instructions and warnings before first application of the device.

A previous knowledge regarding handling and operation is required and essential. Unpack the instrument carefully and examine the devices for any possibility of damage. In case of any damage or missing items contact your distributor immediately.

Handling and Operation

The Endoscopic Drainage Stents are placed in conjunction with an endoscope, a guide wire, a pusher and, when necessary, a guiding catheter

To facilitate the application, an endoscopic sphincterotomy prior to guide wire placement is recommended. (Please observe HF-generator and sphincterotome manufacturer's instructions).

It is recommended NOT to remove the stent by inflating a balloon inside the stent.

The patient and the stent should be regularly followed-up to remove the stent immediately in case of clog, migration or fallout from the papilla.

Procedure

For the Application System, please proceed as follows:

- Use the included positioning sleeve for easier introducing by sliding it to the proximal end of the stent. So the rear flap will fit back into it.
- Ensure that the Luer-Lock connection between pusher and guiding catheter is firmly closed. Loosen the sealing cap at the Y-adapter.
- Place the biliary stent on to the device by sliding it over the distal end of guiding catheter up to the pusher.
- Now introduce the loaded device over the proximal end of the guide wire into the operating channel. Positioning sleeve does not fit to the channel and will stay outside.

- Advance the device under constant fluoroscopic control carefully until the radiopaque metal tip of the guiding catheter passes the obstructed area.
- Closing the sealing cap at the Y-adapter prevents movement of the guide wire inside of the guiding catheter. Contrast injection is possible through the Luer-Lock-Port at the Y-adapter.
- Open up the Luer-Lock connection between guiding catheter and pusher and advance the biliary stent to its final position by moving forward the pusher carefully.

Keep guiding catheters position while placing the stent

- Pull back the guiding catheter and guide wire into the operating channel. Prevent movement of the stent by keeping the pusher in its last position. After that the stent will return to its original shape.
- Finally, remove pusher, guiding catheter and guide wire out of the duodenoscope.

Injection of contrast medium into the biliary duct is possible via the guiding catheter by using the included Luer-Lock adapter.

Contraindication

The contraindications include without being limited to poor general condition such as respiratory failure, risk of aspiration, general blood coagulation and wound healing disorders, asymptomatic stenosis or jaundice, endoscopically impassible stenosis, acute, non-biliary pancreatitis, absent patient consent.

Storage

This device is delivered non sterile and must be sterilised before first use according with the instruction described below.

Do not put any objects on the instrument or its package!
Do not store the instruments near aggressive chemical products!

Do not expose the instruments to direct or indirect sunlight or other ultra-violet rays!

Keep in dry area

Complaints will not be taken under consideration if the instruments have been stored improperly.

Cleaning and Sterilization

G-Flex recommend to follow one of the following sequences when sterilizing the device:

- Manual Cleaning, Ultrasonic Cleaning, Automatic/Thermal Cleaning and Autoclave; or
- Manual Cleaning, Ultrasonic Cleaning, Automatic/Thermal Cleaning and Ethylene Oxide Gas Sterilization

Either sterilize by Autoclave or by ETO. There is no need to do both.

Manual Cleaning

The cleaning of the device should be done immediately after each use:

- Immerse the device in a suitable liquid detergent or disinfectant (Always observe the minimum diameter of 20 cm when winding the instrument and applicability of the cleaning product!).
- Wipe the instrument with a soft towel.

Ultrasonic Cleaning

The ultrasonic cleaning facilitates the dislodging of residual material and must be done before sterilization. The ultrasonic cleaning must last for at least 30 minutes.

- Clean the device right after its use.
- Put the instrument in an ultrasonic cleaning tank always observing the winding minimum diameter 20 cm (Please observe the ultrasonic equipment manufacturer's instructions!).
- Use only tap water and the detergents which dissolve albumen, always according to ultrasonic equipment manufacturer's instructions.
- Rinse the instrument with tap water and wipe dry by using gauze pads.

Automatic / Thermal Cleaning

Alternatively or additionally to the manual cleaning, this instrument may also be cleaned in automatic / thermal cleaning machines by using approved alkaline detergents. Please refer to relevant instructions of the cleaning machine supplier.

Sterilization - Autoclave

Before autoclaving sterilization the instrument should be cleaned carefully as described in Manual Cleaning, Ultrasonic Cleaning and Automatic/Thermal Cleaning.

This symbol means that the product can be sterilized by autoclave according to the below parameters

Parameters:

Temperature: 134 °C / 273°F

Pressure: 3 bar

Exposure: ≥ 18 min

Also in here the minimum 20 cm diameter for winding the instrument must be noted.

The original packages of G-FLEX are not autoclavable!

Sterilization - Ethylene Oxide Gas Sterilization

Please follow manufacturer's warnings regarding to the biological indicators

Parameters:

Temperature:	57 °C
Pressure:	1,7 bar
Exposure:	max. 4 h
Relative Humidity:	50 %
Gas Concentration:	12 %
Aeration Time:	7 days in ambient temperature or 12 hours by 50-57 °C

User

The users of G-FLEX instruments must be specialists in their fields. An appropriate and specific training for preparation, care and maintenance of the flexible instruments is required.

Repairing and complaints

In health protection of our employees only disinfected or sterilized instruments will be accepted for analysis or repairing. In any case disinfection or sterilization dates and validity must be labelled outside the package. If this requirement has not been fulfilled the instruments will be returned to the sender without analysis or repairing.

Product disposal

After use, these products may be a potential biohazard. Handle and dispose of in accordance with accepted medical practice and applicable local, state and federal laws and regulations.

Legal foundation

The Law of the European Union applied.

Product support

In case of questions or difficulties concerning our instruments please contact your local distributor or G-FLEX directly during regular working hours.

Times available: Monday to Friday; 9am to 5pm (MET)

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