**EN - STONE EXTRACTION BALLOON
CATHETERS**

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EN - INSTRUCTION FOR USE

Intended Use

Stone Extraction Balloon Catheters are intended to remove stones from the biliary duct. When necessary, contrast injection is possible through one of the Luer-Lock connectors.

Product Description

The device is made of a teflon tube with a latex balloon at the distal end and a Y connector at the proximal end for the balloon inflation and injection of contrast or the use of a guide wire.

Advices before first application

Please read carefully and follow all safety operating instructions and warnings before first application of the device.

A previous knowledge regarding handling and operation is required and essential. Unpack the instrument carefully and examine the devices for any possibility of damage. In case of any damage or missing items contact your distributor immediately.

If the package present any damaged, the sterility of the device is not guarantee.

Handling and Operation

The device must be inserted deflated into the endoscope operating channel. Please observe compatibility of catheter diameter and the endoscope operating channel diameter (see catalogue).

If possible, the use of a guide wire is recommendable.

Procedure

After locating the stone to be removed, proceed as follows:

- Insert the deflated balloon into the biliary duct until the distal tip is located above the stone to be removed.
- Pull the plunger of the included syringe and connect it to the transparent stopcock at the proximal end of the catheter .
- Open the stopcock.
- Push the plunger of the syringe carefully to inflate the balloon. Air capacity specified for each balloon must be respected to avoid its destruction.
- Close the stopcock to prevent the air to leak out of the balloon.
- Gently withdraw the inflated balloon catheter to pull the stones out of the biliary duct and drop them into the duodenum. They can be expelled naturally out of the body through the intestines.
- Open the stopcock and deflate the balloon.
- Remove the entire device out of the endoscope carefully.

Sharp stones can destroy the balloon and should be removed with an extractor basket or a lithotripsy basket.

Contraindications

Use of this natural latex rubber balloon is contraindicated in patients with a known hypersensitivity to latex.

Complications

Potential complication: pancreatitis, cholangitis, aspiration, perforation, hemorrhage, infection, sepsis, allergic reaction to contrast or medication, hypertension, respiratory depression or arrest, cardiac arrhythmia or arrest.

Additional complications that can occur during endoscopic balloon extraction include: stone impaction, localized inflammation, pressure necrosis.

Storage

This device is delivered sterilized and in order to maintain this state until its first use it must be kept in its original sterilized package.

Do not put any objects on the instrument or its package!
Do not store the instruments near aggressive chemical products!

Do not expose the instruments to direct or indirect sunlight or other ultra-violet rays!

Keep in dry area

Complaints will not be taken under consideration if the instruments have been stored improperly.

Cleaning and Sterilization

This device is designed and warranted for single use!

Risk in case of reuse

This device is for single use only. G-Flex denied all responsibility in case of reuse. The reuse of single use device might present risk to the safety of the user or patient, due to possible uncontrolled contamination and / or unreliability of the device performance.

User

The users of G-FLEX instruments must be specialists in their fields. An appropriate and specific training for preparation, care and maintenance of the flexible instruments is required.

Repairing and complains

In health protection of our employees only disinfected or sterilized instruments will be accepted for analysis or repairing. In any case disinfection or sterilization dates and validity must be labelled outside the package. If this requirement has not been fulfilled the instruments will be returned to the sender without analysis or repairing.

Legal foundation

The Law of the European Union applied.

Product support

In case of questions or difficulties concerning our instruments please contact your local distributor or G-FLEX directly during regular working hours.

Times available: Monday to Friday; 9am to 5pm (MET)

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G-Flex USA, Inc., 74 Southwood Pkwy Suite 600