



1. Intended Use

Stone Extraction Baskets are intended to catch and remove stones from ureter.

2. Product Description

The device is made of a basket at the distal end (nitinol) and a handle at the proximal end. The handle can be removable depending on the model.

3. Advices before first application

Please read carefully and follow all safety operating instructions and warnings before first application of the device.

A previous knowledge regarding handling and operation is required and essential.

Unpack the instrument carefully and examine the devices for any possibility of damage. In case of any damage or missing items contact your distributor immediately.



If the package present any damaged, the sterility of the device is not guarantee.

4. Handling and Operation

The device must be passed through the endoscope channel with its basket withdrawn into the tube. Please observe compatibility between the device and the endoscope operating channel diameters.

5. Procedure

- Once the stone is located by usual means do as following:
- Place the Stone Extraction Basket's distal tip approx. 20 mm in front of the stone to be removed.
- Push the basket completely out of the tube and bring the stone in the basket centre.
- Pull the handle back until the stone is fixed in the basket and pull the device out of the endoscope until the basket touches the distal tip of the endoscope.



Excessive traction should be avoided to prevent damage the endoscope.

- By carefully taking out the endoscope and the device, the stone will be removed from the ureter. Keep the handle of the device in position for fixing the stone.

6. Contraindication



Stone Extraction Basket (Nitinol) are not recommended for Lithotripsy.

7. Storage

This device is delivered sterilized and in order to maintain this state until its first use it must be kept in its original sterilized package.



Do not put any objects on the instrument or its package!

Do not store the instruments near aggressive chemical products!



Do not expose the instruments to direct or indirect sunlight or other ultra-violet rays!



Keep in dry area

Complaints will not be taken under consideration if the instruments have been stored improperly.





8. Cleaning and Sterilization



This device is designed and warranted for single use!

9. Risk in case of reuse

This device is for single use only. G-Flex denied all responsibility in case of reuse. The reuse of single use device might present risk to the safety of the user or patient, due to possible uncontrolled contamination and / or unreliability of the device performance.

10. User

The users of G-FLEX instruments must be specialists in their fields. An appropriate and specific training for preparation, care and maintenance of the flexible instruments is required.

11. Repairing and complains

In health protection of our employees only disinfected or sterilized instruments will be accepted for analysis or repairing. In any case disinfection or sterilization dates and validity must be labelled outside the package. If this requirement has not been fulfilled the instruments will be returned to the sender without analysis or repairing.

12. Legal foundation

The Law of the European Union applied.

13. Product support

In case of questions or difficulties concerning our instruments please contact your local distributor or G-FLEX directly during regular working hours.

Times available: Monday to Friday; 9am to 5pm (MET)

14. Symbols



Date of manufacture



Do not use if package is damaged



Caution



Use-by date



Keep away from sunlight



Do not reuse



Batch or lot code



Keep dry



Consult the instruction for use



Sterilized using ethylene oxide



Quantity of units per box



Temperature limit

