EN - CYSTO-GASTRO SET

FR - SET CYSTO-GASTRO

NL - CYSTE-MAAG-SET

ES - SET DE GASTRO-QUISTE

DE - CYSTO-GASTRO SET

PT - CONJUNTO CISTO-GÁSTRICO

IT - SET CISTO-GASTRO

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
EN - INSTRUCTION FOR USE

Intended Use

Cysto-Gastro-Sets are intended to be used to puncture pancreatic cysts or pseudo-cysts endoscopically (via the transgastric or transduodenal wall) and to dilate the opening to the cyst (6, 8.5 or 10Fr models) as an alternative to surgical or percutaneous treatments. It is also indicated to dilate the access in the procedure of Endoscopic ultrasonography-guided biliary drainage (EUS-BD).

Product Description

The device is made of an outer flexible teflon tube equipped with a metal tip at the distal end and an HF connector in the proximal end (all models) images on page 3:

1. HF dilator
2. Guiding catheter

Handle

3. Luer-lock
4. Guide wire entry port
5. Dilator HF Connector

Advices before first application

Please read carefully and follow all safety operating instructions and warnings before first application of the device.

A previous knowledge regarding handling and operation is required and essential.

Unpack the instrument carefully and examine the devices for any possibility of damage. In case of any damage or missing items contact your distributor immediately.

If the package present any damaged, the sterility of the device is not guarantee.

Handling and Operation

Both parts of the device are supplied with HF connectors that need to be fixed to an appropriate active cords.

It is recommended to use an echo-endoscope.

Cysto-Gastro-Sets require an endoscope's working channel of at least 3.7 mm Ø for 10Fr Set, 3.2 mm Ø for the 8.5 Fr Set or 2.8 mm for the 6 Fr.

Use cut frequency only since the use or excess of coagulation frequency can damage the product. Recommended frequency is between 80 and 120watts.

Procedure

After the size and position of the cyst or pseudo-cyst are determined by endosonography, the device is inserted into the endoscope's operating channel over a guide wire (.025" or .035").

- For the 6 Fr model you can puncture the stomach or duodenal wall and the cyst directly with the catheter of the device.
- For the 8.5 and 10 Fr model the stomach or the duodenal wall and the cyst have to be previously puncture with an EUS Needle or a diathermic needle.
- If necessary, fluids can be aspirated from the cavity through the lateral Luer-Lock adapter at the proximal end of the device.
- Now connect the active cord to the lateral connector and to the HF generator.
- Switch on electricity and dilate the access previously created and cyst's opening by pushing forward the device's metal tip carefully.
- Remove the active cord and leave the device at its position.
- Finally, remove the device out of the working channel, while ensuring that the guide wire stays in position.
- Use the guide wire to place a Double Pigtail or Metallic stent in order to drain the cyst, after the cystostomy.

Possible Complications

Possible complications associated with the gastrointestinal endoscopical puncture of cysts or pseudo-cysts include pain, abscess, hemorrhage/bleeding, pancreatitis, perforation small risk of infection, abdominal pain, coagulation disorders and contraindications related to the use of electrical scalpel.

Storage

This device is delivered sterilized and in order to maintain this state until its first use it must be kept in its original sterilized package.

Do not put any objects on the instrument or its package!

Do not store the instruments near aggressive chemical products!

Do not expose the instruments to direct or indirect sunlight or other ultra-violet rays!

Keep in dry area

Complaints will not be taken under consideration if the instruments have been stored improperly.

Cleaning and Sterilization

The disposable version of this device is designed and warranted for single use!

Risk in case of reuse

This devices is for single use only. G-Flex denied all responsibility in case of reuse. The reuse of single use device might present risk to the safety of the user or patient, due to possible uncontrolled contamination and / or unreliability of the device performance.

User

The users of G-FLEX instruments must be specialists in their fields. An appropriate and specific training for preparation, care and maintenance of the flexible instruments is required.

Repairing and complaints

In health protection of our employees only disinfected or sterilized instruments will be accepted for analysis or repairing. In any case disinfection or sterilization dates and validity must be labelled outside the package. If this requirement has not been fulfilled the instruments will be returned to the sender without analysis or repairing.

Product disposal

After use, these products may be a potential biohazard. Handle and dispose of in accordance with accepted medical practice and applicable local, state and federal laws and regulations.

Legal foundation

The Law of the European Union applied.

Product support

In case of questions or difficulties concerning our instruments please contact your local distributor or G-FLEX directly during regular working hours.

Times available: Monday to Friday; 9am to 5pm (MET)

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