
EN - DISPOSABLE EXTRACTION BASKETS
FR - PANIER D'EXTRACTION JETABLE
NL - WEGWERP-EXTRACTIEMAND
ES - CESTA DE EXTRACCIÓN DESECHABLE
**DE - EXTRAKTIONSKÖRBE FÜR DEN
EINMALGEBRAUCH**
PT - CESTOS DE EXTRAÇÃO DESCARTÁVEIS
IT - CESTELLI DI ESTRAZIONEMONOUSO
EL - ΚΑΛΪΘΗ ΕΞΑΓΩΓΉΣ ΜΪΑΣ ΧΡΉΣΗΣ
AR - سلات الاستخراج أحادية الاستعمال
TR - TEK KULLANIMLIK ÇIKARMA SEPETLERİ

G-Flex Europe S.P.R.L. - 20, Rue de l'industrie 1400
Nivelles - Belgium | Phone: +32 67.88.36.65
Email: info@g-flex.com | Website: www.g-flex.com

EN - INSTRUCTION FOR USE	5
FR - MODE D'EMPLOI	8
NL - GEBRUIKSAANWIJZING	12
ES - INSTRUCCIONES DE USO	15
DE - GEBRAUCHSANWEISUNG	19
PT - INSTRUÇÕES DE USO	22
IT - ISTRUZIONI PER L'USO	26
EL - ΟΔΗΓΊΕΣ ΧΡΉΣΗΣ	29
AR - العربية - تعليمات الاستخدام	33
TR - KULLANIM TALIMATLARI	37

EN - INSTRUCTION FOR USE

Intended Use

Disposable Extraction Baskets are intended to catch and remove stones or to crush stones in the biliary and pancreatic ducts of the human body.

Product Description

The **extraction basket** is made of a basket at the distal end and a handle at the proximal end.

1. Basket
2. Catheter

The **lithotripsy handle** is made of a winding cylinder with a locking system fixed to a wheel, an aluminum body with an opening for the spiral and a handle (see page 4). The lithotripsy handle is delivered non-sterile and is made of the following parts:

8. Spiral entrance
9. Winding cylinder with orifice for tracking wire
10. Locking system
11. Wheel
12. Handle

Advices before first application

Please read carefully and follow all safety operating instructions and warnings before first application of the device.

A previous knowledge regarding handling and operation is required and essential. Unpack the instrument carefully and examine the devices for any possibility of damage. In case of any damage or missing items contact your distributor immediately.

If the package presents any damage, the sterility of the device is not guaranteed.

Handling and Operation

The device has to be withdrawn completely into the Teflon tube to be introduced into the endoscope operating channel. Please observe compatibility between the device and the endoscope operating channel diameters.

Procedure

Extraction of Stone

- Once the stone is located by usual means proceed as following:
- Place the Extraction Basket's distal tip above the stone to be removed.
- Push the basket completely out of the catheter and pull the basket down to the stone.
- If necessary, rotate the basket or move the catheter up and down in short and quick moves in order to facilitate to place the basket wires around the stone.

- Pull the handle back until the stone is fixed in the basket and pull the device toward the endoscope until the stone is extracted.

Excessive traction should be avoided to prevent damage the endoscope.

Emergency lithotripsy

- In case a stone gets stuck inside of an Extraction Basket proceed as following:
- With the help of a cutting plier, cut the catheter and tracking wire of the basket as close as possible to the handle of the device, just before the shrinking tube.
- Dispose the handle and pull gently the catheter over the cut wire 1-2cm to ensure the wire is protected and can not damage the endoscope.
- While keeping both tracking wire and catheter inside of the patient, gently remove the endoscope. Ensure not to move the catheter up or down while removing the endoscope.
- Once the endoscope is completely removed, remove the plastic catheter and replace it by the short and thick metallic spiral (Ref. GF641).
- Once the spiral reaches the stone, push the tracking wire through the lithotripsy handle spiral entrance until a slight tension can be felt.
- Lock the lithotripsy handle at the locking system.
- Gently begin to rotate the wheel until the stone gets crushed. When applying rotation on the wheel, allow some delays so that the mechanical force is perfectly transferred to the distal end of the spiral before repeating a new rotation (at least 40 seconds)
- Remove the device from the patient.

It is essential to keep the spiral as straight as possible to avoid friction between the tracking wire and the spiral. Any friction could lead to break the tracking wire and cause the procedure to fail.

Do not rotate the wheel quickly as this may lead to break the tracking wire and cause the procedure to fail.

Handle Assembling (See page 4)

- Keep the locking system in the open position
- Push the traction cable tip through the spiral entrance to the winding cylinder orifice.
- Push the handle over the traction cable until it couples to the spiral tip and insert the spiral tip in the spiral entrance until resistance is felt (excessive force should be avoided).
- Push the spiral straighten to the basket.
- To release the tracking wire, place the locking system in the close position.

Crushing the Stone

- Turn the manual wheel at the handle, until the stone is crushed. When applying rotation on the wheel, allow some delays so that the mechanical force is perfectly transferred to the distal end of the spiral before repeating a new rotation (at least 40 seconds)

A lithotripsy has no guarantee of total success! If the stone to be crushed is very hard the basket might be torn because of the force applied on it. In this case the stone along with the torn basket should be surgically removed.

The lithotripsy should be interrupted immediately if the stone crushing could not be achieved.

Possible Complications

Possible complications associated with perforation, bleeding, pancreatitis, cholangitis and cholecystitis.

Storage

This device is delivered sterilized and in order to maintain this state until its first use it must be kept in its original sterilized package.

Do not put any objects on the instrument or its package!

Do not store the instruments near aggressive chemical products!

Do not expose the instruments to direct or indirect sunlight or other ultra-violet rays!

Keep in dry area

Complaints will not be taken under consideration if the instruments have been stored improperly.

Cleaning and Sterilization

The Disposable Extraction Basket is designed and warranted for single use!

Risk in case of reuse

The Disposable extraction basket is for single use only. G-Flex denied all responsibility in case of reuse. The reuse of single use device might present risk to the safety of the user or patient, due to possible uncontrolled contamination and / or unreliability of the device performance.

User

The users of G-FLEX instruments must be specialists in their fields. An appropriate and specific training for preparation, care and maintenance of the flexible instruments is required.

Repairing and complaints

In health protection of our employees only disinfected or sterilized instruments will be accepted for analysis or repairing. In any case disinfection or sterilization dates

and validity must be labelled outside the package. If this requirement has not been fulfilled the instruments will be returned to the sender without analysis or repairing.

Product disposal

After use, these products may be a potential biohazard. Handle and dispose of in accordance with accepted medical practice and applicable local, state and federal laws and regulations.

Legal foundation

The Law of the European Union applied.

Product support

In case of questions or difficulties concerning our instruments please contact your local distributor or G-FLEX directly during regular working hours.

Times available: Monday to Friday; 9am to 5pm (MET)

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