EN - EXTRACTION BASKETS & LITHOTRIPSY SYSTEM

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


Current Name: IFU-Basket & Lithotripsy System

Creation Date: 11/12/2019


Version: 1

This document is only controlled in electronic format by the QM.

EN - INSTRUCTION FOR USE	6
FR - MODE D'EMPLOI.....	11
NL - GEBRUIKSAANWIJZING	16
ES - INSTRUCCIONES DE USO	21
DE - GEBRAUCHSANWEISUNG	27
PT - INSTRUÇÕES DE USO	32
IT - ISTRUZIONI PER L'USO	38
EL - ΟΔΗΓΊΕΣ ΧΡΉΣΗΣ	43
AR العربية - تعليمات الاستخدام.....	49
TR - KULLANIM TALIMATLARI	54

EN - INSTRUCTION FOR USE

Intended Use

Extraction Basket & Lithotripsy System are intended to catch and remove stones or to crush stones in the biliary and pancreatic ducts of the human body.

Product Description

The **extraction basket** is made of a basket at the distal end and a handle at the proximal end. Some models have a removable handle (see page 3)

1. Basket
2. Catheter
3. Luer-lock connection
4. Locking system of the removable handle
5. Removable handle

The **spiral** is made of:

6. Rigid shaft
7. Spiral

The **lithotripsy handle** is made of a winding cylinder with a locking system fixed to a wheel, an aluminum body with an opening for the spiral and a handle (see page 4). The lithotripsy handle is delivered non-sterile and is made of the following parts:

8. Spiral entrance
9. Winding cylinder with orifice for tracking wire
10. Locking system
11. Wheel
12. Handle

Advices before first application

Please read carefully and follow all safety operating instructions and warnings before first application of the device.

A previous knowledge regarding handling and operation is required and essential.

Unpack the instrument carefully and examine the devices for any possibility of damage. In case of any damage or missing items contact your distributor immediately.

Handling and Operation

The device has to be withdrawn completely into the Teflon tube to be introduced into the endoscope operating channel. Please observe compatibility between the device and the endoscope operating channel diameters.

Procedure

Extraction of Stone

- Once the stone is located by usual means proceed as following:
- Place the Extraction Basket's distal tip above the stone to be removed.
- Push the basket completely out of the catheter and pull the basket down to the stone.

- If necessary, rotate the basket or move the catheter up and down in short and quick moves in order to facilitate to place the basket wires around the stone.
- Pull the handle back until the stone is fixed in the basket and pull the device toward the endoscope until the stone is extracted.

Excessive traction should be avoided to prevent damage the endoscope.

Emergency lithotripsy

- In case a stone gets stuck inside of an Extraction Basket proceed as following:
- With the help of a cutting plier, cut the catheter and tracking wire of the basket as close as possible to the handle of the device, just before the shrinking tube.
- Dispose the handle and pull gently the catheter over the cut wire 1-2cm to ensure the wire is protected and can not damage the endoscope.
- While keeping both tracking wire and catheter inside of the patient, gently remove the endoscope. Ensure not to move the catheter up or down while removing the endoscope.
- Once the endoscope is completely removed, remove the plastic catheter and replace it by the short and thick metallic spiral (Ref. GF641).
- Once the spiral reaches the stone, push the tracking wire through the lithotripsy handle spiral entrance until a slight tension can be felt.
- Lock the lithotripsy handle at the locking system.
- Gently begin to rotate the wheel until the stone gets crushed. When applying rotation on the wheel, allow some delays so that the mechanical force is perfectly transferred to the distal end of the spiral before repeating a new rotation (at least 40 seconds)
- Remove the device from the patient.

It is essential to keep the spiral as straight as possible to avoid friction between the tracking wire and the spiral. Any friction could lead to break the tracking wire and cause the procedure to fail.

Do not rotate the wheel quickly as this may lead to break the tracking wire and cause the procedure to fail.

Lithotripsy through the working channel (for baskets GF264GW, GF265GW, GF266GW, GF267GW only)

How to use the lithotripsy basket with a guide wire

This model of basket can be introduced into the bile duct directly or over a previously positioned guide wire. To do so, please follow the instructions below:

- Loosen the sealing cap at the Y-adapter.
- Remove the basket pulling it from its handle (without unscrewing it) from the distal outer teflon tube.
- Insert the teflon tube over the guide wire through the working channel of the endoscope until the metallic tip of the tube reaches a position behind the stone to be crushed.

- With the teflon tube in position, remove the guide wire from the proximal end of the teflon tube.
- With two fingers, squeeze the basket and introduce it through the proximal end of the teflon tube (see page 5) until the rigid part of the wire along with the handle reach the teflon tube.
- Slowly expose the basket and catch the aimed stone to be crushed.

Changing the Teflon tube for the metal spiral

For crushing the stone it is necessary to change the Teflon tube for the fitted metal spiral. These action makes a stable counter-pressure during the crushing procedure possible.

- Loosen the sealing cap at the Y-adapter.
- Remove the Teflon tube completely and carefully out of the endoscope, sliding it over the traction cable.
- Insert the metal spiral's distal end by sliding it over the traction cable into the operating channel until it reaches the basket with the stone inside.

Handle Assembling (See page 4)

- Keep the locking system in the open position
- Push the traction cable tip through the spiral entrance to the winding cylinder orifice.
- Push the handle over the traction cable until it couples to the spiral tip and insert the spiral tip in the spiral entrance until resistance is felt (excessive force should be avoided).
- Push the spiral straighten to the basket.
- To release the tracking wire, place the locking system in the close position.

Crushing the Stone

- Turn the manual wheel at the handle, until the stone is crushed. When applying rotation on the wheel, allow some delays so that the mechanical force is perfectly transferred to the distal end of the spiral before repeating a new rotation (at least 40 seconds)

A lithotripsy has no guarantee of total success! If the stone to be crushed is very hard the basket might be torn because of the force applied on it. In this case the stone along with the torn basket should be surgically removed.

The lithotripsy should be interrupted immediately if the stone crushing could not be achieved.

Possible Complications

Possible complications associated with perforation, bleeding, pancreatitis, cholangitis and cholecystitis.

Storage

This device is delivered non sterile and must be sterilised before first use according with the instruction described below.

Do not put any objects on the instrument or its package!
Do not store the instruments near aggressive chemical products!

Do not expose the instruments to direct or indirect sunlight or other ultra-violet rays!

Keep in dry area

Complaints will not be taken under consideration if the instruments have been stored improperly.

Cleaning and Sterilization

G-Flex recommends to follow one of the following sequences when sterilizing the device:

- Manual Cleaning, Ultrasonic Cleaning, Automatic/Thermal Cleaning and Autoclave; or
- Manual Cleaning, Ultrasonic Cleaning, Automatic/Thermal Cleaning and Ethylene Oxide Gas Sterilization

Either sterilize by Autoclave or by ETO. There is no need to do both.

Manual Cleaning

The cleaning of the device should be done immediately after each use:

- Immerse the device in a suitable liquid detergent or disinfectant (Always observe the minimum diameter of 20 cm when winding the instrument and applicability of the cleaning product!).
- Wipe the instrument with a soft towel.

Ultrasonic Cleaning

The ultrasonic cleaning facilitates the dislodging of residual material and must be done before sterilization. The ultrasonic cleaning must last for at least 30 minutes.

- Clean the device right after its use.
- Put the instrument in an ultrasonic cleaning tank always observing the winding minimum diameter 20 cm (Please observe the ultrasonic equipment manufacturer's instructions!).
- Use only tap water and the detergents which dissolve albumen, always according to ultrasonic equipment manufacturer's instructions.
- Rinse the instrument with tap water and wipe dry by using gauze pads.

Automatic / Thermal Cleaning

Alternatively or additionally to the manual cleaning, this instrument may also be cleaned in automatic / thermal cleaning machines by using approved alkaline detergents. Please refer to relevant instructions of the cleaning machine supplier.

Sterilization - Autoclave

Before autoclaving sterilization the instrument should be cleaned carefully as described in Manual Cleaning, Ultrasonic Cleaning and Automatic/Thermal Cleaning.

This symbol means that the product can be sterilized by autoclave according to the below parameters

Parameters:

Temperature: 134 °C / 273°F

Pressure: 3 bar

Exposure: ≥ 18 min

Also in here the minimum 20 cm diameter for winding the instrument must be noted.

The original packages of G-FLEX are not autoclavable!

Sterilization - Ethylene Oxide Gas Sterilization

Please follow manufacturer's warnings regarding to the biological indicators

Parameters:

Temperature: 57 °C

Pressure: 1,7 bar

Exposure: max. 4 h

Relative Humidity: 50 %

Gas Concentration: 12 %

Aeration Time: 7 days in ambient temperature or 12 hours by 50-57 °C

User

The users of G-FLEX instruments must be specialists in their fields. An appropriate and specific training for preparation, care and maintenance of the flexible instruments is required.

Repairing and complaints

In health protection of our employees only disinfected or sterilized instruments will be accepted for analysis or repairing. In any case disinfection or sterilization dates and validity must be labelled outside the package. If this requirement has not been fulfilled the instruments will be returned to the sender without analysis or repairing.

Product disposal

After use, these products may be a potential biohazard. Handle and dispose of in accordance with accepted medical practice and applicable local, state and federal laws and regulations.

Legal foundation

The Law of the European Union applied.

Product support

In case of questions or difficulties concerning our instruments please contact your local distributor or G-FLEX directly during regular working hours.

Times available: Monday to Friday; 9am to 5pm (MET)

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