**EN - SCLEROTHERAPY & ENDOSCOPIC
NEEDLES**

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Current Name: IFU-Sclerotherapy & Endoscopic Needles **Creation Date:** 03/11/2020 **Version** 10_

This document is only controlled in electronic format by the QM.

EN - INSTRUCTION FOR USE	3
FR - MODE D'EMPLOI	5
NL - GEBRUIKSAANWIJZING	8
ES - INSTRUCCIONES DE USO	11
DE - GEBRAUCHSANWEISUNG.....	14
PT - INSTRUÇÕES DE USO.....	17
IT - ISTRUZIONI PER L'USO	20
EL - ΟΔΗΓΪΕΣ ΧΡΪΗΣΗΣ	23
AR العربية - تعليمات الاستخدام.....	26
TR - KULLANIM TALIMATLARI	29

EN - INSTRUCTION FOR USE

Intended Use

Sclerotherapy & Endoscopic Needles are intended to introduce a sclerosing agent (any therapeutic vasoconstrictor agents, Histoacryl®, Cyanoacrylate, epinephrine) into selected sites to control actual or potential bleeding lesions in the digestive system or to inject liquid to aid in polypectomy procedures.

Product Description

Sclerotherapy & Endoscopic Needles are made of a flexible catheter with an inner tube connected to the needle at the distal end and a handle with a piston and a luer-lock connection at the proximal end.

Advices before first application

Please read carefully and follow all safety operating instructions and warnings before first application of the device.

A previous knowledge regarding handling and operation is required and essential.

Unpack the instrument carefully and examine the devices for any possibility of damage. In case of any damage or missing items contact your distributor immediately.

If the package present any damaged, the sterility of the device is not guarantee.

Handling and Operation

Check to ensure that the needle is operating properly by pushing the needle piston in and out while holding the handle or catheter. The handle is equipped with two stops. You should feel a “click” when the needle is in the fully deployed position and when the needle is in the fully retracted position.

Sclerotherapy & endoscopic needles are composed by a catheter and should never be manipulated with its catheter wined as this may damage the device and make its usage impossible!

For ethanolamine and cyanoacrylate only use a 19GA needle (GF262303)

Procedure

For Sclerotherapy and Endoscopic Needles, please proceed as follows

- Prepare a Luer-Lock adaptable syringe with the agent to be injected.
- Make sure that the needle is in its fully exposed position. Then attach the syringe to the needle via the Luer-Lock connection at the top of the handle and flush out the air in the piston using the liquid agent. Visually verify that the liquid agent reaches the distal tip of needle and all air has been evacuated from the inner catheter.

- Now retract the needle into its catheter to full withdrawn position. (The inner piston part of the handle is retracted to the position UP)
- Then introduce the catheter into the working channel of an appropriately sized endoscope. Advance in small increments until the catheter emerges from the scope.
- Fully deploy the needle until you feel the piston “click” into place. (The inner piston part of the handle is fully pushed into the outer part of the handle - Position DOWN)
- Proceed with the injection at an oblique angle to the tissue site
- Following the injection(s), retract the needle until the needle piston “clicks” into the fully retracted position. (Position UP)

When using the irrigation flush port, the needle must be fully retracted until the piston “click” into the retracted position. Prepare a luer-lock adaptable syringe with the irrigation agent to be injected and attach the syringe to the irrigation flush port via the luer-lock protruding from the side of the lateral connection port.

When inserting or withdrawing device from scope, ensure that the distal end of the scope is not in the retroflexed position.

- Upon completion of the procedure, remove the instrument from the scope.

When the device is being advanced or withdrawn from the endoscope, it is important to ensure that the needle is in its fully retracted position to avoid damaging the working channel of the endoscope. (The inner piston part of the handle is retracted to the position UP)

Contraindication

- Fever, allergic reaction to medication, perforation, sepsis
- Respiratory depression or respiratory arrest, cardiac arrhythmias, or halt
- Patients with lesions inappropriate for injection therapy

Storage

This device is delivered sterilized and in order to maintain this state until its first use it must be kept in its original sterilized package.

Do not put any objects on the instrument or its package!
Do not store the instruments near aggressive chemical products!

Do not expose the instruments to direct or indirect sunlight or other ultra-violet rays!

Keep in dry area

Complaints will not be taken under consideration if the instruments have been stored improperly.

Cleaning and Sterilization

This device is designed and warranted for single use!

Risk in case of reuse

This device is for single use only. G-Flex denied all responsibility in case of reuse. The reuse of single use device might present risk to the safety of the user or patient, due to possible uncontrolled contamination and / or unreliability of the device performance.

User

The users of G-FLEX instruments must be specialists in their fields. An appropriate and specific training for preparation, care and maintenance of the flexible instruments is required.

Repairing and complaints

In health protection of our employees only disinfected or sterilized instruments will be accepted for analysis or repairing. In any case disinfection or sterilization dates and validity must be labelled outside the package. If this requirement has not been fulfilled the instruments will be returned to the sender without analysis or repairing.

Product disposal

After use, these products may be a potential biohazard. Handle and dispose of in accordance with accepted medical practice and applicable local, state and federal laws and regulations.

Legal foundation

The Law of the European Union applied.

Product support

In case of questions or difficulties concerning our instruments please contact your local distributor or G-FLEX directly during regular working hours.

Times available: Monday to Friday; 9am to 5pm (MET)

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