EN - DISPOSABLE HEMOCLIP

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NL - WEGWERP-HEMOCLIP

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DESECHABLE

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IT - CLIP PER EMOSTASI MONOUSO

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

EN - INSTRUCTION FOR USE

Intended Use

Disposable Hemoclip are intended to place clip within the Gastrointestinal (GI) tract for the following purposes:

- Endoscopic marking;
- Hemostasis:
 - Mucosal/sub-mucosal defects (< 3 cm)
 - Bleeding ulcers
 - Arteries (< 2 mm)
 - Polyps (< 1.5 cm in diameter)
 - Diverticula in the colon
- Anchor to attached jejunal feeding tubes to the wall of the small bowel;
- Anchor to self expandable metallic stents (SEMS) to prevent migration;
- As a supplementary method to close luminal perforations (< 20 mm) in the GI tract that can be treated conservatively

Product Description

The U2525 series device is made of a pre-loaded, rotatable, radiopaque, single-use clip (See picture on page 3):

1. Clip jaws
2. Clip body

And a flexible delivery system made of:

3. Spring tube
4. Hole in the slider
5. Slider
6. Locking system in the
7. Handle
8. Thumb ring

The U25(11-16) series device is made of a pre-loaded, rotatable, radiopaque, single-use clip and a flexible delivery system made of (See picture on page 3):

1. Thumb ring
2. Handle
3. Slider
4. Rotation Wheel
5. Spring tube
6. Clip seat
7. Clip body
8. Clip jaws

The handle in the proximal part to manipulate the rotation and the deployment, a long flexible tube and the clip at the distal part. The handle is also used for opening and closing of the clip. It is designed to work in gastroscope and colonoscopes (see model working length) with a working channel equal to or greater than 2.8 mm diameter.

Intended Population

Our device is intended for both adult and paediatric patient according to the physician recommendation taking into account the contraindications.

Advices before first application

Please read carefully and follow all safety operating instructions and warnings before first application of the device.

A previous knowledge regarding handling and operation is required and essential. Unpack the instrument carefully and examine the devices for any possibility of damage. In case of any damage or missing items contact your distributor immediately.

If the package present any damaged, the sterility of the device is not guarantee.

Handling and Operation

Ensure the device and its packaging are intact before use. Remove the device from it's packaging before use.

Please observe compatibility between the operating channel and device diameter. Before inserting the device into the biopsy channel ensure that the clip is closed.


Procedure

Insertion into the Biopsy Channel

- Insert the device into the biopsy channel until it comes out of the distal part of the endoscope.
- Once under direct visualization, you are ready for the next step

Deployment of the clip


- The clip is deployed closed. To open the clip, move the slider distally (away from the thumb ring)(see picture below).

Rotation of the clip


- **U2525 Series** - To rotate the clip, in order to adjust the angle between the clip and the lesion, just rotate the complete handle in any direction while holding the blue sheet still. (See picture below)

- **U25(11-16) Series** - To rotate the clip, in order to adjust the angle between the clip and the lesion, just rotate the complete handle in any direction. (See picture below)


Closing the clip


- To close the clip jaws on the desired location move the slider proximally (close to the thumb ring) until the first tactile resistance is felt in the handle (see picture below).

- **U2525 Series** - To permanently deploy the clip, continue moving the slider proximally (toward the thumb ring) beyond the tactile resistance point until the Hole in the slider reaches the Locking System. At this point, click may be heard or felt. The Locking System ensures that the tracking wire stays inside the spring tube (spiral) and does not damage the endoscope biopsy channel. (see picture below).

- **U25(11-16) Series** - To permanently deploy the clip, continue moving the slider proximally (toward the thumb ring) beyond the tactile resistance point. At this point, click may be heard or felt. (see picture below).

Removal of the Delivery System

- Remove the whole delivery system from the biopsy channel and check the result with the endoscope, if needed repeat the whole procedure with other clips.

Contraindications

- Do not use this device when hemostasis cannot be verified visually with an endoscopic field of view.
- Do not perform MRI procedures on patients who have had clips placed within their gastrointestinal tract, as this could be harmful to the patient.
- Arteries greater than 2 mm.
- Polyps greater than 1.5 cm in diameter.
- Mucosal / Sub- Mucosal defects greater than 3 cm.

Possible Complications

A possible complication is re-bleeding.

Storage

This device is delivered sterilized and in order to maintain this state until its first use it must be kept in its original sterilized package.

Do not put any objects on the instrument or its package!
Do not store the instruments near aggressive chemical products!

Do not expose the instruments to direct or indirect sunlight or other ultra-violet rays!

Keep in dry area

Complaints will not be taken under consideration if the instruments have been stored improperly.

Cleaning and Sterilization

This device is designed and warranted for single use!

Risk in case of reuse

This device is for single use only. G-Flex denied all responsibility in case of reuse. The reuse of single use device might present risk to the safety of the user or patient, due to possible uncontrolled contamination and / or unreliability of the device performance.

User

The users of G-FLEX instruments must be specialists in their fields. An appropriate and specific training for preparation, care and maintenance of the flexible instruments is required.

Serious Incident, Complaints and repairs

In case of a serious incident, please report to the manufacturer and to the competent authority of the European Member State where the user is established. Shall the incident occur outside the Europe, please report to the manufacturer and follow the local applicable regulation.

In health protection of our employees only disinfected or sterilised instruments will be accepted for analysis or repairing. In any case disinfection or sterilisation dates and validity must be labelled outside the package. If this requirement has not been fulfilled the instruments will be returned to the sender without analysis or repairing.

Declarations

We, hereby, declare that the medical device neither contain nor is produced with tissue of animal origin.

We, hereby, declare that the medical device neither contain nor is produced with tissue or extracted material of human origin.

We, hereby, declare that the medical device do not contain medicinal substance.

We, hereby, declare that the medical device do not contain human blood derivative.

We, hereby, declare that the medical device do not contain phthalates or other products that could leak from the device that may be carcinogenic, mutagenic or toxic.

We, hereby, declare that the medical device is not composed of substances or combination of substances intended to be introduced into the human body.

Product disposal

After use, these products may be a potential biohazard. Handle and dispose of in accordance with accepted medical practice and applicable local, state and federal laws and regulations.

Legal foundation

The Law of the European Union applied.

Product support

In case of questions or difficulties concerning our instruments please contact your local distributor or G-FLEX directly during regular working hours.

Times available: Monday to Friday; 9am to 5pm (MET)

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