EN - MULTIBAND LIGATOR
FR - LIGATEUR MULTIBANDES
NL - MULTIBAND-LIGATOR
ES - LIGADOR DE MULTIBANDAS
DE - MULTIBAND-LIGATUR-SET
PT - KIT DE LIGADURA
IT - DISPOSITIVO PER LEGATURA
MULTIPLA
EL - ΣΥΣΚΕΥΉ ΑΠΟΛΊΝΩΣΗΣ
ΠΟΛΛΑΠΛΊΩΝ ΔΑΚΤΥΛΊΩΝ
AR - رابط متعدد الأربطة
TR - ÇOK BANTLI LIGATÖR
SV - MULTIBAND LIGATOR





G-Flex Europe S.P.R.L. - 20, Rue de l'industrie 1400
Nivelles - Belgium | Phone: +32 67.88.36.65
Email: info@g-flex.com | Website: www.g-flex.com



EN - Versions	Standard & Reloadable
FR - Versions	Standard & Rechargeable
NL - Versies	Standaard & Reloadable
ES - Versiones	Estándar & Recargable
DE - Versionen	Standard & Nachladbar
PT - Versões	Padrão & Recarregável
IT - Versioni	Standard e Ricaricabile
EL - Εκδόσεις	Κανονική και Με επαναφόρτωση
AR - الإصدارات	قياسي وقابل لإعادة التحميل
TR - Modeller	Standart ve Tekrar Yüklenebilir
SV - Versioner	Standard och omladdningsbar

EN - INSTRUCTION FOR USE	6
FR - MODE D'EMPLOI.....	12
NL - GEBRUIKSAANWIJZING.....	18
ES - INSTRUCCIONES DE USO	25
DE - GEBRAUCHSANWEISUNG.....	31
PT - INSTRUÇÕES DE USO	38
IT - ISTRUZIONI PER L'USO	44
EL - ΟΔΗΓΊΕΣ ΧΡΉΣΗΣ.....	50
AR العربية - تعليمات الاستخدام.....	57
TR - KULLANIM TALIMATLARI.....	64
SV - BRUKSANVISNING	70

AUTOCLAVABLE
134°C
273°F

EN - INSTRUCTION FOR USE

Intended Use

The Multiband Ligator is used via endoscope to ligate oesophageal varices at and above the oesophageal junction. The Multiband Ligator is supplied non-sterile and is intended for single use only.

The Ligating Ring is intended to ligate oesophageal varices at and above the oesophageal junction or to ligate anorectal haemorrhoids. The ligating ring is supplied non-sterile and is intended for single use only.

Please read carefully the below explanation for the reusable model that is composed of reusable handle which is used in combination with a disposable part

Product Description

All versions are composed of the following parts (see pictures on page 3):

1. Big wheel
2. **Multiband Ligator Handle**
3. Blue Luer-lock cap
4. Silicone sealing ring
5. **Multiband Ligator tip**
6. Alignment arrows
7. Small wheel
8. Tracking wire
9. Plastic part of the introduction catheter
10. **Introduction catheter**
11. Tracking wire inside the introduction catheter
12. Latex or Latex-Free ring
13. Penultimate ring in another color.
14. **Band Barrel (Cap)**
15. Soft part of the cap
16. Loop of the band barrel
17. Hook
18. **Dismounted handle (Reusable)**
19. Hole in the handle
20. Locking clip
21. **Reloading Unit (Disposable)**

The Ligating Ring is composed of a latex or latex-free ring

Intended Population

Our device is intended for only adult patient according to the physician recommendation taking into account the contraindications.

Advices before first application

Please read carefully and follow all safety operating instructions and warnings before first application of the device.

A previous knowledge regarding handling and operation is required and essential.

Unpack the instrument carefully and examine the devices for any possibility of damage. In case of any damage or missing items contact your distributor immediately.

Handling and Operation

This device is delivered non sterile and must **NOT** be sterilised before use!

This device **CANNOT** be used with an over tube.

For the ligating ring please refer to the application system's manufacturer instruction.

Mounting the Multiband Ligator:

1. Unpack and check all parts of the system carefully. Ensure that all parts are in the package.
2. The tracking wire of the Multiband Ligator is already pre-mounted into an introduction catheter.
3. Straighten the endoscope. We recommend the use of a biopsy valve. Keep its lid open.
4. Insert the catheter into the biopsy channel until it comes out of the distal tip of the scope.
5. Fix the handle inside the biopsy channel port.

Hold the **distal tip** of the introduction catheter and **remove the plastic part by slowly pulling it out of the endoscope until only tracking wire is left for the next fixing step.**

Do not roll up the catheter on your hands to remove as it will block the inner tracking wire and then prevent the catheter to go out from the biopsy channel

6. Take the band barrel loop and connect it with the tracking wire loop properly as shown on page 3, making a flat knot.

In order to tighten the knot, pull only the wires instead of pulling the whole cap as this could misfire a ring before the procedure.

7. Ensure that the handle has the alignment arrows facing each other.
8. Using the small wheel, **carefully and slowly** pull the band barrel close to the tip of the scope.
9. Mount the band barrel on the tip of the scope and ensure that the soft part is properly fixed on the scope in a way that it cannot fall off.
10. Rotate the band barrel in order to prevent the tracking wire being in front of the endoscope camera.
11. Now bend the endoscope in a J shape as to reproduce the position of introduction of the endoscope into the patient.
12. Slowly tension the tracking wire by turning the small wheel until you feel a small tension.
13. Verify the tension is enough by either looking at the monitor if the endoscope is already connected to the tower and lit up, or looking directly into the cap through the distal tip of the endoscope if the same is not connected to the tower. The wires must be in a straight position, tensioned.

Do not tension beyond this slight tension point to avoid shooting a ring.

If after inserting the handle of the device into the working channel entry port the handle is not stable enough or can move with soft movements of the scope, please remove the biopsy valve along with the handle of the device by 2cm up from the endoscope. Then press the handle further into the biopsy valve and place both handle and biopsy valve together into the endoscope working channel entry port again.

Procedure

Ligation of oesophageal varices

1. Lubricate the scope and exterior part of the band barrel ensuring that no lubricant has entered the barrel.

- Roll back the big wheel 180° to permit that the tracking wires is not fully straightened so that no band is shot while inserting the endoscope into the patient.
- When in place, tension again the tracking wire using the small wheel to prepare for the shooting (be sure that the alignment arrows on the handle are facing each other).
- Locate the selected varix and aspirate it into the band barrel.
- To fire the band, slowly rotate 180° forward the big wheel. More than one ligation band for each varix might be needed to control acute bleeding.
- Release the suction button of the scope, insufflate air and then withdraw the scope slightly to release the ligated varice.
- Repeat the procedure until the last varix is banded.

 The luer-lock connection is available to ensure irrigation of the scope biopsy channel if required. After removing the blue cap, attach the luer-lock connector to a syringe filled with sterile water and irrigate.

 After irrigation put back the blue cap to ensure good suction.

 In case additional bands are needed, remove the scope and attach a new Multiband Ligator starting back from the mounting steps.



Removing the Multiband Ligator

- Once the ligation procedure is completed, remove the scope from the patient.
- Dismantle the the Multiband Ligator as follows:
 - If any unused bands are still on the barrel shoot all remaining bands
 - Remove the Multiband Ligator handle from the biopsy channel by pulling it, the remaining tracking wire attached with the handle should be removed also.
 - Remove the band barrel from the scope tip.
 - Dispose all parts per institutional guidelines for biohazardous medical waste.

 If all bands have been used, you can remove the Multiband Ligator from the biopsy channel and continue using the scope and its biopsy channel while still inside the patient.

Reloading the Multiband Ligator (Reloadable version ONLY)

- Once the complete system is out of the patient and removed as described above, you may start the reloading process.
- Loosen completely the tracking wire by turning backwards the big wheel
- Remove the blue luer-lock cap
- Press the locking clip towards the center of the handle.
- Remove the used L-shape reloading unit
- Remove the blue luer-lock cap from the new reloading units
- Insert the new L-shape reloading unit through the hole of the handle (see picture aside). You should hear a "click" once it is correctly in place.
- Put back the blue luer-lock cap
- Attach the tracking wire loop through the shaft hook as indicated aside
- While holding with your thumb the tracking wire inside the hook, turn the small wheel forward until the tracking wire is slightly tensioned
- The reloaded unit is now ready for use and you may follow the instructions from "Mounting the Multiband Ligator"

Ligating Ring

 Please follow the ligator manufacturer instruction for the usage of the ligating ring.

Contraindication

The contraindications include without being limited to cricopharyngeal or oesophageal narrowing/ stricture, tortuous oesophageal diverticula, known or suspected oesophageal perforation, asymptomatic rings or webs, coagulopathy. Use of ligation bands is contraindicated in patients with a known hypersensitivity to latex.

Ligation of anorectal hemorrhoids is contraindicated in severe inflammatory bowel disease and portal hypertension. Hemorrhoidal banding may result in severe pain if the procedure is performed below the dental line.

Possible Complications

Possible complications associated with the gastrointestinal endoscopy include and not limited to perforations, fever, aspirations, haemorrhage, infections, hypo-tension, allergic reactions to medication, cardiac arrhythmia or arrest, respiratory depression or arrest. Furthermore, complications that may occur with oesophageal banding include and not limited to chest tightness, retrosternal pain, dysphagia or odynophagia, and post-EVL ulcer with or without bleeding, increase of portal hypertensive gastropathy, bacteremia, nausea, laryngeal, retrosternal pain, laceration, stricture formation, oesophageal perforation, obstructions.

Also, complications that may occur with anorectal haemorrhoids and not limited to severe pain, bleeding, urinary symptoms, swelling and edam, tissue ulceration, bands slipping off. Hemorrhoidal banding may result in severe pain if the procedure is performed below the dental line.

Precautions

It is vital the biopsy channel size of the endoscope is co-ordinated with compatible devices in order to obtain the best results during the procedure. The minimum biopsy channel diameter required is 2.8 mm and length of up to 150 cm. The outer diameter of the scope to be used must be from 8,5 mm up to 11,5 mm. Band ligation may not be effective when applied to small varices.

It is vital to ensure the understanding of the technical principles, clinical applications and associated risks with GI endoscopy as well as endoscopic banding before using this product. This product must only be used by or under direct supervision of physician thoroughly trained in therapeutic endoscopy as well as vessel banding.

Oesophageal ligators are not intended for ligation of varices below the gastro-oesophageal junction.

Banding should begin at the gastro-oesophageal junction and proceed up the oesophagus in an helicoidal manner, respecting a minimum of 2cm in between 2 ligatures. Passing the endoscope over an already ligated varices may dislodge the band from its position.

Routine endoscopic examination by the physician is recommended to confirm the diagnosis requiring treatment of oesophageal varices prior to assembling the Multiband Ligator.

Storage

Do not put any objects on the instrument or its package!
Do not store the instruments near aggressive chemical products!

Do not expose the instruments to direct or indirect sunlight or other ultra-violet rays!

Keep in dry area.

Complaints will not be taken under consideration if the instruments have been stored improperly.

Cleaning and Sterilization

This device is designed and warranted for single use!

For the reusable version, only the handle of this device can be reprocessed, please consult the following:

G-Flex recommend to follow one of the following sequences when sterilizing the device:

- Manual Cleaning, Ultrasonic Cleaning, Automatic/Thermal Cleaning and Autoclave; or
- Manual Cleaning, Ultrasonic Cleaning, Automatic/Thermal Cleaning and Ethylene Oxide Gas Sterilization

Either sterilize by Autoclave or by ETO. There is no need to do both.

Manual Cleaning

The cleaning of the device should be done immediately after each use:

- Immerse the device in a suitable liquid detergent or disinfectant (Always observe the minimum diameter of 20 cm when winding the instrument and applicability of the cleaning product!).
- Wipe the instrument with a soft towel.

Ultrasonic Cleaning

The ultrasonic cleaning facilitates the dislodging of residual material and must be done before sterilization. The ultrasonic cleaning must last for at least 30 minutes.

- Clean the device right after its use.
- Put the instrument in an ultrasonic cleaning tank always observing the winding minimum diameter 20 cm (Please observe the ultrasonic equipment manufacturer's instructions!).
- Use only tap water and the detergents which dissolve albumen, always according to ultrasonic equipment manufacturer's instructions.
- Rinse the instrument with tap water and wipe dry by using gauze pads.

Automatic / Thermal Cleaning

Alternatively or additionally to the manual cleaning, this instrument may also be cleaned in automatic / thermal cleaning machines by using approved alkaline detergents. Please refer to relevant instructions of the cleaning machine supplier.

Sterilization - Autoclave

Before autoclaving sterilization the instrument should be cleaned carefully as described in Manual Cleaning, Ultrasonic Cleaning and Automatic/Thermal Cleaning.

This symbol means that the product can be sterilized by autoclave according to the below parameters

Parameters:

Temperature: 134 °C / 273°F

Pressure: 3 bar

Exposure: ≥ 18 min

Also in here the minimum 20 cm diameter for winding the instrument must be noted.

The original packages of G-FLEX are not autoclavable!

Sterilization - Ethylene Oxide Gas Sterilization

Please follow manufacturer's warnings regarding to the biological indicators

Parameters:

Temperature: 57 °C

Pressure:	1,7 bar
Exposure:	max. 4 h
Relative Humidity:	50 %
Gas Concentration:	12 %
Aeration Time:	7 days in ambient temperature or 12 hours by 50-57 °C

Risk in case of reuse

This device is for single use only. G-Flex denied all responsibility in case of reuse. The reuse of single use device might present risk to the safety of the user or patient, due to possible uncontrolled contamination and / or unreliability of the device performance.

User

The users of G-FLEX instruments must be specialists in their fields. An appropriate and specific training for preparation, care and maintenance of the flexible instruments is required.

Serious Incident, Complaints and repairs

In case of a serious incident, please report to the manufacturer and to the competent authority of the European Member State where the user is established. Shall the incident occur outside the Europe, please report to the manufacturer and follow the local applicable regulation.

In health protection of our employees only disinfected or sterilised instruments will be accepted for analysis or repairing. In any case disinfection or sterilisation dates and validity must be labelled outside the package. If this requirement has not been fulfilled the instruments will be returned to the sender without analysis or repairing.

Declarations

We, hereby, declare that the medical device neither contain nor is produced with tissue of animal origin.

We, hereby, declare that the medical device neither contain nor is produced with tissue or extracted material of human origin.

We, hereby, declare that the medical device do not contain medicinal substance.

We, hereby, declare that the medical device do not contain human blood derivative.

We, hereby, declare that the medical device do not contain phthalates or other products that could leak from the device that may be carcinogenic, mutagenic or toxic.

We, hereby, declare that the medical device is not composed of substances or combination of substances intended to be introduced into the human body.

Product disposal

After use, these products may be a potential biohazard. Handle and dispose of in accordance with accepted medical practice and applicable local, state and federal laws and regulations.

Legal foundation

The Law of the European Union applied.

Product support

In case of questions or difficulties concerning our instruments please contact your local distributor or G-FLEX directly during regular working hours.

Times available: Monday to Friday, 9am to 5pm (MET)

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