EN - STENTS APPLICATION SYSTEM
FR - SYSTÈME D'APPLICATION D'ENDOPROTHÈSE
NL - STENTPLAATSINGS SYSTEEM
ES - INTRODUCTORES DE STENT
DE - STENTS-APPLIKATIONSSYSTEM
PT - SISTEMA DE APLICAÇÃO DE PRÓTESE
IT - SISTEMA DI APPLICAZIONE DE STENT
EL - ΣΥΣΤΗΜΑ ΕΦΑΡΜΟΓΉΣ ΤΩΝ
ΕΝΔΟΠΡΟΘΈΣΕΩΝ
AR - نظام تطبيق الدعامات
TR - STENTI UYGULAMA SISTEMI
SV - STENTAPPLICERINGSSYSTEM

G-Flex Europe S.P.R.L. - 20, Rue de l'industrie 1400
Nivelles - Belgium | **Phone:** +32 67.88.36.65
Email: info@g-flex.com | **Website:** www.g-flex.com

EN - INSTRUCTION FOR USE	3
FR - MODE D'EMPLOI	8
NL - GEBRUIKSAANWIJZING	13
ES - INSTRUCCIONES DE USO	19
DE - GEBRAUCHSANWEISUNG.....	25
PT - INSTRUÇÕES DE USO.....	31
IT - ISTRUZIONI PER L'USO	36
EL - ΟΔΗΓΪΕΣ ΧΡΪΗΣΗΣ	42
AR - العربية - تعليمات الاستخدام.....	48
TR - KULLANIM TALIMATLARI	53
SV - BRUKSANVISNING	58

EN - INSTRUCTION FOR USE

Intended Use

Stents Application System are intended to place stent via endoscopes.

Product Description

The application system in 8.5 Fr (color code: orange) or 10 Fr (color code: red) is made of an outer pusher over a guiding catheter that has a radiopaque metal ring at the distal extremity and a Y connector at the proximal end.

The application system in 5 Fr (color code: none) or 7 Fr (color code: yellow) is made of an outer pusher only with a radiopaque metal ring at the distal extremity.

If the application system is not presented with your stent, please consider having one on hands always respecting the French size of the stent. (Refer to the current catalogue for more information)

Intended Population

Our device is intended for only adult patient according to the physician recommendation taking into account the contraindications.

Advices before first application

Please read carefully and follow all safety operating instructions and warnings before first application of the device.

A previous knowledge regarding handling and operation is required and essential.

Unpack the instrument carefully and examine the devices for any possibility of damage. In case of any damage or missing items contact your distributor immediately.

Handling and Operation

The Endoscopic Drainage Stents are placed in conjunction with an endoscope, a guide wire, a pusher and, when necessary, a guiding catheter

To facilitate the application, an endoscopic sphincterotomy prior to guide wire placement is recommended. (Please observe HF-generator and sphincterotome manufacturer's instructions).

It is recommended NOT to remove the stent by inflating a balloon inside the stent.

The patient and the stent should be regularly followed-up to remove the stent immediately in case of clog, migration or fallout from the papilla.

Procedure

For stents in 5 Fr or 7 Fr, please proceed as follows:

- Choose the shape, diameter and length of the stent to be placed.

 Note: make sure to have on hands an application system with the same diameter as the stent selected (5 Fr or 7 Fr).

- Remove the stent with its positioning sleeve and the application system out of their packaging.
- Over the guide wire, place the tapered distal extremity of the stent, followed by the application system and the positioning sleeve over the application system until the stent reach the entrance of the biopsy channel.
- Push the positioning sleeve forward over the proximal flap to lower it in order to protect it while inserting the stent and the rest of the device inside the biopsy channel.

 Note the positioning sleeve will stay out of the biopsy channel.

- To place the stent in the desired location, hold the guide wire with one hand while pushing the application system forward with the other hand until the stent is in place.
- Once the stent is placed, pull the guide wire back into the biopsy channel while holding the application system in place. After the complete removal of the guide wire, the stent will regain its original shape.
- Finally, remove application system out of the scope.

For stents in 8.5 Fr or 10 Fr, please proceed as follows:

- Choose the shape, diameter and length of the stent to be placed.

 Note: make sure to have on hands an application system with the same diameter as the stent selected (8.5 Fr or 10 Fr).

- Remove the stent with its positioning sleeve and the application system out of their packaging.
- Insert the positioning sleeve from the distal end of the guiding catheter all the way up the proximal Y connector.
- Place the stent from its proximal extremity (NOT tapered side) over the guiding catheter until it reaches the distal extremity of the pusher.
- Insert the assembled stent and application system over the guide wire by inserting the guide wire inside the distal end of the guiding catheter until the tapered distal extremity stent reach the entrance of the biopsy channel.
- Push the positioning sleeve forward over the proximal flap to lower it in order to protect it while inserting the stent and the rest of the device inside the biopsy channel.

 Note the positioning sleeve will stay out of the biopsy channel.

- Advance the guiding catheter until the radiopaque ring passes the obstructed area.
- Once the guiding catheter is in the desired place, unscrew the luer-lock connection of the Y-connector to release the pusher from the guiding catheter.
- Slowly advance the pusher until it reaches the proximal extremity of the stent.

To hold the guide wire in place with the guiding catheter, use the sealing cap on the Y-connector. Injection of contrast is possible via the guiding catheter Y-connector lateral connection. Do this with the sealing cap close.

- To place the stent into its final position, hold the guiding catheter with one hand while slowly advancing the pusher.

Keep guiding catheters position while placing the stent

- Once the stent is placed, pull the guiding catheter and the guide wire back into the biopsy channel while keeping the pusher in its position. After the complete removal of those devices, the stent will regain its original shape.
- Finally, remove pusher, guiding catheter and guide wire out of the scope.

For stents removal, please proceed as follows:

- As the stent has been placed in the desired location with the proximal extremity appearing into the digestive system cavity, it can be removed by gently pulling it with a foreign body retriever forceps or a snare or even a basket.
- Once the stent has been pulled out remove it from the patient along with the scope and the removing device used.

Contraindication

The contraindications include without being limited to poor general condition such as respiratory failure, risk of aspiration, general blood coagulation and wound healing disorders, asymptomatic stenosis or jaundice, endoscopically impassible stenosis, acute, non-biliary pancreatitis, absent patient consent. The device is not indicating for paediatric procedures

Storage

This device is delivered non sterile and must be sterilised before first use according with the instruction described below.

Do not put any objects on the instrument or its package!
Do not store the instruments near aggressive chemical products!

Do not expose the instruments to direct or indirect sunlight or other ultra-violet rays!

Keep in dry area

Complaints will not be taken under consideration if the instruments have been stored improperly.

Cleaning and Sterilization

G-Flex recommend to follow one of the following sequences when sterilizing the device:

- Manual Cleaning, Ultrasonic Cleaning, Automatic/Thermal Cleaning and Autoclave; or
- Manual Cleaning, Ultrasonic Cleaning, Automatic/Thermal Cleaning and Ethylene Oxide Gas Sterilization

Either sterilize by Autoclave or by ETO. There is no need to do both.

Manual Cleaning

The cleaning of the device should be done immediately after each use:

- Immerse the device in a suitable liquid detergent or disinfectant (Always observe the minimum diameter of 20 cm when winding the instrument and applicability of the cleaning product!).
- Wipe the instrument with a soft towel.

Ultrasonic Cleaning

The ultrasonic cleaning facilitates the dislodging of residual material and must be done before sterilization. The ultrasonic cleaning must last for at least 30 minutes.

- Clean the device right after its use.
- Put the instrument in an ultrasonic cleaning tank always observing the winding minimum diameter 20 cm (Please observe the ultrasonic equipment manufacturer's instructions!).
- Use only tap water and the detergents which dissolve albumen, always according to ultrasonic equipment manufacturer's instructions.
- Rinse the instrument with tap water and wipe dry by using gauze pads.

Automatic / Thermal Cleaning

Alternatively or additionally to the manual cleaning, this instrument may also be cleaned in automatic / thermal cleaning machines by using approved alkaline detergents. Please refer to relevant instructions of the cleaning machine supplier.

Sterilization - Autoclave

Before autoclaving sterilization the instrument should be cleaned carefully as described in Manual Cleaning, Ultrasonic Cleaning and Automatic/Thermal Cleaning.

This symbol means that the product can be sterilized by autoclave according to the below parameters

Parameters:

Temperature: 134 °C / 273 °F
 Pressure: 3 bar
 Exposure: ≥ 18 min

Also in here the minimum 20 cm diameter for winding the instrument must be noted.

The original packages of G-FLEX are not autoclavable!

Sterilization - Ethylene Oxide Gas Sterilization

Please follow manufacturer's warnings regarding to the biological indicators

Parameters:

Temperature:	57 °C
Pressure:	1,7 bar
Exposure:	max. 4 h
Relative Humidity:	50 %
Gas Concentration:	12 %
Aeration Time:	7 days in ambient temperature or 12 hours by 50-57 °C

User

The users of G-FLEX instruments must be specialists in their fields. An appropriate and specific training for preparation, care and maintenance of the flexible instruments is required.

Serious Incident, Complaints and repairs

In case of a serious incident, please report to the manufacturer and to the competent authority of the European Member State where the user is established. Shall the incident occur outside the Europe, please report to the manufacturer and follow the local applicable regulation.

In health protection of our employees only disinfected or sterilised instruments will be accepted for analysis or repairing. In any case disinfection or sterilisation dates and validity must be labelled outside the package. If this requirement has not been fulfilled the instruments will be returned to the sender without analysis or repairing.

Declarations

We, hereby, declare that the medical device neither contain nor is produced with tissue of animal origin.

We, hereby, declare that the medical device neither contain nor is produced with tissue or extracted material of human origin.

We, hereby, declare that the medical device do not contain medicinal substance.

We, hereby, declare that the medical device do not contain human blood derivative.

We, hereby, declare that the medical device do not contain phthalates or other products that could leak from the device that may be carcinogenic, mutagenic or toxic.

We, hereby, declare that the medical device is not composed of substances or combination of substances intended to be introduced into the human body.

Product disposal

After use, these products may be a potential biohazard. Handle and dispose of in accordance with accepted medical practice and applicable local, state and federal laws and regulations.

Legal foundation

The Law of the European Union applied.

Product support

In case of questions or difficulties concerning our instruments please contact your local distributor or G-FLEX directly during regular working hours.

Times available: Monday to Friday; 9am to 5pm (MET)

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