EN - BIOPSY FORCEPS
FR - PINÇE À BIOPSIE
NL - BIOPSIETANG
ES - FORCEPS DE BIOPSIA
DE - BIOPSIEZANGE
PT - FÓRCEPS DE BIOPSIA
IT - PINZA PER BIOPSIA
EL - ΛΑΒΉΔΑ ΒΙΟΨΪΑΣ
AR - ملقط الخزعة
TR - BIYOPSI FORSEPSI
SV - BIOPSITÅNG
CS – BIOPTICKÉ KLEŠTĚ
DA - BIOPSITANG
HU – BIOPSZIÁS FOGÓ
RO - PENSĂ PENTRU BIOPSIE

G-Flex Europe S.R.L. - 20, Rue de l'industrie 1400
Nivelles - Belgium | Phone: +32 67.88.36.65
Email: info@g-flex.com | Website: www.g-flex.com

EN - Catalog number	GF120
FR - Numéro de catalogue	GF120-E
NL - Catalogusnummer	GF160
ES - Número de catálogo	GF16034
DE - Katalog-Nummer	GF160-E
PT - Número do catalogo	GF160P
IT - Numero di catalogo	GF160P-E
EL - Αριθμός καταλόγου	GF160PS
AR - العربية- رقم الكتالوج	GF160R-E
TR - Katalog numarası	GF160S
SV - Katalognummer	GF160S-E
CS - Katalogové číslo	GF180
DA - Katalognummer	GF180-E
HU – Katalógusszám	GF180S
RO - Număr de catalog	GF180S-E
	GF230
	GF23034
	GF230-E
	GF230P
	GF230P-E
	GF230S
	GF230S-E
	GU100
	GU160

EN - Instruction for use.....	4
FR - Mode d'emploi	9
NL - Gebruiksaanwijzing	15
ES - Instrucciones de uso	20
DE - Gebrauchsanweisung	25
PT - Instruções de Uso.....	31
IT - Istruzioni per l'uso	36
EL - Οδηγίες χρήσης.....	42
AR العربية - تعليمات الاستخدام.....	47
TR - Kullanım talimatları.....	53
SV - Bruksanvisning	58
CS – Návod k použití.....	63
DA - Brugsanvisning.....	68
HU – Használati útmutató	74
RO - Instrucțiuni de utilizare	80

EN - INSTRUCTION FOR USE

1. Intended Use

Biopsy Forceps are used for taking out biopsy specimens (tissues samples) from the human body, through endoscopes.

2. Intended Population

Our device is intended for both adult and paediatric patient according to the physician recommendation taking into account the contraindications.

3. Contraindications

All the usual contraindications for the use of biopsy forceps apply to the use of G-Flex disposable biopsy forceps. These include:

Related to the cardiovascular system:

Shock, Hypoxia, Comorbidity (e.g. severe cardiopulmonary diseases and decompensation), Coagulopathy, Uncontrollable hemorrhagic diathesis, Acute myocardial infarction, Tissue samples from vessel injury in ulceration, Use on the central cardiovascular system.

Related to the gastrointestinal system:

Peritonitis, Fragility of the intestinal wall (e.g. ulcerative colitis, diverticulitis, colitis ulcerosa, toxic megacolon), Fulminant colitis, Ileus, Not prepped bowel, Recent gastrointestinal anastomosis, Acute abdomen.

General contraindications:

Non fasted patient, Severe respiratory insufficiency, Acute perforation of the tracts, Primary endoscopy, colonoscopy, rectoscopy; Sepsis, Pregnancy. Moreover, an analysis of the risk/benefit ratio should be established by endoscopist and cardiologist for patients using anticoagulants or anti-aggregants.

4. Possible Complications

Possible complications are perforation, bleeding, difficulties to remove/insert, mechanical issues, bent device issues, breakage issues.

5. User

The users of G-FLEX instruments must be specialists in their fields. An appropriate and specific training for preparation, care and maintenance of the flexible instruments is required.

6. Product Description

The device is composed by a jaw (in different shapes and sizes), a spiral spring hose (in different lengths and diameters) and an operating handle.

7. Advices before first application

Please read carefully and follow all safety operating instructions and warnings before first application of the device.

A previous knowledge regarding handling and operation is required and essential. Unpack the instrument carefully and examine the devices for any possibility of damage. In case of any damage or missing items contact your distributor immediately.

Considering that device are not delivered sterilized the user must do it before the first application (See Cleaning and Sterilization)

8. Handling and Operation

The device must be passed through the endoscope operating channel with its cups closed. Please observe compatibility between the biopsy forceps and the endoscope operating channel diameter.

Unwind the biopsy forceps, manipulate backward - forward smoothly the plastic spool of the handle to verify the opening / closing of the cups.

9. Procedure

- Approach the biopsy forceps near to the tissue (lesion) to be taken from the patient's body.
- Open the cups by pushing the plastic spool at the handle.
- Place the forceps on the tissue and take the sample gently by pulling back the plastic spool (this manipulation will close the forceps cups).
- Keep the forceps cups closed and remove biopsy forceps out of the endoscope carefully.
- This tissue specimen taken will be ready for further histological examination.

10. Precautions/ Warnings

Biopsy Forceps are composed by a spiral and should never be manipulated with its catheter winded as this may damage the device and make its usage impossible!

Never clean instruments with metal brushes or wire wool

11. Cleaning and Sterilization

Before the first use and prior to all subsequent uses, the device must be cleaned, disinfected and sterilised. Effective cleaning and disinfection are essential for an effective sterilisation.

The reprocessing of reusable endoscopic equipment must remain under the responsibility of individuals with thorough training in reprocessing methods and be aware of the following points:

- Local hospital's reprocessing guidelines
- National health and safety guidelines
- The instructions in this instruction for use

- Additional national specifications related to effective inactivation of prions (i.e: use of highly alkaline cleaning agent and prolonged sterilisation time)
 - That the methods used to clean, disinfected and sterilise are adequately validated for this type of instruments
 - That the equipment used is regularly maintained and checked
 - That the validated parameters are maintained during each reprocessing cycles
- The reprocessing of the instrument must be initiated immediately after use, since dried organic residue can impact effective cleaning, disinfection and sterilisation.

Reprocessing methods

A mechanical cleaning/disinfection process should be used whenever possible over manual cleaning.

Therefore, one of the following 2 options must be used for the reprocessing:

1. Perform the pre-treatment and manual cleaning before the mechanical cleaning/disinfection followed by the sterilization (3 steps).
2. In case a mechanical cleaning/disinfection is not available the device must then be pre-treated and manually cleaned followed by a sterilization (2 steps).

Material Resistance

Ensure that the cleaning and disinfection agent do not contain the following:

- Aldehydes (this will fix blood-based soiling)
- Strong organic, mineral or oxidising acids (lowest permissible pH value 5.5)
- Strong alkalis (highest permissible pH 11, neutral / enzymatic or slightly alkaline cleaners are recommended)
- Organic solvents (e.g. ether, ketones, benzines), fluorinated alcohols
- Oxidising agents (e.g. hydrogen peroxide)
- Halogens (chlorine, iodine, bromine)
- Aromatic / halogenated hydrocarbons
- Formamide
- Trichloroethylene / perchloroethylene

Pre-treatment / Manual Cleaning

1. Immersed the instrument in cleaning solution within maximum 2 hours after use
2. Remove any soiling visible on the jaws using a clean soft brush (or a clean, soft and lint-free cloth) which is used solely for this purpose. Never use metal brushes or wire wool.
3. Open and close the jaws, immersed in cleaning solution at least 3 times
4. Immediately place the instrument wined with a diameter not less than 20 cm into an ultrasound bath filled with cleaning solution for 30 minutes (time and concentration may vary according to the solution manufacturer's specifications. Cleaning method has been validated with ultrasound bath using 0.5% neodisher® Mediclean forte (Dr. Weigert GmbH) at 35 kHz)
5. Rinse the instrument for at least 1 minute under running water (temperature < 35°C/85°F)

Mechanical Cleaning / Disinfection

Proceed as followed:

1. Place the instrument wined with a diameter not less than 20 cm into the cleaning and disinfection equipment that fulfils the following requirement:
 - Equipment is effective and suitable for the cleaning and disinfection of medical instruments (i.e: FDA approved, CE marked)
 - Only sterile or low-germ (max. 10 bacteria/ml) and low-endotoxin (max. 0.25 endotoxin units /ml) water is used (e.g. purified water/highly purified water)
 - The air used for drying is filtered (oil free, low-germ and low-particle)
 - The device is maintained and checked regularly
 - The cleaning and disinfection solutions used are compatible with the instrument (See "Material Resistance")
 - The concentration, temperatures and exposure times specified by the cleaning / disinfectant agent manufacturer is followed at all time, this includes the rinsing specifications
2. Start the cleaning and disinfection program:
 - The programme must be suitable for the instrument
 - Thermal disinfection programme should be used (90°C / 194°F for at least 5 minutes)
 - The program contains at least 3 rinsing cycles (depleting steps) after cleaning (or neutralization if applied) or conductance control recommended to effectively prevent detergent residues
3. Remove the instruments from the equipment when the program is finished
4. Check the functionality of the device by opening and closing the jaws several times
5. Check the device for signs of corrosion, damaged surfaces, cracks, soiling or discolouration
6. Remove any damaged device from use
7. Any soiled device must be cleaned and disinfected again
8. Package the device into a single-use sterilisation pack or a sterilisation container (suitable for steam sterilisation of at least 134°C / 270 °F) as soon as possible after removal from the equipment

Sterilisation

Only steam sterilisation must be used.

Flash sterilisation, hot air sterilisation, radiation sterilisation, formaldehyde, ethylene oxide sterilisation or plasma sterilisation are not permitted.

The following aspects shall be taken into account:

- Fractionated vacuum method (at least 3 vacuum steps and with adequate product drying not less than 20 minutes)
- Steam steriliser compliant with DIN EN 13060/DIN EN 285 or ANSI AAMI ST79
- Validated in accordance with DIN EN ISO 17665 (valid IQ/OQ and product-specific performance assessment (PQ))

- Maximum sterilisation temperature 134°C / 273°F plus tolerance as defined in DIN EN ISO 17665
- Sterilisation time (exposure time at sterilisation temperature) at least 4 min at 134°C/270°F).
- For prion inactivation (depending on the local recommendation) at least 18 minutes at 134°C/273°F

12. Storage

This device is delivered non sterile and in order to maintain this state until its first use it must be kept in its original package.

Do not put any objects on the instrument or its package!
Do not store the instruments near aggressive chemical products!

Do not expose the instruments to direct or indirect sunlight or other ultra-violet rays!

Keep in dry area

Complaints will not be taken under consideration if the instruments have been stored improperly.

13. Risk in case of reuse

The quantity of reprocessing cycles the device can withstand is closely related to the care given to the device during its intended use and the care taken during the cleaning and sterilization. We have validated 200 reprocessing cycles but the user can use for longer looking, in between each use, for signs of performance degradation such as kinks on the spiral, damaged to the jaws and stain. Shall any of these be present, it is recommended to segregate the device and use another one.

14. Serious Incident, Complaints and repairs

In case of a serious incident, please report to the manufacturer and to the competent authority of the European Member State where the user is established. Shall the incident occur outside the Europe, please report to the manufacturer and follow the local applicable regulation.

In health protection of our employees only disinfected or sterilised instruments will be accepted for analysis or repairing. In any case disinfection or sterilisation dates and validity must be labelled outside the package. If this requirement has not been fulfilled the instruments will be returned to the sender without analysis or repairing.

15. Declarations

We, hereby, declare that the medical device neither contain nor is produced with tissue of animal origin.

We, hereby, declare that the medical device neither contain nor is produced with tissue or extracted material of human origin.

We, hereby, declare that the medical device do not contain medicinal substance.

We, hereby, declare that the medical device do not contain human blood derivative.
We, hereby, declare that the medical device do not contain phthalates or other products that could leak from the device that may be carcinogenic, mutagenic or toxic.

We, hereby, declare that the medical device is not composed of substances or combination of substances intended to be introduced into the human body.

16. Product disposal

After use, these products may be a potential biohazard. Handle and dispose of in accordance with accepted medical practice and applicable local, state and federal laws and regulations.

17. Legal foundation

The Law of the European Union applied.

18. Product support

In case of questions or difficulties concerning our instruments please contact your local distributor or G-FLEX directly during regular working hours.

Times available: Monday to Friday; 9am to 5pm (MET)

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