EN - CYSTO-GASTRO SET
FR - SET CYSTO-GASTRO
NL - CYSTE-MAAG-SET
ES - SET DE GASTRO-QUISTE
DE - CYSTO-GASTRO SET
PT - CONJUNTO CISTO-GÁSTRICO
IT - SET CISTO-GASTRO
EL - ΣΕΤ ΓΙΑ Κ΄ΥΣΤΕΙΣ ΚΑΙ
ΓΑΣΤΡΟΣΚ΄ΟΠΗΣΗ
AR - مجموعة CYSTO-GASTRO
TR - SISTO-GASTRO SET
SV – CYSTO-GASTRO-SET
CS – CYSTO-GASTRO SETY
DA - CYSTO-GASTRO-SÆT
HU – CYSTO-GASTRO KÉSZLET
RO - CYSTO-GASTRO SET

G-Flex Europe S.R.L. - 20, Rue de l'industrie 1400
Nivelles - Belgium | Phone: +32 67.88.36.65
Email: info@g-flex.com | Website: www.g-flex.com

Current Name: IFU-Cysto-Gastro Set



Creation Date: 2024-12-08 Version 22

This document is only controlled in electronic format by the QM.

EN - Catalog number
FR - Numéro de catalogue
NL - Catalogusnummer
ES - Número de catálogo
DE - Katalog-Nummer
PT - Número do catalogo
IT - Numero di catalogo
EL - Αριθμός καταλόγου
AR - العربية- رقم الكتالوج
TR - Katalog numarası
SV - Katalognummer
CS - Katalogové číslo
DA - Katalognummer
HU – Katalógusszám
RO - Număr de catalog

CYSTO06U
CYSTO08UK
CYSTO10U
CYSTO10UK
CYSTO85U

EN - INSTRUCTION FOR USE	5
FR - MODE D'EMPLOI	8
NL - GEBRUIKSAANWIJZING	13
ES - INSTRUCCIONES DE USO	17
DE - GEBRAUCHSANWEISUNG.....	21
PT - INSTRUÇÕES DE USO.....	25
IT - ISTRUZIONI PER L'USO	29
EL - ΟΔΗΓΊΕΣ ΧΡΉΣΗΣ	33
AR - العربية - تعليمات الاستخدام.....	37
TR - KULLANIM TALIMATLARI	41
SV - BRUKSANVISNING	45
CS - NÁVOD K POUŽITÍ.....	49
DA - BRUGSANVISNING	53
HU – HASZNÁLATI ÚTMUTATÓ	57
RO - INSTRUCȚIUNI DE UTILIZARE	61

EN - INSTRUCTION FOR USE

1. Intended Use

Cysto-Gastro-Sets are intended to be used to puncture pancreatic cysts or pseudo-cysts endoscopically (via the transgastric or transduodenal wall) and to dilate the opening to the cyst (6, 8.5 or 10Fr models) as an alternative to surgical or percutaneous treatments.

The claimed clinical performance is that the Cysto-Gastro Sets allow an effective cyst puncture and dilation prior the placement of drainage devices (stents) - (success rate >90%).

2. Intended Population

Our device is intended for only adult patient according to the physician recommendation taking into account the contraindications.

3. Contraindications

Contraindications include blood coagulation diseases, interposing vessels, nickel sensitivity and contraindications related to the use of electrical scalpel.

4. Possible Complications

Possible complications associated with the gastrointestinal endoscopic puncture of cysts or pseudo-cysts include pain, abscess, hemorrhage/bleeding, pancreatitis, perforation, small risk of infection, abdominal pain, coagulation disorders and contraindications related to the use of electrical scalpel, allergic reaction to medication, allergic reaction to nickel, hypotension, respiratory depression or arrest, cardiac arrhythmia or arrest.

5. User & Environment

The users of G-FLEX instruments must be specialists in their fields. An appropriate and specific training for preparation, care and maintenance of the flexible instruments is required.

This device is intended to be used in operating room in hospitals or clinics.

6. Product Description

The device is made of an outer flexible teflon tube equipped with a metal tip at the distal end and an HF connector in the proximal end (all models) images on page 4:

1. HF dilator
2. Guiding catheter with diathermic needle

Handle

3. Luer-lock
4. Guide wire entry port
5. Dilator HF Connector

7. Advices before first application

Please read carefully and follow all safety operating instructions and warnings before first application of the device.

A previous knowledge regarding handling and operation is required and essential. Unpack the instrument carefully and examine the devices for any possibility of damage. In case of any damage or missing items contact your distributor immediately.

If the package present any damaged, the sterility of the device is not guarantee.

8. Handling and Operation

Both parts of the device are supplied with HF connectors that need to be fixed to an appropriate active cords.

It is recommended to use an echo-endoscope.

Cysto-Gastro-Sets require an endoscope's working channel of at least 3.7 mm Ø for 10Fr Set, 3.2 mm Ø for the 8.5 Fr Set or 2.8 mm for the 6 Fr.

The devices are compatible with ERBE ICC generators.

Please refer to the respective manufacturer's manual.

If HF from other manufacturers are used, check their compatibility and follow their instructions. Applicable limit values:

- Max. Power: 300W
- Max. HF voltage: 1600 Vpeak

9. Procedure

After the size and position of the cyst or pseudo-cyst are determined by endosonography, the device is inserted into the endoscope's operating channel over a guide wire (.025" or .035").

- For the 6 Fr model you can puncture the stomach or duodenal wall and the cyst directly with the catheter of the device.
- For the 8.5 and 10 Fr model the stomach or the duodenal wall and the cyst have to be previously puncture with an EUS Needle or a diathermic needle.
- If necessary, fluids can be aspirated from the cavity through the lateral Luer-Lock adapter at the proximal end of the device.
- Now connect the active cord to the lateral connector and to the HF generator.
- Switch on electricity and dilate the access previously created and cyst's opening by pushing forward the device's metal tip carefully.
- Remove the active cord and leave the device at its position.
- Finally, remove the device out of the working channel, while ensuring that the guide wire stays in position.
- Use the guide wire to place a Double Pigtail or Metallic stent in order to drain the cyst, after the cystostomy.

10. Precaution/ Warnings

The surgery room should contain non flammable substances.

-  The device when applied to a patient with a pacemaker implanted may cause malfunction or failure of the pacemaker, seriously affecting the patient. Before proceeding, always confirm with a cardiologist or the manufacturer of the pacemaker that it is safe.
-  After the device is connected with the equipment, the cutting wire cannot touch any other equipment. Spark discharge may cause clinician burns.
-  When using an electrocardiograph or other physiological monitoring equipment simultaneously with the device on a patient, any monitoring equipment should be placed as far away as possible from the electrodes used with the device. Needle monitoring electrodes should not be used, as they may cause patient burns.
-  Electromagnetic interferences can occur when the HF-generators are used in close proximity of diagnostic or monitoring equipment leading to inappropriate assessment of the intraoperative state of the patient.
-  When using the device in the vicinity of the heart, be sure to use it with the minimum necessary output. Spark discharge during operation may affect the heart.

11. Cleaning and Sterilization

The disposable version of this device is designed and warranted for single use!

12. Storage

This device is delivered sterilized and in order to maintain this state until its first use it must be kept in its original sterilized package.

-  Do not put any objects on the instrument or its package!
Do not store the instruments near aggressive chemical products!

-  Do not expose the instruments to direct or indirect sunlight or other ultra-violet rays!

-  Keep in dry area

Complaints will not be taken under consideration if the instruments have been stored improperly.

13. Risk in case of reuse

This devices is for single use only. G-Flex denied all responsibility in case of reuse. The reuse of single use device might present risk to the safety of the user or patient, due to possible uncontrolled contamination and / or unreliability of the device performance.

14. Serious Incident, Complaints and repairs

In case of a serious incident, please report to the manufacturer and to the competent authority of the European Member State where the user is established. Shall the incident occur outside the Europe, please report to the manufacturer and follow the local applicable regulation.

In health protection of our employees only disinfected or sterilised instruments will be accepted for analysis or repairing. In any case disinfection or sterilisation dates and validity must be labelled outside the package. If this requirement has not been fulfilled the instruments will be returned to the sender without analysis or repairing.

15. Declarations

We, hereby, declare that the medical device neither contain nor is produced with tissue of animal origin.

We, hereby, declare that the medical device neither contain nor is produced with tissue or extracted material of human origin.

We, hereby, declare that the medical device do not contain medicinal substance.

We, hereby, declare that the medical device do not contain human blood derivative.

We, hereby, declare that the medical device do not contain phthalates or other products that could leak from the device that may be carcinogenic, mutagenic or toxic.

We, hereby, declare that the medical device is not composed of substances or combination of substances intended to be introduced into the human body.

16. Product disposal

After use, these products may be a potential biohazard. Handle and dispose of in accordance with accepted medical practice and applicable local, state and federal laws and regulations.

17. Legal foundation

The Law of the European Union applied.

18. Product support

In case of questions or difficulties concerning our instruments please contact your local distributor or G-FLEX directly during regular working hours.

Times available: Monday to Friday; 9am to 5pm (MET)

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