EN - DISPOSABLE HEMOCLIP
FR - HÉMOCLIP À USAGE UNIQUE
NL - WEGWERP-HEMOCLIP
ES - APLICADOR DE CLIP DESECHABLE
DE - HÄMOKLEMME FÜR DEN
EINMALGEBRAUCH
PT - HEMOCLIP DESCARTÁVEL
IT - CLIP PER EMOSTASI MONOUSO
EL - ΑΙΜΟΣΤΑΤΙΚΌ ΚΛΙΠ Μ΄ΙΑΣ ΧΡ΄ΗΣΗΣ
AR - مشبك الأوعية أحادي الاستعمال
TR - TEK KULLANIMLIK HEMOKLIPS
SV - HEMOCLIP FÖR ENGÅNGSBRUK
CS – JEDNORÁZOVÁ HEMOSVORKA
DA - HÆMOKLIPS TIL ENGANGSBRUG
HU – EGYSZER HASZNÁLATOS HEMOCLIP
RO - CLEMĂ HEMOSTATICĂ DE UNICĂ
FOLOSINȚĂ



G-Flex Europe S.R.L. - 20, Rue de l'industrie 1400
Nivelles - Belgium | Phone: +32 67.88.36.65
Email: info@g-flex.com | Website: www.g-flex.com

EN - Catalog number	
FR - Numéro de catalogue	
NL - Catalogusnummer	
ES - Número de catálogo	
DE - Katalog-Nummer	
PT - Número do catalogo	
IT - Numero di catalogo	U2511-180
EL - Αριθμός καταλόγου	U2511-230
AR - العربية- رقم الكتالوج	U2516-180
TR - Katalog numarası	U2525-230
SV - Katalognummer	
CS - Katalogové číslo	
DA - Katalognummer	
HU – Katalógusszám	
RO - Număr de catalog	

EN - INSTRUCTION FOR USE	4
FR - MODE D'EMPLOI	9
NL - GEBRUIKSAANWIJZING	15
ES - INSTRUCCIONES DE USO	21
DE - GEBRAUCHSANWEISUNG.....	27
PT - INSTRUÇÕES DE USO.....	33
IT - ISTRUZIONI PER L'USO	39
EL - ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ	45
AR العربية - تعليمات الاستخدام.....	51
TR - KULLANIM TALIMATLARI	57
SV - BRUKSANVISNING	63
CS – NÁVOD K POUŽITÍ	69
DA - BRUGSANVISNING	74
HU – HASZNÁLATI ÚTMUTATÓ	80
RO - INTRUCȚIUNI DE UTILIZARE	86

U25(11-16) Series

EN - INSTRUCTION FOR USE

1. Intended Use

Disposable Hemoclip are intended to place clip within the Gastrointestinal (GI) tract for the following purposes:

- Endoscopic marking;
- Hemostasis:
- Anchor to attached jejunal feeding tubes to the wall of the small bowel;
- Anchor to self expandable metallic stents (SEMS) to prevent migration;
- As a supplementary method to close luminal perforations (< 10 mm) in the GI tract that can be treated conservatively

2. Intended Population

Our device is intended for only adult patient according to the physician recommendation taking into account the contraindications.

Pregnant woman are not contraindicated for the use of disposable hemoclip.

3. Contraindications

- Do not use this device when hemostasis cannot be verified visually with an endoscopic field of view.
- Arteries greater than 2 mm.
- Polyps greater than 20 mm in diameter.
- Mucosal / Sub- Mucosal defects up to 2.5 cm.

4. Possible Complications

Possible complications are Recurrent bleeding; Mechanical issues / Break issues; Migration issues; Perforations

5. User

The users of G-FLEX instruments must be specialists in their fields. An appropriate and specific training for preparation, care and maintenance of the flexible instruments is required.

6. Product Description

The U2525 series device is made of a pre-loaded, rotatable, radiopaque, single-use clip (See picture on page 3):

1. Clip jaws

2. Clip body

And a flexible delivery system made of:

3. Spring tube

4. Hole in the slider

5. Slider

6. Locking system in the

7. Handle

8. Thumb ring

The U25(11-16) series device is made of a pre-loaded, rotatable, radiopaque, single-use clip and a flexible delivery system made of (See picture on page 3):

1. Thumb ring
2. Handle
3. Slider
4. Rotation Wheel
5. Spring tube
6. Clip seat
7. Clip body
8. Clip jaws

The handle in the proximal part to manipulate the rotation and the deployment, a long flexible tube and the clip at the distal part. The handle is also used for opening and closing of the clip. It is designed to work in gastroscope and colonoscopes (see model working length) with a working channel equal to or greater than 2.8 mm diameter.

7. Advices before first application

Please read carefully and follow all safety operating instructions and warnings before first application of the device.

A previous knowledge regarding handling and operation is required and essential. Unpack the instrument carefully and examine the devices for any possibility of damage. In case of any damage or missing items contact your distributor immediately.

If the package present any damaged, the sterility of the device is not guarantee.

8. Handling and Operation

Ensure the device and its packaging are intact before use. Remove the device from it's packaging before use.

Please observe compatibility between the operating channel and device diameter.

Before inserting the device into the biopsy channel ensure that the clip is closed.

The Disposable Hemoclip can be opened and closed no more than 5 times prior the final deployment.


9. Procedure

Insertion into the Biopsy Channel

- Insert the device into the biopsy channel until it comes out of the distal part of the endoscope.
- Once under direct visualization, you are ready for the next step


Deployment of the clip

- The clip is deployed closed. To open the clip, move the slider distally (away from the thumb ring)(see picture below).


Rotation of the clip

- **U2525 Series** - To rotate the clip, in order to adjust the angle between the clip and the lesion, just rotate the complete handle in any direction while holding the blue sheet still. (See picture below)
- **U25(11-16) Series** - To rotate the clip, in order to adjust the angle between the clip and the lesion, just rotate the complete handle in any direction. (See picture below)


Closing the clip

- To close the clip jaws on the desired location move the slider proximally (close to the thumb ring) until the first tactile resistance is felt in the handle (see picture below).


- **U2525 Series** - To permanently deploy the clip, continue moving the slider proximally (toward the thumb ring) beyond the tactile resistance point until the Hole in the slider reaches the Locking System. At this point, click may be heard or felt. The Locking System ensures that the tracking wire stays inside the spring

tube (spiral) and does not damage the endoscope biopsy channel. (see picture below).

- **U25(11-16) Series** - To permanently deploy the clip, continue moving the slider proximally (toward the thumb ring) beyond the tactile resistance point. At this point, click may be heard or felt. (see picture below).

Removal of the Delivery System

- Remove the whole delivery system from the biopsy channel and check the result with the endoscope, if needed repeat the whole procedure with other clips.

10. Precautions/ Warnings

MR Conditional:

Non-clinical testing, a patient with this clip can be scanned safely immediately after placement under the following conditions:

- Static magnetic field of 3.0 Tesla or less
- Maximum spatial magnetic gradient field of 1600 gauss/cm or less
- Normal operating mode: Maximum MR System reported, whole-body- average specific absorption rate (SAR) of 2.0W/Kg for 15 minutes of scanning or less (i.e., per scanning sequence).

MRI- Related Heating:

In non clinical testing, the clip produced maximum temperature rise of no more than 5°C, during 15 minutes of MR imaging performed in a 3.0 Tesla MR System.

Image Artefacts:

MR image quality may be compromised if the area of interest is within approximately 30mm of the position of the clip as found during non-clinical testing, in a 3.0 Tesla MR system.

WARNING : There is a potential risk of rebleeding of small arteries due to torque place on the clip during the MRI.

11. Cleaning and Sterilization

This device is designed and warranted for single use!

12. Storage

This device is delivered sterilized and in order to maintain this state until its first use it must be kept in its original sterilized package.

Do not put any objects on the instrument or its package!
Do not store the instruments near aggressive chemical products!

Do not expose the instruments to direct or indirect sunlight or other ultra-violet rays!

Keep in dry area

Complaints will not be taken under consideration if the instruments have been stored improperly.

13. Risk in case of reuse

This device is for single use only. G-Flex denied all responsibility in case of reuse. The reuse of single use device might present risk to the safety of the user or patient, due to possible uncontrolled contamination and / or unreliability of the device performance.

14. Serious Incident, Complaints and repairs

In case of a serious incident, please report to the manufacturer and to the competent authority of the European Member State where the user is established. Shall the incident occur outside the Europe, please report to the manufacturer and follow the local applicable regulation.

In health protection of our employees only disinfected or sterilised instruments will be accepted for analysis or repairing. In any case disinfection or sterilisation dates and validity must be labelled outside the package. If this requirement has not been fulfilled the instruments will be returned to the sender without analysis or repairing.

15. Declarations

We, hereby, declare that the medical device neither contain nor is produced with tissue of animal origin.

We, hereby, declare that the medical device neither contain nor is produced with tissue or extracted material of human origin.

We, hereby, declare that the medical device do not contain medicinal substance.

We, hereby, declare that the medical device do not contain human blood derivative.

We, hereby, declare that the medical device do not contain phthalates or other products that could leak from the device that may be carcinogenic, mutagenic or toxic.

We, hereby, declare that the medical device is not composed of substances or combination of substances intended to be introduced into the human body.

16. Product disposal

After use, these products may be a potential biohazard. Handle and dispose of in accordance with accepted medical practice and applicable local, state and federal laws and regulations.

17. Legal foundation

The Law of the European Union applied.

18. Product support

In case of questions or difficulties concerning our instruments please contact your local distributor or G-FLEX directly during regular working hours.

Times available: Monday to Friday; 9am to 5pm (MET)

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