EN - DISPOSABLE POLYPECTOMY SNARES

FR - ANSE À POLYPECTOMIE JETABLE

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Current Name: IFU-Disposable Polypectomy Snare Creation Date: 2025-03-26 Version 22

This document is only controlled in electronic format by the QM.

EN - Catalog number	
FR - Numéro de catalogue	
NL - Catalogusnummer	GF4300 GF4615
ES - Número de catálogo	GF4301 GF4615R
DE - Katalog-Nummer	GF4406 GF461525CR
PT - Número do catalogo	GF4410 GF4625
IT - Numero di catalogo	GF4532R GF4625R
EL - Αριθμός καταλόγου	GF4564R GF4635
AR - العربية- رقم الكاتالوج	GF4575R GF4635R
TR - Katalog numarası	GF4606 GF4710R
SV - Katalognummer	GF4606R GF4815R
CS - Katalogové číslo	GF4610 GF481525CR
DA - Katalognummer	GF4610R GF4825R
HU – Katalógusszám	
RO - Număr de catalog	

EN - INSTRUCTION FOR USE	4
FR - MODE D'EMPLOI	8
NL - GEBRUIKSAANWIJZING	14
ES - INSTRUCCIONES DE USO	19
DE - GEBRAUCHSANWEISUNG.....	24
PT - INSTRUÇÕES DE USO.....	30
IT - ISTRUZIONI PER L'USO	35
EL - ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ	40
AR - العربية - تعليمات الاستخدام	46
TR - KULLANIM TALIMATLARI	51
SV - BRUKSANVISNING	56
CS – NÁVOD K POUŽITÍ	61
DA - BRUGSANVISNING.....	66
HU – HASZNÁLATI ÚTMUTATÓ	70
RO - INSTRUCTIUNI DE UTILIZARE	75

EN - INSTRUCTION FOR USE

1. Intended Use

Disposable Polypectomy Snares are intended to be used to remove and/or cauterize any kind of polyps in the gastrointestinal tract.

2. Intended Population

Our device is intended for both adult and paediatric patient according to the physician recommendation taking into account the contraindications.

3. Contraindications

- Missing accordance and failing cooperation of the patients
- Non-fasting patients,
- Intestine wall fragility e.g: High-florid inflammation of the colon, toxic megacolon
- Peritonitis, acute abdomen
- Intestine perforation
- Ileus
- Sepsis
- Severe cardio-pulmonary diseases and decompensation
- Non-controllable hemorrhagic diatheses (prothrombin time below 50% partial thromboplastin time PTT, prolonged by more than two-fold, thrombocytes below 5000
- Insufficient preparation (cleaning)
- Pacemakers that cannot be set to a fixed frequency due to parasympathetic rhythm,
- Hip endoprothesis or metal implants,
- Pregnancy
- Recently created gastrointestinal anastomosis.

4. Possible complications

Possible complications associated with current application include perforation and bleeding.

5. User

The users of G-FLEX instruments must be specialists in their fields. An appropriate and specific training for preparation, care and maintenance of the flexible instruments is required.

6. Product Description

Disposable Polypectomy Snares are made of a flexible tube with a loop in different sizes and shapes at the distal end and a handle to expose and retract the loop into the catheter. The handle has a HF connector to be attached to a high frequency (HF) generator and an active cord.

7. Advices before first application

Please read carefully and follow all safety operating instructions and warnings before first application of the device.

A previous knowledge regarding handling and operation is required and essential. Unpack the instrument carefully and examine the devices for any possibility of damage. In case of any damage or missing items contact your distributor immediately.

 If the package present any damaged, the sterility of the device is not guarantee.

8. Handling and Operation

Insertion into the Endoscope

The Disposable Polypectomy Snare must be passed through the endoscope operating channel with its snare withdrawn into the tube. Please observe compatibility between the Disposable Polypectomy Snare and the operating channel diameters. Minimum operating channel diameter of 2.8 mm.

HF active cord Connection

Connect a compatible active cord to electricity connector at the handle and to the HF-generator.

HF Generators

Disposable Polypectomy Snares can be connected to all licensed HF monopolar generators by using the active cord.

Always check compatibility between HF generator and active cord. In doubt contact the manufacturer of generator or your local distributor.

The family of products has been tested in combination with the HF generator model ERBE ICC and the impedance has been controlled using the HF generator protocol.

Any HF generators:

- able to deliver the Max voltage of 1500Vp,
- compliant to the IEC 60601-1 standard at the time of its acquisition and
- able to perform an impedance control of the external electrical circuit prior delivering the HF energy

can be used in combination with the device.

Please refer to the respective manufacturer's manual.

If HF from other manufacturers are used, check their compatibility and follow their instructions. Applicable limit values:

- Max. Power: 300W
- Max. HF voltage CUT mode: 1000 Vp
- Max. HF voltage COAG mode: 2500 Vp

9. Procedure

 Clear the area from mucus or other fluids around the polyp as it may decrease the current density reducing the cutting power.

- Place the distal tip of the tube right in front of the polyp to be cut.
- By pushing the handle spool expose the snare fully out of the tube and capture the polyp into the open snare.

- Pull back the spool until the snare closes on the polyp stalk and stays tighten around the polyp pedicle.
- Verify, with the HF generator and the attached HF-cord, that the electrical circuit has been correctly created.
- Apply HF-electric current on the snare (Please observe HF generator manufacturer's instructions!).
- Cut the polyp by pulling the snare completely into the tube.
- Stop the HF-electric current and remove carefully the Disposable Polypectomy Snare out of the endoscope.

 Do not touch other tissue areas with the snare while using HF-electric current (high hurt risk for the patient)!

10. Precautions / Warnings

 The instrument and active cord need to be checked for cracks and/or damages on its insulation.

 Polypectomy snares are composed by a catheter and should never be manipulated with its catheter winded as this may damage the device and make its usage impossible!

 Make sure to have on hand a grasping forceps in case the snare get stuck in the tissue and needs to be removed by other means.

The surgery room should contain non flammable substances.

 After the device is connected with the equipment, the cutting wire cannot touch any other equipment. Spark discharge may cause clinician burns.

 When using an electrocardiograph or other physiological monitoring equipment simultaneously with the device on a patient, any monitoring equipment should be placed as far away as possible from the electrodes used with the device. Needle monitoring electrodes should not be used, as they may cause patient burns.

 Disposable polypectomy snares successfully pass radiated disturbance electrical testing. When the device is operating with the HF generator delivering HF current, electromagnetic interferences can occur when the HF-generators are used in proximity of diagnostic or monitoring equipment leading to potential inappropriate assessment of the intraoperative state of the patient.

 When using the device in the vicinity of the heart, be sure to use it with the minimum necessary output. Spark discharge during operation may affect the heart.

 The device when applied to a patient with an implantable pacemaker or any active implantable device may cause malfunction or failure of the pacemaker or the active implantable device, seriously affecting the patient. Before proceeding, always confirm with a cardiologist or the manufacturer of the pacemaker or the active implantable device that it is safe.

 If the patient required the use of an external defibrillator during the procedure involving the Disposable Polypectomy Snares, it must be totally removed or disconnected from the HF surgical generator before delivering the electrical shock.

 Do not operate the device if the entire length of the surgical accessories is exceeding the maximum permissible length defined for the HF surgical generator to be used

 Before to operate the HF surgical generator, verify that all the electrical connections have been properly established and no visible metallic parts are accessible on the snare's handle avoiding contact with wet tissues.

 The snare is a surgery accessory intended to deliver high-frequency electrical power to the body tissue thanks to the embedded electrical wire. Due to the very low electrical resistance, the snare is not dissipating heat. It's a media to deliver the energy to the interface between the tip of the snare and the body tissues to be treated. The significant electrical resistant of the body tissues to be treated (in contact with the tip of the snare) induce heat to the tissues. During the short duration of power delivery by the HF generator, electromagnetic interferences can occur which could affect diagnostic or monitoring equipment leading to inappropriate assessment of the intraoperative state of the patient. Such potential interferences are inherent to the use of HF surgical generator in combination with surgery tool such as the snare.

 The snare must be used with HF cord, HF generator and grounding pad that shall be CE marked and compliant with the electrical safety standards (IEC 60601), to ensure patient safety.

11. Cleaning and Sterilization

 This device is designed and warranted for single use!

12. Storage

This device is delivered sterilized and in order to maintain this state until its first use it must be kept in its original sterilized package.

 Do not put any objects on the instrument or its package!

Do not store the instruments near aggressive chemical products!

 Do not expose the instruments to direct or indirect sunlight or other ultra-violet rays!

 Keep in dry area

Complaints will not be taken under consideration if the instruments have been stored improperly.

13. Risk in case of reuse

This devices is for single use only. G-Flex denied all responsibility in case of reuse. The reuse of single use device might present risk to the safety of the user or patient, due to possible uncontrolled contamination and / or unreliability of the device performance.

14. Serious Incident, Complaints and repairs

In case of a serious incident, please report to the manufacturer and to the competent authority of the European Member State where the user is established. Shall the incident occur outside the Europe, please report to the manufacturer and follow the local applicable regulation.

In health protection of our employees only disinfected or sterilised instruments will be accepted for analysis or repairing. In any case disinfection or sterilisation dates and validity must be labelled outside the package. If this requirement has not been fulfilled the instruments will be returned to the sender without analysis or repairing.

15. Declarations

We, hereby, declare that the medical device neither contain nor is produced with tissue of animal origin.

We, hereby, declare that the medical device neither contain nor is produced with tissue or extracted material of human origin.

We, hereby, declare that the medical device do not contain medicinal substance.

We, hereby, declare that the medical device do not contain human blood derivative.

We, hereby, declare that the medical device do not contain phthalates or other products that could leak from the device that may be carcinogenic, mutagenic or toxic.

We, hereby, declare that the medical device is not composed of substances or combination of substances intended to be introduced into the human body.

16. Product disposal

After use, these products may be a potential biohazard. Handle and dispose of in accordance with accepted medical practice and applicable local, state and federal laws and regulations.

17. Legal foundation

The Law of the European Union applied.

18. Product support

In case of questions or difficulties concerning our instruments please contact your local distributor or G-FLEX directly during regular working hours.

Times available: Monday to Friday; 9am to 5pm (MET)

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