EN - BILIARY & PANCREATIC STENT

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EN - INSTRUCTION FOR USE

1. Intended Use

The Biliary & Pancreatic Stent are used to drain biliary or pancreatic obstructions. Biliary & Pancreatic Stents enable effective drainage of biliary and pancreatic obstructions, improving clinical outcomes.

The claimed clinical performance are:

- Successful drainage of biliary ducts: Technical success by insertion of at least one stent in the biliary or pancreatic ducts (79-100%) and clinical success by the disappearance of symptoms or clinical improvement by passage of fluid (55.6-100%)
- Successful drainage of pancreatic ducts: Technical success by insertion of at least one stent in the biliary or pancreatic ducts (80-100%) and clinical success by the disappearance of symptoms or clinical improvement by passage of fluid (68.4-100%)

2. Intended Population

Our device is intended for only adult patient according to the physician recommendation taking into account the contraindications.

3. Contraindication

The contraindications include without being limited to poor general condition such as respiratory failure, risk of aspiration, general blood coagulation and wound healing disorders, asymptomatic stenosis or jaundice, endoscopically impassible stenosis, acute, non-biliary pancreatitis, absent patient consent. The device is not indicating for paediatric procedures

4. Possible Complications

Stent migration, bleeding and perforation due to the stent migration, abdominal pain, stent occlusion, infection and cholangitis due to the stent occlusion/ obstruction, jaundice, pancreatitis.

5. User

The users of G-FLEX instruments must be specialists in their fields. An appropriate and specific training for preparation, care and maintenance of the flexible instruments is required.

6. Product Description

Biliary & Pancreatic Stent are made of radiopaque PE or PU plastic with or without flaps, with 1 or 2 or without pigtails. In order to facilitate the mounting orientation of the stent, they all have always and only the distal end tapered. The 7 Fr stents are designed to go over a 0.035" guide wire, please test the compatibility before the procedure. The 8.5 Fr and 10 Fr stents are compatible with common application systems designed for endoscopic plastic stent placement if they respect the stent's diameter and compatibility should be tested prior to use.

G-Flex offers kits in 7 Fr, composed of a stent, a positioning sleeve and an outer pusher. We also have kits in 8.5 Fr and 10 Fr, composed by a stent **pre-mounted** on the guiding catheter with a positioning sleeve and an outer pusher. This guiding catheter has a standard luer-lock Y-connector for contrast injection through the guiding catheter. Both kits can be used with a 0.035" guide wire.

If the application system is not presented with your stent, please consider having one on hands always respecting the French size of the stent. (Refer to the current catalogue for more information)

The table below details the composition of the implantable part of the device:

Product Series	Raw material	Percentage
Bxxy series	Polyethylene	80%
	Bismuth subcarbonate	20%
Pxxx series	Polyurethane	80%
	Barium Sulfate	20%
EZYFLEXxx series	Polyethylene	76.15%
	Bismuth subcarbonate	20%
	FIBAPLAST BLAU PE	3.85%

7. Advices before first application

Please read carefully and follow all safety operating instructions and warnings before first application of the device.

A previous knowledge regarding handling and operation is required and essential.

Unpack the instrument carefully and examine the devices for any possibility of damage. In case of any damage or missing items contact your distributor immediately.

If the package present any damaged, the sterility of the device is not guarantee.

8. Handling and Operation

Our Biliary & Pancreatic Stent are packed with an implant card that should be completed by the healthcare institution/provider and hand over to the patient.

The card comes with one side pre-printed with the manufacturer and stent traceability information and the side that must be completed according to the example shown on page 4. A website with more information is available at the following address:

www.g-flex.com/implantcard

 The Biliary & Pancreatic Stent are placed in conjunction with an endoscope, a guide wire, a pusher and, when necessary, a guiding catheter

 To facilitate the application, an endoscopic sphincterotomy prior to guide wire placement is recommended. (Please observe HF-generator and sphincterotome manufacturer's instructions).

 It is recommended NOT to remove the stent by inflating a balloon inside the stent.

 The patient and the stent should be regularly followed-up to remove the stent immediately in case of clog, migration or fallout from the papilla.

9. Procedure

For stent pre-mounted over the application system in 8.5 Fr or 10 Fr, please proceed as follows:

- Choose the diameter and length of the stent to be placed.
- Remove the kit from its packaging and remove the yellow stopper at the distal extremity of the stent.
- Introduce slowly the complete device over the guide wire by inserting the guide wire inside the distal end of the guiding catheter until the tapered distal extremity stent reach the entrance of the biopsy channel.
- Push the positioning sleeve forward over the proximal flap to lower it in order to protect it while inserting the stent and the rest of the device inside the biopsy channel.

 Note the positioning sleeve will stay out of the biopsy channel.

- Advance the guiding catheter until the radiopaque ring passes the obstructed area.
- Once the guiding catheter is in the desired place, unscrew the luer-lock connection of the Y-connector to release the pusher from the guiding catheter.
- Slowly advance the pusher until it reaches the proximal extremity of the stent.

 To hold the guide wire in place with the guiding catheter, use the sealing cap on the Y-connector. Injection of contrast is possible via the guiding catheter Y-connector lateral connection. Do this with the sealing cap closed.

- To place the stent into its final position, hold the guiding catheter with one hand while slowly advancing the pusher.

 Keep guiding catheters position while placing the stent

- Once the stent is placed, pull the guiding catheter and the guide wire back into the biopsy channel while keeping the pusher in its position. After the complete removal of those devices, the stent will regain its original shape.
- Finally, remove pusher, guiding catheter and guide wire out of the scope.

For 7 Fr stents in kit, please proceed as follows:

- Choose the length of the stent to be placed
- Remove the kit (stent, positioning sleeve and application system) out of the packaging.
- Over the guide wire, place the tapered distal extremity of the stent, followed by the application system and the positioning sleeve over the application system until the stent reach the entrance of the biopsy channel.
- Push the positioning sleeve forward over the proximal flap to lower it in order to protect it while inserting the stent and the rest of the device inside the biopsy channel.

 Note the positioning sleeve will stay out of the biopsy channel.

- To place the stent in the desired location, hold the guide wire with one hand while pushing the application system forward with the other hand until the stent is in place.
- Once the stent is placed, pull the guide wire back into the biopsy channel while holding the application system in place. After the complete removal of the guide wire, the stent will regain its original shape.
- Finally, remove application system out of the scope

For stents in 5 Fr or 7 Fr, please proceed as follows:

- Choose the shape, diameter and length of the stent to be placed.

 Note: make sure to have on hands an application system with the same diameter as the stent selected (5 Fr or 7 Fr).

- Remove the stent with its positioning sleeve and the application system out of their packaging.
- Over the guide wire, place the tapered distal extremity of the stent, followed by the application system and the positioning sleeve over the application system until the stent reach the entrance of the biopsy channel.
- Push the positioning sleeve forward over the proximal flap to lower it in order to protect it while inserting the stent and the rest of the device inside the biopsy channel.

 Note the positioning sleeve will stay out of the biopsy channel.

- To place the stent in the desired location, hold the guide wire with one hand while pushing the application system forward with the other hand until the stent is in place.
- Once the stent is placed, pull the guide wire back into the biopsy channel while holding the application system in place. After the complete removal of the guide wire, the stent will regain its original shape.
- Finally, remove application system out of the scope.

For stents in 8.5 Fr or 10 Fr, please proceed as follows:

- Choose the shape, diameter and length of the stent to be placed.

 Note: make sure to have on hands an application system with the same diameter as the stent selected (8.5 Fr or 10 Fr).

- Remove the stent with its positioning sleeve and the application system out of their packaging.
- Insert the positioning sleeve from the distal end of the guiding catheter all the way up the proximal Y connector.
- Place the stent from its proximal extremity (NOT tapered side) over the guiding catheter until it reaches the distal extremity of the pusher.
- Insert the assembled stent and application system over the guide wire by inserting the guide wire inside the distal end of the guiding catheter until the tapered distal extremity stent reach the entrance of the biopsy channel.
- Push the positioning sleeve forward over the proximal flap to lower it in order to protect it while inserting the stent and the rest of the device inside the biopsy channel.

 Note the positioning sleeve will stay out of the biopsy channel.

- Advance the guiding catheter until the radiopaque ring passes the obstructed area.
- Once the guiding catheter is in the desired place, unscrew the luer-lock connection of the Y-connector to release the pusher from the guiding catheter.
- Slowly advance the pusher until it reaches the proximal extremity of the stent.

 To hold the guide wire in place with the guiding catheter, use the sealing cap on the Y-connector. Injection of contrast is possible via the guiding catheter Y-connector lateral connection. Do this with the sealing cap close.

- To place the stent into its final position, hold the guiding catheter with one hand while slowly advancing the pusher.

 Keep guiding catheters position while placing the stent

- Once the stent is placed, pull the guiding catheter and the guide wire back into the biopsy channel while keeping the pusher in its position. After the complete removal of those devices, the stent will regain its original shape.
- Finally, remove pusher, guiding catheter and guide wire out of the scope.

For stents removal, please proceed as follows:

- As the stent has been placed in the desired location with the proximal extremity appearing into the digestive system cavity, it can be removed by gently pulling it with a foreign body retriever forceps or a snare or even a basket.
- Once the stent has been pulled out remove it from the patient along with the scope and the removing device used.

10. Precaution/ Warnings

The stent is radiopaque, MRI compatibility has not been established.
In case of long-term use, the stent should be monitored for replacement every three months intervals as these stents are not intended for permanent implant use.
The device shall be used with licensed iodine-based or barium-based contrast agent.

11. Cleaning and Sterilization

 These products are designed and warranted for single use!

12. Storage

This device is delivered sterilized and in order to maintain this state until its first use it must be kept in its original sterilized package.

 Do not put any objects on the instrument or its package!
Do not store the instruments near aggressive chemical products!

 Do not expose the instruments to direct or indirect sunlight or other ultra-violet rays!

 Keep in dry area

Complaints will not be taken under consideration if the instruments have been stored improperly.

13. Risk in case of reuse

These devices are for single use only. G-Flex denied all responsibility in case of reuse. The reuse of single use device might present risk to the safety of the user or patient, due to possible uncontrolled contamination and / or unreliability of the device performance.

14. Serious Incident, Complaints and repairs

In case of a serious incident, please report to the manufacturer and to the competent authority of the European Member State where the user is established. Shall the incident occur outside the Europe, please report to the manufacturer and follow the local applicable regulation.

In health protection of our employees only disinfected or sterilised instruments will be accepted for analysis or repairing. In any case disinfection or sterilisation dates and validity must be labelled outside the package. If this requirement has not been fulfilled the instruments will be returned to the sender without analysis or repairing.

15. Declarations

We, hereby, declare that the medical device neither contain nor is produced with tissue of animal origin.

We, hereby, declare that the medical device neither contain nor is produced with tissue or extracted material of human origin.

We, hereby, declare that the medical device do not contain medicinal substance.

We, hereby, declare that the medical device do not contain human blood derivative.

We, hereby, declare that the medical device do not contain phthalates or other products that could leak from the device that may be carcinogenic, mutagenic or toxic.

We, hereby, declare that the medical device is not composed of substances or combination of substances intended to be introduced into the human body.

16. Product disposal

After use, these products may be a potential biohazard. Handle and dispose of in accordance with accepted medical practice and applicable local, state and federal laws and regulations.

17. Legal foundation

The Law of the European Union applied.

18. Product support

In case of questions or difficulties concerning our instruments please contact your local distributor or G-FLEX directly during regular working hours.

Times available: Monday to Friday; 9am to 5pm (MET)

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