- EN - MULTIBAND LIGATOR (EZYSHOOT®)**  
**FR - LIGATEUR MULTIBANDES (EZYSHOOT®)**  
**NL - MULTIBAND-LIGATOR (EZYSHOOT®)**  
**ES - LIGADOR DE MULTIBANDAS (EZYSHOOT®)**  
**DE - MULTIBAND-LIGATUR-SET (EZYSHOOT®)**  
**PT - KIT DE LIGADURA (EZYSHOOT®)**  
**IT - DISPOSITIVO PER LEGATURA MULTIPLA (EZYSHOOT®)**  
**EL - ΣΥΣΚΕΥΉ ΑΠΟΛΊΝΩΣΗΣ ΠΟΛΛΑΠΛΊΩΝ ΔΑΚΤΥΛΊΩΝ (EZYSHOOT®)**  
**AR - رابط متعدد الأربطة (EZYSHOOT®)**  
**TR - ÇOK BANTLI LIGATÖR (EZYSHOOT®)**  
**SV - MULTIBAND LIGATOR (EZYSHOOT®)**  
**CS – VÍCEGUMIČKOVÝ LIGÁTOR (EZYSHOOT®)**  
**DA - MULTIBÅND-LIGATOR (EZYSHOOT®)**  
**HU – TÖBBGYŰRŰS LIGÁTOR (EZYSHOOT®)**  
**RO - LIGATOR CU BENZI MULTIPLE (EZYSHOOT®)**  
**BG – МНОГОЛЕНТОВ ЛИГАТОР (EZYSHOOT®)**  
**SK – VIACGUMIČKOVÝ LIGÁTOR (EZYSHOOT®)**  
**LT – DAUGIAJUOSTIS LIGATORIUS (EZYSHOOT®)**  
**NO - MULTIBÅNDS LIGATOR (EZYSHOOT®)**




G-Flex Europe S.R.L. - 20, Rue de l'industrie 1400  
Nivelles - Belgium | Phone: +32 67.88.36.65  
Email: [info@g-flex.com](mailto:info@g-flex.com) | Website: [www.g-flex.com](http://www.g-flex.com)

EN - Versions	Standard	
FR - Versions	Standard	
NL - Versies	Standaard	
ES - Versiones	Estándar	
DE - Versionen	Standard	
PT - Versões	Padrão	
IT - Versioni	Standard	
EL - Εκδόσεις	Κανονική	
AR - الإصدارات	قياسي	
TR - Modeller	Standart	
SV - Versioner	Standard	
CS – Verze	Standard	
DA - Versioner	Standard	
HU - Versions	Szabványos	
RO - Versiuni	Standard	
SK – Verzie	Štandardný	
BG – Версии	Стандартен	
LT - Versijos	Standartas	
NO – Versjoner	Standard	
		GF-OVL100
		GF-OVL300
		GF-OVL501
		GF-OVL510

EN - INSTRUCTION FOR USE .....	6
FR - MODE D'EMPLOI.....	10
NL - GEBRUIKSAANWIJZING.....	15
ES - INSTRUCCIONES DE USO .....	19
DE - GEBRAUCHSANWEISUNG.....	24
PT - INSTRUÇÕES DE USO .....	28
IT - ISTRUZIONI PER L'USO .....	33
EL - ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ.....	37
AR - العربية - تعليمات الاستخدام.....	43
TR - KULLANIM TALIMATLARI.....	47
SV - BRUKSANVISNING .....	52
CS – NÁVOD K POUŽITÍ .....	56
DA - BRUGSANVISNING .....	60
HU – HASZNÁLATI ÚTMUTATÓ.....	64
RO - INSTRUCȚIUNI DE UTILIZARE.....	69
BG – ИНСТРУКЦИЯ ЗА УПОТРЕБА .....	73
SK – NÁVOD NA POUŽITIE.....	79
LT – NAUDOJIMO INSTRUKCIJA .....	83
NO – BRUKSANVISNING.....	87

# EN - INSTRUCTION FOR USE

## 1. Intended Use

The Multiband Ligator is used via endoscope to ligate oesophageal varices at and above the oesophageal junction. The G-Flex Multiband Ligator generates a haemostasis effect by providing mechanical compression in the clamped tissues and vessels.

### Performances

The claimed clinical performances are primary and secondary prophylaxis of oesophageal variceal bleeding: eradication of the varices (60-80%). Several sessions (2-6) may be required until the total eradication of the varices.

Acute oesophageal variceal bleeding management: stop bleeding 80-100%

## 2. Intended Population

Our device is intended for both adult and paediatric patient according to the physician recommendation taking into account the contraindications.

## 3. Contraindication

The contraindications include without being limited to cricopharyngeal or oesophageal narrowing/stricture, tortuous oesophageal diverticula, known or suspected oesophageal perforation, asymptomatic rings or webs, coagulopathy, active infection in the esophageal area such as esophagitis, distal esophageal tumors, pregnancy. Use of ligation bands is contraindicated in patients with a known hypersensitivity to latex and for any other condition which would otherwise contraindicate gastrointestinal endoscopy.

## 4. Possible Complications

Possible complications associated with the gastrointestinal endoscopy include and not limited to perforations, fever, aspirations, haemorrhage, infections, hypo-tension, allergic reactions to medication, cardiac arrhythmia or arrest, respiratory depression or arrest. Furthermore, complications that may occur with oesophageal banding include and not limited to band dislodgment, chest tightness, retrosternal pain, dysphagia or odynophagia, and post-EVL ulcer with or without bleeding, increase of portal hypertensive gastropathy, bacteremia, nausea, laryngeal, retrosternal pain, laceration, stricture formation, oesophageal perforation, obstructions, variceal recurrence and rebleeding.

## 5. User

The users of G-FLEX instruments must be specialists in their fields. An appropriate and specific training for preparation, care and maintenance of the flexible instruments is required.

## 6. Product Description

All versions are composed of the following parts (see pictures on page 4):

1. Big wheel
2. **Multiband Ligator Handle**
3. Blue Luer-lock cap
4. Silicone sealing ring
5. **Multiband Ligator tip**
6. Alignment arrows
7. Small wheel
8. Tracking wire
9. Plastic part of the introduction catheter

## 10. Introduction catheter

11. Tracking wire inside the introduction catheter
12. Latex or Latex-Free ring
13. Penultimate ring in another color.
14. **Band Barrel (Cap)**
15. Soft part of the cap
16. Loop of the band barrel

## 7. Advices before first application

Please read carefully and follow all safety operating instructions and warnings before first application of the device.

A previous knowledge regarding handling and operation is required and essential.

Unpack the instrument carefully and examine the devices for any possibility of damage. In case of any damage or missing items contact your distributor immediately.

## 8. Handling and Operation

This device is delivered non sterile and must **NOT** be sterilised before use!

This device **CANNOT** be used with an over tube.

### Mounting the Multiband Ligator:

1. Unpack and check all parts of the system carefully. Ensure that all parts are in the package.
2. The tracking wire of the Multiband Ligator is already pre-mounted into an introduction catheter.
3. Straighten the endoscope. We recommend the use of a biopsy valve. Keep its lid open.
4. Insert the catheter into the biopsy channel until it comes out of the distal tip of the scope.
5. Fix the handle inside the biopsy channel port.

Hold the **distal tip** of the introduction catheter and **remove the plastic part by slowly** pulling it out of the endoscope until only tracking wire is left for the next fixing step.

**Do not roll up the catheter on your hands to remove as it will block the inner tracking wire and then prevent the catheter to go out from the biopsy channel**

6. Take the band barrel loop and connect it with the tracking wire loop properly as shown on page 5, making a flat knot.

In order to tighten the knot, pull only the wires instead of pulling the whole cap as this could misfire a ring before the procedure.

7. Ensure that the handle has the alignment arrows facing each other.
8. Using the small wheel, **carefully and slowly** pull the band barrel close to the tip of the scope.
9. Mount the band barrel on the tip of the scope and ensure that the soft part is properly fixed on the scope in a way that it cannot fall off.
10. Rotate the band barrel in order to prevent the tracking wire being in front of the endoscope camera.
11. Now bend the endoscope in a J shape as to reproduce the position of introduction of the endoscope into the patient.
12. Slowly tension the tracking wire by turning the small wheel until you feel a small tension.
13. Verify the tension is enough by either looking at the monitor if the endoscope is already connected to the tower and lit up, or looking directly into the cap through the distal tip of the endoscope if the same is not connected to the tower. The wires must be in a straight position, tensioned.

 **Do not tension beyond** this slight tension point to avoid shooting a ring.

 If after inserting the handle of the device into the working channel entry port the handle is not stable enough or can move with soft movements of the scope, please remove the biopsy valve along with the handle of the device by 2cm up from the endoscope. Then press the handle further into the biopsy valve and place both handle and biopsy valve together into the endoscope working channel entry port again.

## 9. Procedure

### Ligation of oesophageal varices

1. Lubricate the scope and exterior part of the band barrel ensuring that no lubricant has entered the barrel.
2. Roll back the big wheel 180° to permit that the tracking wires is not fully straightened so that no band is shot while inserting the endoscope into the patient.
3. When in place, tension again the tracking wire using the small wheel to prepare for the shooting (be sure that the alignment arrows on the handle are facing each other).
4. Locate the selected varix and aspirate it into the band barrel.
5. To fire the band, slowly rotate 180° forward the big wheel. More than one ligation band for each varix might be needed to control acute bleeding.
6. Release the suction button of the scope, insufflate air and then withdraw the scope slightly to release the ligated varice.
7. Repeat the procedure until the last varix is banded.

 The luer-lock connection is available to ensure irrigation of the scope biopsy channel if required. After removing the blue cap, attach the luer-lock connector to a syringe filled with sterile water and irrigate.

 After irrigation put back the blue cap to ensure good suction.

 In case additional bands are needed, remove the scope and attach a new Multiband Ligator starting back from the mounting steps.

### Removing the Multiband Ligator

1. Once the ligation procedure is completed, remove the scope from the patient.
2. Dismantle the the Multiband Ligator as follows:
  - a. If any unused bands are still on the barrel shoot all remaining bands
  - b. Remove the Multiband Ligator handle from the biopsy channel by pulling it, the remaining tracking wire attached with the handle should be removed also.
  - c. Remove the band barrel from the scope tip.
  - d. Dispose all parts per institutional guidelines for biohazardous medical waste.

 If all bands have been used, you can remove the Multiband Ligator from the biopsy channel and continue using the scope and its biopsy channel while still inside the patient.

## 10. Precautions/Warnings

It is vital the biopsy channel size of the endoscope is co-ordinated with compatible devices in order to obtain the best results during the procedure. The minimum biopsy channel diameter required is 2.8 mm and length of up to 150 cm. The outer diameter of the scope to be used must be from 8,5 mm up to 11,5 mm. Band ligation may not be effective when applied to small varices.

It is vital to ensure the understanding of the technical principles, clinical applications and associated risks with GI endoscopy as well as endoscopic banding before using this product. This product must only be used by or under direct supervision of physician thoroughly trained in therapeutic endoscopy as well as vessel banding.

Oesophageal ligators are not intended for ligation of varices below the gastro-oesophageal junction.

Banding should begin at the gastro-oesophageal junction and proceed up the oesophagus in an helicoidal manner, respecting a minimum of 2cm in between 2 ligatures. Passing the endoscope over an already ligated varices may dislodge the band from its position.

Routine endoscopic examination by the physician is recommended to confirm the diagnosis requiring treatment of oesophageal varices prior to assembling the Multiband Ligator.

## 11. Cleaning and Sterilization

This device is designed and warranted for single use!

## 12. Storage

Do not put any objects on the instrument or its package!

Do not store the instruments near aggressive chemical products!

Do not expose the instruments to direct or indirect sunlight or other ultra-violet rays!

Keep in dry area.

Complaints will not be taken under consideration if the instruments have been stored improperly.

## 13. Risk in case of reuse

This device is for single use only. G-Flex denied all responsibility in case of reuse. The reuse of single use device might present risk to the safety of the user or patient, due to possible uncontrolled contamination and / or unreliability of the device performance.

## 14. Serious Incident, Complaints and repairs

In case of a serious incident, please report to the manufacturer and to the competent authority of the European Member State where the user is established. Shall the incident occur outside the Europe, please report to the manufacturer and follow the local applicable regulation.

In health protection of our employees only disinfected or sterilised instruments will be accepted for analysis or repairing. In any case disinfection or sterilisation dates and validity must be labelled outside the package. If this requirement has not been fulfilled the instruments will be returned to the sender without analysis or repairing.

## 15. Declarations

We, hereby, declare that the medical device neither contain nor is produced with tissue of animal origin.

We, hereby, declare that the medical device neither contain nor is produced with tissue or extracted material of human origin.

We, hereby, declare that the medical device do not contain medicinal substance.

We, hereby, declare that the medical device do not contain human blood derivative.

We, hereby, declare that the medical device do not contain phthalates or other products that could leak from the device that may be carcinogenic, mutagenic or toxic.

We, hereby, declare that the medical device is not composed of substances or combination of substances intended to be introduced into the human body.

## 16. Product disposal

After use, these products may be a potential biohazard. Handle and dispose of in accordance with accepted medical practice and applicable local, state and federal laws and regulations.

## 17. Legal foundation

The Law of the European Union applied.

## 18. Product support

In case of questions or difficulties concerning our instruments please contact your local distributor or G-FLEX directly during regular working hours.

**Times available:** Monday to Friday, 9am to 5pm (MET)

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