**EN - SCLEROTHERAPY & ENDOSCOPIC NEEDLES**  
**FR - SCLÉROTHÉRAPIE & AIGUILLE ENDOSCOPIQUE**  
**NL - SCLEROTHERAPIE- &**  
**ENDOSCOPISCHE NAALDEN**  
**ES - AGUJA ENDOSCOPICA & DE ESCLEROTERAPIA**  
**DE - SKLEROThERAPIE- UND ENDOSKOPIENADELN**  
**PT - AGULHAS DE ESCLEROTERAPIA E**  
**ENDOSCOPIA**  
**IT - AGHI PER SCLEROTERAPIA ED ENDOSCOPICI**  
**EL - ΒΕΛ'ΟΝΕΣ ΣΚΛΗΡΟΘΕΡΑΠΕΙΑΣ ΚΑΙ**  
**ΕΝΔΟΣΚΟΠΗΣΗΣ**  
**أبر المعالجة بالتصليب والتنظير الداخلي - AR**  
**TR - SKLEROTERAPİ İGNELERİ**  
**SV - SKLEROTERAPI VE ENDOSKOPI İGNELERİ**  
**CS – SKLEROTERAPEUTICKÉ A ENDOSKOPICKÉ**  
**JEHLY**  
**DA - SKLEROTERAPI- OG ENDOSKOPINÅLE**  
**HU – SZKLEROTERÁPIÁS ÉS ENDOSZKÓPOS TŰK**  
**RO - ACE PENTRU SCLEROTERAPIE ȘI**  
**ENDOSCOPICE**

G-Flex Europe S.R.L. - 20, Rue de l'industrie 1400  
Nivelles - Belgium | Phone: +32 67.88.36.65  
Email: [info@g-flex.com](mailto:info@g-flex.com) | Website: [www.g-flex.com](http://www.g-flex.com)

EN - Catalog number	
FR - Numéro de catalogue	
NL - Catalogusnummer	
ES - Número de catálogo	
DE - Katalog-Nummer	
PT - Número do catalogo	GF252202
IT - Numero di catalogo	GF252203
EL - Αριθμός καταλόγου	GF252203-25
العربية - رقم الكatalog	GF262302
TR - Katalog numarası	
SV - Katalognummer	
CS - Katalogové číslo	
DA - Katalognummer	
HU - Katalógusszám	
RO - Număr de catalog	

EN - INSTRUCTION FOR USE .....	3
FR - MODE D'EMPLOI .....	6
NL - GEBRUIKSAANWIJZING .....	10
ES - INSTRUCCIONES DE USO .....	13
DE - GEBRAUCHSANWEISUNG.....	17
PT - INSTRUÇÕES DE Uso.....	21
IT - ISTRUZIONI PER L'USO.....	24
EL - ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ .....	28
AR - العربية - تعليمات الاستخدام.....	32
TR - KULLANIM TALIMATLARI .....	36
SV - BRUKSANVISNING .....	39
CS – NÁVOD K POUŽITÍ .....	42
DA - BRUGSANVISNING.....	46
HU – HASZNÁLATI ÚTMUTATÓ .....	49
RO - INSTRUCTIUNI DE UTILIZARE .....	53

# **EN - INSTRUCTION FOR USE**

## **1. Intended Use**

Sclerotherapy & Endoscopic Needles are intended to introduce a sclerosing agent (any therapeutic vasoconstrictor agents, Histoacryl®, Cyanoacrylate, epinephrine) into selected sites to control actual or potential bleeding lesions in the digestive system or to inject liquid to aid in polypectomy procedures.

 Off-label use of the device or sclerosing agents is not promoted by G-flex.  
Healthcare professional shall ensure that the medicinal substance used is licensed in the territory.

## **2. Intended Population**

Our device is intended for both adult and paediatric patient according to the physician recommendation taking into account the contraindications

## **3. Contraindication**

- Patients with lesions inappropriate for injection therapy

## **4. Possible Complications**

- Fever, allergic reaction to medication, perforation, sepsis
- Respiratory depression or respiratory arrest, cardiac arrhythmias, or halt

## **5. User**

The users of G-FLEX instruments must be specialists in their fields. An appropriate and specific training for preparation, care and maintenance of the flexible instruments is required.

## **6. Product Description**

Sclerotherapy & Endoscopic Needles are made of a flexible catheter with an inner tube connected to the needle at the distal end and a handle with a piston and a luer-lock connection at the proximal end.

## **7. Advices before first application**

Please read carefully and follow all safety operating instructions and warnings before first application of the device.

A previous knowledge regarding handling and operation is required and essential. Unpack the instrument carefully and examine the devices for any possibility of damage. In case of any damage or missing items contact your distributor immediately.

 If the package present any damaged, the sterility of the device is not guarantee.

## **8. Handling and Operation**

Check to ensure that the needle is operating properly by pushing the needle piston in and out while holding the handle or catheter. The handle is equipped with two

stops. You should feel a “click” when the needle is in the fully deployed position and when the needle is in the fully retracted position.

## 9. Procedure

### For Sclerotherapy and Endoscopic Needles, please proceed as follows

- Prepare a Luer-Lock adaptable syringe with the agent to be injected.
- Make sure that the needle is in its fully exposed position. Then attach the syringe to the needle via the Luer-Lock connection at the top of the handle and flush out the air in the piston using the liquid agent. Visually verify that the liquid agent reaches the distal tip of needle and all air has been evacuated from the inner catheter.
- Now retract the needle into its catheter to full withdrawn position. (The inner piston part of the handle is retracted to the position UP)
- Then introduce the catheter into the working channel of an appropriately sized endoscope. Advance in small increments until the catheter emerges from the scope.
- Fully deploy the needle until you feel the piston “click” into place. (The inner piston part of the handle is fully pushed into the outer part of the handle - Position DOWN)
- Proceed with the injection at an oblique angle to the tissue site
- Following the injection(s), retract the needle until the needle piston “clicks” into the fully retracted position. (Position UP)

 When using the irrigation flush port, the needle must be fully retracted until the piston “click” into the retracted position. Prepare a luer-lock adaptable syringe with the irrigation agent to be injected and attach the syringe to the irrigation flush port via the luer-lock protruding from the side of the lateral connection port.

 When inserting or withdrawing device from scope, ensure that the distal end of the scope is not in the retroflexed position.

- Upon completion of the procedure, remove the instrument from the scope.

 When the device is being advanced or withdrawn from the endoscope, it is important to ensure that the needle is in its fully retracted position to avoid damaging the working channel of the endoscope. (The inner piston part of the handle is retracted to the position UP)

## 10. Precautions/ Warnings

 Sclerotherapy & endoscopic needles are composed by a catheter and should never be manipulated with its catheter winded as this may damage the device and make its usage impossible!

 For ethanolamine and cyanoacrylate only use a 19GA needle (GF262303)

## **11. Cleaning and Sterilization**

This device is designed and warranted for single use!

## **12. Storage**

This device is delivered sterilized and in order to maintain this state until its first use it must be kept in its original sterilized package.

Do not put any objects on the instrument or its package!

Do not store the instruments near aggressive chemical products!

Do not expose the instruments to direct or indirect sunlight or other ultra-violet rays!

Keep in dry area

Complaints will not be taken under consideration if the instruments have been stored improperly.

## **13. Risk in case of reuse**

This devices is for single use only. G-Flex denied all responsibility in case of reuse. The reuse of single use device might present risk to the safety of the user or patient, due to possible uncontrolled contamination and / or unreliability of the device performance.

## **14. Serious Incident, Complaints and repairs**

In case of a serious incident, please report to the manufacturer and to the competent authority of the European Member State where the user is established. Shall the incident occur outside the Europe, please report to the manufacturer and follow the local applicable regulation.

In health protection of our employees only disinfected or sterilised instruments will be accepted for analysis or repairing. In any case disinfection or sterilisation dates and validity must be labelled outside the package. If this requirement has not been fulfilled the instruments will be returned to the sender without analysis or repairing.

## **15. Declarations**

We, hereby, declare that the medical device neither contain nor is produced with tissue of animal origin.

We, hereby, declare that the medical device neither contain nor is produced with tissue or extracted material of human origin.

We, hereby, declare that the medical device do not contain medicinal substance.

We, hereby, declare that the medical device do not contain human blood derivative.

We, hereby, declare that the medical device do not contain phthalates or other products that could leak from the device that may be carcinogenic, mutagenic or toxic.

We, hereby, declare that the medical device is not composed of substances or combination of substances intended to be introduced into the human body.

## 16. Product disposal

After use, these products may be a potential biohazard. Handle and dispose of in accordance with accepted medical practice and applicable local, state and federal laws and regulations.

## 17. Legal foundation

The Law of the European Union applied.

## 18. Product support

In case of questions or difficulties concerning our instruments please contact your local distributor or G-FLEX directly during regular working hours.

**Times available:**      Monday to Friday;      9am to 5pm (MET)

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