EN - SPHINCTEROTOME/PAPILLOTOME

FR - SPHINCTÉROTOME/PAPILLOTOME

NL - SFINCTEROTOOM/PAPILLOTOOM

ES - ESFINTERÓTOMO/PAPILÓTOMO

DE - SPHINKTEROTOM/PAPILLOTOM

PT - ESFINCTERÓTOMO/PAPILÓTOMO

IT - SFINTEROTOMA/PAPILLOTOMA

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EN - INSTRUCTION FOR USE

1. Intended Use

Sphincterotomes/Papillotomes are intended for performing sphincterotomy (papilotomy) through the use of high frequency electrical current.

The Sphincterotomes are indicated for use in :

- The endoscopic sphincterotomy of the Papilla of Vater and/or the Sphincter of Oddi.
- The selective cannulation of the Common Bile Ducts (Selective cannulation is a prerequisite for successful pancreaticobiliary endoscopic interventions such as stone removal or stent placement)

The Sphincterotomes can also be used to inject contrast medium.

The Needle-Knife models K18NTL and GW08-000HF are intended for precut sphincterotomy in case of difficult cannulation.

The claimed performance are:

- Incision rate >99.5%
- Successful canulation / catheterization >80%
- Complete stone removal >90%

2. Intended Population

Our device is intended for adult patient according to the physician recommendation taking into account the contraindications.

3. Contraindications

- Patients with severe cardiopulmonary or renal dysfunction.
- Patients with acute pancreatitis or acute exacerbation of chronic pancreatitis.
- Patients with severe biliary tract infection.
- Patients who are allergic to iodine contrast agents.
- Patients who have upper gastrointestinal stenosis or obstruction.
- Contraindications for this device are those specific to endoscopic retrograde cholangiopancreatography(ERCP) and endoscopic sphincterotomy(ES).
- Hemodynamic instability in critically ill patients;

4. Possible Complications

Possible complications associated with current application include perforation and bleeding. Permanent and/or serious side effects of Sphincterotome Papillotome are related to the most frequent complications of endoscopic retrograde cholangiopancreatography (ERCP) and endoscopic biliary sphincterotomy, which include: acute pancreatitis, post-sphincterotomy bleeding/hemorrhage, sepsis of biliary origin (cholangitis and cholecystitis) and duodenal perforation.

5. User

The users of G-FLEX instruments must be specialists in their fields. An appropriate and specific training for preparation, care and maintenance of the flexible instruments is required.

6. Product Description

Sphincterotomes or also called Papillotomes are made of a teflon flexible tube with a cutting wire at the distal end in different shapes and sizes. At the proximal end, there is a handle with an HF connector. Most models have 3 lumen, one that can be used to inject contrast through a lateral port and the other that can be used to insert a CE certified guide wire. Please refer to the catalogue for more information

7. Advices before first application

Please read carefully and follow all safety operating instructions and warnings before first application of the device.

A previous knowledge regarding handling and operation is required and essential.

Unpack the instrument carefully and examine the devices for any possibility of damage. In case of any damage or missing items contact your distributor immediately.

If the package present any damaged, the sterility of the device is not guarantee.

8. Handling and Operation

Insertion into the Endoscope

The device must be passed through the endoscope operating channel with its cutting wire free from tension. Please ensure compatibility between the operating channel and device diameter. Any CE certified endoscope brand with a minimum operating channel diameter of 2.8 mm is compatible with the Sphincterotomes/ Papillotomes.

For the Knife-Snake® (GW08-000HF), any lumen that accepts a 0.035" guidewire is compatible.

HF active cord Connection

Connect a compatible active cord to electricity connector at the handle and to the HF-generator.

HF Generators

Sphincterotomes/Papillotomes can be connected to all licensed HF monopolar generators by using the active cord.

Always check compatibility between HF generator and active cord. In doubt contact the manufacturer of generator or your local distributor.

The family of products has been tested in combination with the HF generator model ERBE ICC and the impedance has been controlled using the HF generator protocol.

Any HF generators:

- able to deliver the Max voltage of 1500Vp,
- compliant to the IEC 60601-1 standard at the time of its acquisition and
- able to perform an impedance control of the external electrical circuit prior delivering the HF energy

can be used in combination with the device.

Please refer to the respective manufacturer's manual.

If HF from other manufacturers are used, check their compatibility and follow their instructions. Applicable limit values:

- Max. Power: 300W
- Max. HF voltage CUT mode: 1000 Vp

9. Procedure

- Connect the active cord to HF connector at the handle and to the HF-generator. (refer to the HF generator manufacturer instruction for proper connection)
- After connecting the active cord to the HF-generator the sphincterotomy (papillotomy) could be initiated.
- Verify, with the HF generator and the attached HF-cord, that the electrical circuit has been correctly created.
- Apply HF-electric current on the papillotomes (Please observe HF generator manufacturer's instructions!).

Do not touch other tissue areas with the cutting wire when using HF electric current (high risk of hurt for the patient).

Do not inject contrast when using HF electric current.

For devices that accept guide wire, it is strongly recommended to pull the guide wire back into the tube while using HF current in order to avoid a possible induction of voltage on the guide wire

Procedure for Knife Snake®

Whenever a precut papillotome is needed after a regular (lateral cut) wire papillotome is already in place, instead of exchanging the current papillotome for a needle-knife papillotome, an alternative solution can be the Knife-Snake® precut papillotome. For its use, please follow the instructions below:

- Remove completely the guidewire from the guidewire lumen of your double/triple lumen papillotome.
- With the papillotome and duodenoscope elevator untensioned, through the guide wire lumen, insert the Knife Snake® until its end can be seen through endoscopic view.
- Attach the yellow connector of the device to the Luer-Lock connector from the guide wire lumen entry of the papillotome and ensure it is well fixed.
- Now bend the papillotome to the desired precut position.
- Adjust the length of the Knife Snake® to the desired length of cutting wire and then fix the length at the yellow connector by pushing it distally.
- Connect the active cord to HF connector at the handle and to the HF-generator. (refer to the HF generator manufacturer instruction for proper connection)
- After connecting the active cord to the HF-generator the sphincterotomy (papillotomy) could be initiated.

- Once the desired cut opening has been reached, disconnect the HF active cord from the device and remove it from the papillotome replacing it again by a guide wire.

 Ensure that your papillotome guide wire lumen accepts a 0,035" guide wire

 Do not bend the papillotome or endoscope while inserting the device as this may block the device.

10. Precautions/ Warnings

 Remove the stiffening wire before use.

 After the device is connected with the equipment, the cutting wire cannot touch any other equipment. Spark discharge may cause clinician burns.

 When using an electrocardiograph or other physiological monitoring equipment simultaneously with the device on a patient, any monitoring equipment should be placed as far away as possible from the electrodes used with the device. Needle monitoring electrodes should not be used, as they may cause patient burns.

 Sphincterotome Papillotome successfully pass radiated disturbance electrical testing. When the device is operating with the HF generator delivering HF current, electromagnetic interferences can occur when the HF-generators are used in proximity of diagnostic or monitoring equipment leading to potential inappropriate assessment of the intraoperative state of the patient

 Do not bend the tube (< 15mm). It may cause tube to be deformed permanently, thus, the device cannot be used.

 When using the device in the vicinity of the heart, be sure to use it with the minimum necessary output. Spark discharge during operation may affect the heart.

 During the preparation, do not quickly push and pull the thumb ring causing the distal end of the tube to be bent less than 90°. Do not do anything that might cause tube to kink. Rotating the handle while inserting the device into endoscope channel would cause tube to kink.

 Sphincterotomes/Papillotomes are composed by a catheter and for that should always be straighten when used to avoid device damage.

 **For the K18NTL model, the handle need to be operated forward. Operating it backward may damage the device.**

 **The operating room must have good ventilation to avoid the accumulation of potentially flammable gases, as the use of electric current in the presence of these gases can cause an explosion.**

 The Sphincterotomes/Papillotomes must be used with HF cord, HF generator and grounding pad shall be CE marked and compliant with the electrical safety standards (IEC 60601), to ensure patient safety.

 The contrast medium used with the device must be licensed in the EU. Off-label use is not promoted and may cause serious harm to the patient.

 The device when applied to a patient with an implantable pacemaker or any active implantable device may cause malfunction or failure of the pacemaker or the active implantable device, seriously affecting the patient. Before proceeding, always confirm with a cardiologist or the manufacturer of the pacemaker or the active implantable device that it is safe.

 If the patient required the use of an external defibrillator during the procedure involving the papillotome, it must be totally removed or disconnected from the HF surgical generator before delivering the electrical shock.

 Do not operate the device if the entire length of the surgical accessories is exceeding the maximum permissible length defined for the HF surgical generator to be used

 Before to operate the HF surgical generator, verify that all the electrical connections have been properly established and no visible metallic parts are accessible on the papillotome's handle avoiding contact with wet tissues.

 The papillotome is a surgery accessory intended to deliver high-frequency electrical power to the body tissue thanks to the embedded electrical wire. Due to the very low electrical resistance, the device is not dissipating heat. It's a media to deliver the energy to the interface between the tip of the papillotome and the body tissues to be treated. The significant electrical resistant of the body tissues to be treated (in contact with the tip of the papillotome) induce heat to the tissues. During the short duration of power delivery by the HF generator, electromagnetic interferences can occur which could affect diagnostic or monitoring equipment leading to inappropriate assessment of the intraoperative state of the patient. Such potential interferences are inherent to the use of HF surgical generator in combination with surgery tool such as the papillotome.

11. Cleaning and Sterilization

 This device is designed and warranted for single use!

12. Storage

This device is delivered sterilized and in order to maintain this state until its first use it must be kept in its original sterilized package.

 Do not put any objects on the instrument or its package!
Do not store the instruments near aggressive chemical products!

 Do not expose the instruments to direct or indirect sunlight or other ultra-violet rays!

 Keep in dry area

Complaints will not be taken under consideration if the instruments have been stored improperly.

13. Risk in case of reuse

This devices is for single use only. G-Flex denied all responsibility in case of reuse. The reuse of single use device might present risk to the safety of the user or patient,

due to possible uncontrolled contamination and / or unreliability of the device performance.

14. Serious Incident, Complaints and repairs

In case of a serious incident, please report to the manufacturer and to the competent authority of the European Member State where the user is established. Shall the incident occur outside the Europe, please report to the manufacturer and follow the local applicable regulation.

In health protection of our employees only disinfected or sterilised instruments will be accepted for analysis or repairing. In any case disinfection or sterilisation dates and validity must be labelled outside the package. If this requirement has not been fulfilled the instruments will be returned to the sender without analysis or repairing.

15. Declarations

We, hereby, declare that the medical device neither contain nor is produced with tissue of animal origin.

We, hereby, declare that the medical device neither contain nor is produced with tissue or extracted material of human origin.

We, hereby, declare that the medical device do not contain medicinal substance.

We, hereby, declare that the medical device do not contain human blood derivative.

We, hereby, declare that the medical device do not contain phthalates or other products that could leak from the device that may be carcinogenic, mutagenic or toxic.

We, hereby, declare that the medical device is not composed of substances or combination of substances intended to be introduced into the human body.

16. Product disposal

After use, these products may be a potential biohazard. Handle and dispose of in accordance with accepted medical practice and applicable local, state and federal laws and regulations.

17. Legal foundation

The Law of the European Union applied.

18. Product support

In case of questions or difficulties concerning our instruments please contact your local distributor or G-FLEX directly during regular working hours.

Times available: Monday to Friday; 9am to 5pm (MET)

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