**EN - STONE EXTRACTION BALLOON CATHETERS**

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**Current Name:** IFU-StoneExtractionBalloonCatheter **Creation Date:** 2025-09-29 **Version** 15

*This document is only controlled in electronic format by the QM.*

|                            |           |
|----------------------------|-----------|
| EN - Catalog number        |           |
| FR - Numéro de catalogue   |           |
| NL - Catalogusnummer       |           |
| ES - Número de catálogo    |           |
| DE - Katalog-Nummer        |           |
| PT - Número do catalogo    |           |
| IT - Numero di catalogo    |           |
| EL - Αριθμός καταλόγου     | GF135716  |
| AR - العربية- رقم الكتالوج | GF135716T |
| TR - Katalog numarası      |           |
| SV - Katalognummer         |           |
| CS - Katalogové číslo      |           |
| DA - Katalognummer         |           |
| HU – Katalógusszám         |           |
| RO - Număr de catalog      |           |

|                                      |    |
|--------------------------------------|----|
| EN - INSTRUCTION FOR USE .....       | 4  |
| FR - MODE D'EMPLOI .....             | 7  |
| NL - GEBRUIKSAANWIJZING .....        | 10 |
| ES - INSTRUCCIONES DE USO .....      | 13 |
| DE - GEBRAUCHSANWEISUNG.....         | 16 |
| PT - INSTRUÇÕES DE USO.....          | 20 |
| IT - ISTRUZIONI PER L'USO .....      | 23 |
| EL - ΟΔΗΓΊΕΣ ΧΡΉΣΗΣ .....            | 26 |
| AR - العربية- تعليمات الاستخدام..... | 30 |
| TR - KULLANIM TALIMATLARI .....      | 33 |
| SV - BRUKSANVISNING .....            | 36 |
| CS – NÁVOD K POUŽITÍ .....           | 39 |

|                                     |    |
|-------------------------------------|----|
| DA - BRUGSANVISNING.....            | 42 |
| HU – HASZNÁLATI ÚTMUTATÓ .....      | 45 |
| RO - INSTRUȚIUNI DE UTILIZARE ..... | 48 |

# **EN - INSTRUCTION FOR USE**

## **1. Intended Use**

Stone Extraction Balloon Catheters are intended to remove stones from the biliary duct.

## **2. Intended Population**

Our device is intended for both adult and paediatric patient according to the physician recommendation taking into account the contraindications.

## **3. Contraindication**

Use to this natural latex rubber balloon is contraindicated in patients with a known hypersensitivity to latex and those specific to any other procedures to be performed in conjunction with extraction stone balloon.

## **4. Possible complication**

Potential complication: pancreatitis, cholangitis, aspiration, perforation, hemorrhage, infection, sepsis, allergic reaction to contrast or medication hypertension, respiratory depression or arrest, cardiac arrhythmia or arrest.

Additional complications that can occur during endoscopic balloon extraction include: stone impaction localized inflammation, pressure necrosis.

## **5. User**

The users of G-FLEX instruments must be specialists in their fields. An appropriate and specific training for preparation, care and maintenance of the flexible instruments is required.

## **6. Product Description**

The device is composed of a triple lumen pebax tube with a 3-way splitting connector at the proximal end. It divide the 3 lumen catheter into one luer-lock connector with a 2-way valve to inflate the balloon, one luer-lock connector for a 0.035" guide wire and one luer-lock connector for injection. The Stone Extraction Balloon is on the distal end. The multistage diameters of the balloon are reached using the different syringes supplied.

## **7. Advices before first application**

Please read carefully and follow all safety operating instructions and warnings before first application of the device.

A previous knowledge regarding handling and operation is required and essential.

Unpack the instrument carefully and examine the devices for any possibility of damage. In case of any damage or missing items contact your distributor immediately.

If the package present any damaged, the sterility of the device is not guarantee.

## 8. Handling and Operation

It's important to irrigate the injection and guide wire lumen before introducing the device into the endoscope.

The device must be inserted deflated into the endoscope operating channel. Please observe compatibility of catheter diameter and the endoscope operating channel diameter (see catalogue).

If possible, the use of a guide wire is recommendable.

## 9. Procedure

After locating the stone to be removed, proceed as follows:

- Insert the deflated balloon into the biliary duct through the guide wire until the distal tip is located above the stone to be removed.
- Choose one of the syringes supplied with the balloon according to desired balloon diameter.
- Fill the selected syringe with air and connect it to the proximal luer-lock connector labeled "Balloon" along with the stopcock in the open position
- Push the plunger of the syringe carefully to inflate the balloon and close the stopcock.
- Gently withdraw the inflated balloon catheter to pull the stones out of the biliary duct and drop them into the duodenum. They can be expelled naturally out of the body through the intestines.
- Open the stopcock and deflate the balloon by pulling the plunger of the syringe.
- Remove the entire device out of the endoscope carefully.

## 10. Precautions/ Warnings

At anytime if injection is required, use the connector labeled "Injection"

Sharp stones can destroy the balloon and should be removed with an extractor basket or a lithotripsy basket.

## 11. Cleaning and Sterilization

This device is designed and warranted for single use!

## 12. Storage

This device is delivered sterilized and in order to maintain this state until its first use it must be kept in its original sterilized package.

Do not put any objects on the instrument or its package!  
Do not store the instruments near aggressive chemical products!

Do not expose the instruments to direct or indirect sunlight or other ultra-violet rays!

Keep in dry area

Complaints will not be taken under consideration if the instruments have been stored improperly.

### **13. Risk in case of reuse**

This device is for single use only. G-Flex denied all responsibility in case of reuse. The reuse of single use device might present risk to the safety of the user or patient, due to possible uncontrolled contamination and / or unreliability of the device performance.

### **14. Serious Incident, Complaints and repairs**

In case of a serious incident, please report to the manufacturer and to the competent authority of the European Member State where the user is established. Shall the incident occur outside the Europe, please report to the manufacturer and follow the local applicable regulation.

In health protection of our employees only disinfected or sterilised instruments will be accepted for analysis or repairing. In any case disinfection or sterilisation dates and validity must be labelled outside the package. If this requirement has not been fulfilled the instruments will be returned to the sender without analysis or repairing.

### **15. Declarations**

We, hereby, declare that the medical device neither contain nor is produced with tissue of animal origin.

We, hereby, declare that the medical device neither contain nor is produced with tissue or extracted material of human origin.

We, hereby, declare that the medical device do not contain medicinal substance.

We, hereby, declare that the medical device do not contain human blood derivative.

We, hereby, declare that the medical device do not contain phthalates or other products that could leak from the device that may be carcinogenic, mutagenic or toxic.

We, hereby, declare that the medical device is not composed of substances or combination of substances intended to be introduced into the human body.

### **16. Product disposal**

After use, these products may be a potential biohazard. Handle and dispose of in accordance with accepted medical practice and applicable local, state and federal laws and regulations.

### **17. Legal foundation**

The Law of the European Union applied.

### **18. Product support**

In case of questions or difficulties concerning our instruments please contact your local distributor or G-FLEX directly during regular working hours.

**Times available:** Monday to Friday; 9am to 5pm (MET)

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