EN - BITE BLOCK
FR - CALE-DENTS
NL - BIJTBLOK
ES - PROTECTOR BUCAL
DE - BEIßBLOCK
PT - ABRE-BOCA
IT - BLOCCA MORSO
EL - ΣΤΟΜΑΤΟΔΙΑΣΤΟΛΕΑΣ
AR - قالب العضة
TR - ISIRMA BLOĞU
SV - BITBLOCK
CS – SKUSOVÝ BLOK
DA - BIDEBLOK
HU – HARAPÁSGÁTLÓ
RO - PIESĂ BUCALĂ

G-Flex Europe S.R.L. - 20, Rue de l'industrie 1400
Nivelles - Belgium | **Phone:** +32 67.88.36.65
Email: info@g-flex.com | **Website:** www.g-flex.com

EN - Catalog number	
FR - Numéro de catalogue	
NL - Catalogusnummer	
ES - Número de catálogo	
DE - Katalog-Nummer	
PT - Número do catalogo	
IT - Numero di catalogo	
EL - Αριθμός καταλόγου	U600005
AR - العربية- رقم الكتالوج	U600006
TR - Katalog numarası	
SV - Katalognummer	
CS - Katalogové číslo	
DA - Katalognummer	
HU – Katalógusszám	
RO - Număr de catalog	

EN - INSTRUCTION FOR USE	4
FR - MODE D'EMPLOI	6
NL - GEBRUIKSAANWIJZING	9
ES - INSTRUCCIONES DE USO	11
DE - GEBRAUCHSANWEISUNG.....	14
PT - INSTRUÇÕES DE USO.....	16
IT - ISTRUZIONI PER L'USO	19
EL - ΟΔΗΓΊΕΣ ΧΡΉΣΗΣ	21
AR - العربية - تعليمات الاستخدام.....	24
TR - KULLANIM TALIMATLARI	27
SV - BRUKSANVISNING	29
CS – NÁVOD K POUŽITÍ	32
DA - BRUGSANVISNING.....	34
HU – HASZNÁLATI ÚTMUTATÓ	37
RO - INSTRUCTIUNI DE UTILIZARE	39

EN - INSTRUCTION FOR USE

1. Intended Use

The bite block is used to protect the endoscope insertion tube and esophageal dilators from being bitten by the patient.

2. Intended Population

Our device is intended for both adult and paediatric patient according to the physician recommendation taking into account the contraindications.

3. Contraindication

There are no known contraindications.

4. Possible complication

There are no known complications.

5. User

The users of G-FLEX instruments must be specialists in their fields. An appropriate and specific training for preparation, care and maintenance of the flexible instruments is required.

6. Product Description

The device is made of a plastic part to be placed in the mouth and a strap be fix it around the patient neck.

7. Advices before first application

Please read carefully and follow all safety operating instructions and warnings before first application of the device.

A previous knowledge regarding handling and operation is required and essential. Unpack the instrument carefully and examine the devices for any possibility of damage. In case of any damage or missing items contact your distributor immediately.

8. Handling and Operation

Avoid placing fingers between the bite block and the patient's teeth when placing the bite block in the patient's mouth.

9. Procedure

- Detach one side of the strap from the bite block before placing the bite block into the patient's mouth.
- Place the bite block into the patient's mouth
- Pass the strap around the patient's neck and attach it again to the bite block

10. Precautions/ Warnings

There are no specific warnings for this device

11. Cleaning and Sterilization

This device is designed and warranted for single use!

12. Storage

This device is delivered non sterile and in order to maintain this state until its first use it must be kept in its original package.

Do not put any objects on the instrument or its package!

Do not store the instruments near aggressive chemical products!

Do not expose the instruments to direct or indirect sunlight or other ultra-violet rays!

Keep in dry area

Complaints will not be taken under consideration if the instruments have been stored improperly.

13. Risk in case of reuse

This device is for single use only. G-Flex denied all responsibility in case of reuse. The reuse of single use device might present risk to the safety of the user or patient, due to possible uncontrolled contamination and / or unreliability of the device performance.

14. Serious Incident, Complaints and repairs

In case of a serious incident, please report to the manufacturer and to the competent authority of the European Member State where the user is established. Shall the incident occur outside the Europe, please report to the manufacturer and follow the local applicable regulation.

In health protection of our employees only disinfected or sterilised instruments will be accepted for analysis or repairing. In any case disinfection or sterilisation dates and validity must be labelled outside the package. If this requirement has not been fulfilled the instruments will be returned to the sender without analysis or repairing.

15. Declarations

We, hereby, declare that the medical device neither contain nor is produced with tissue of animal origin.

We, hereby, declare that the medical device neither contain nor is produced with tissue or extracted material of human origin.

We, hereby, declare that the medical device do not contain medicinal substance.

We, hereby, declare that the medical device do not contain human blood derivative.

We, hereby, declare that the medical device do not contain phthalates or other products that could leak from the device that may be carcinogenic, mutagenic or toxic.

We, hereby, declare that the medical device is not composed of substances or combination of substances intended to be introduced into the human body.

16. Product disposal

After use, these products may be a potential biohazard. Handle and dispose of in accordance with accepted medical practice and applicable local, state and federal laws and regulations.

17. Legal foundation

The Law of the European Union applied.

18. Product support

In case of questions or difficulties concerning our instruments please contact your local distributor or G-FLEX directly during regular working hours.

Times available: Monday to Friday; 9am to 5pm (MET)

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