EN - CLEANING BRUSH
FR - BROSSE DE NETTOYAGE
NL - REINIGINGSBORSTEL
ES - CEPILLO DE LIMPIEZA
DE - REINIGUNGSBÜRSTE
PT - ESCOVA DE LIMPEZA
IT - SPAZZOLA PER LA PULIZIA
EL - ΒΟΥΡΤΣΑ ΚΑΘΑΡΙΣΜΟΥ
AR - فرشاة التنظيف
TR - TEMİZLİK FIRÇASI
SV - RENGÖRINGSBORSTE
CS – ČISTICÍ KARTÁČEK
DA - RENGØRINGSBØRSTE
HU – TISZTÍTÓKEFE
RO - PERIE DE CURĂȚARE

G-Flex Europe S.R.L. - 20, Rue de l'industrie 1400
Nivelles - Belgium | **Phone:** +32 67.88.36.65
Email: info@g-flex.com | **Website:** www.g-flex.com

EN - Catalog number	
FR - Numéro de catalogue	
NL - Catalogusnummer	
ES - Número de catálogo	
DE - Katalog-Nummer	
PT - Número do catalogo	
IT - Numero di catalogo	
EL - Αριθμός καταλόγου	GF4110
AR - العربية- رقم الكتالوج	GF4113
TR - Katalog numarası	
SV - Katalognummer	
CS - Katalogové číslo	
DA - Katalognummer	
HU – Katalógusszám	
RO - Număr de catalog	

EN - INSTRUCTION FOR USE	4
FR - MODE D'EMPLOI	8
NL - GEBRUIKSAANWIJZING	13
ES - INSTRUCCIONES DE USO	18
DE - GEBRAUCHSANWEISUNG.....	22
PT - INSTRUÇÕES DE USO.....	27
IT - ISTRUZIONI PER L'USO	32
EL - ΟΔΗΓΊΕΣ ΧΡΉΣΗΣ	36
AR - العربية - تعليمات الاستخدام.....	42
TR - KULLANIM TALIMATLARI	46
SV - BRUKSANVISNING	50
CS – NÁVOD K POUŽITÍ	55
DA - BRUGSANVISNING.....	59
HU – HASZNÁLATI ÚTMUTATÓ	64
RO - INTRUCȚIUNI DE UTILIZARE	68

EN - INSTRUCTION FOR USE

1. Intended Use

The Cleaning Brush is meant to clean the biopsy channel of endoscopes or the endoscope valves. The Sponge is intended to clean the outer part of the endoscope or protect the endoscope during transportation.

2. Intended Population

Not applicable

3. Contraindications

There are no contraindications currently for this device family.

4. Possible Complications

Not applicable

5. User

The users of G-FLEX instruments must be specialists in their fields. An appropriate and specific training for preparation, care and maintenance of the flexible instruments is required.

6. Product Description

The device is made of a long catheter (or spiral) with a brush on each side.

7. Advices before first application

Please read carefully and follow all safety operating instructions and warnings before first application of the device.

A previous knowledge regarding handling and operation is required and essential. Unpack the instrument carefully and examine the devices for any possibility of damage. In case of any damage or missing items contact your distributor immediately.

8. Handling and Operation

The Cleaning Brush has to be passed through the full length of the endoscope biopsy channel several times with the solution described in your endoscope manufacturer notice.

9. Procedure

1. Introduce the brush through one side of the biopsy channel and pull it out from the other side of the same.
2. Repeat this process until the working channel is free from dirt. (Recommended a minimum of 3 times)
3. Ensure to clean the brush before introducing it in the channel each time.

10. Precautions/Warnings

Please note that some models have the distal and proximal ends different from each other. In such models, make sure to introduce the brush with the thinnest diameter tip into the channel.

11. Cleaning and Sterilization

Before the first use and prior to all subsequent uses, the device must be cleaned, disinfected and sterilised. Effective cleaning and disinfection are essential for an effective sterilisation.

The reprocessing of reusable endoscopic equipment must remain under the responsibility of individuals with thorough training in reprocessing methods and be aware of the following points:

- Local hospital's reprocessing guidelines
- National health and safety guidelines
- The instructions in this instruction for use
- Additional national specifications related to effective inactivation of prions (i.e: use of highly alkaline cleaning agent and prolonged sterilisation time)
- That the methods used to clean, disinfected and sterilise are adequately validated for this type of instruments
- That the equipment used is regularly maintained and checked
- That the validated parameters are maintained during each reprocessing cycles

The reprocessing of the instrument must be initiated immediately after use, since dried organic residue can impact effective cleaning, disinfection and sterilisation.

Reprocessing methods

A mechanical cleaning/disinfection process should be used whenever possible over manual cleaning.

Therefore, one of the following 2 options must be used for the reprocessing:

1. Perform the pre-treatment and manual cleaning before the mechanical cleaning/disinfection followed by the sterilization (3 steps).
2. In case a mechanical cleaning/disinfection is not available the device must then be pre-treated and manually cleaned followed by a sterilization (2 steps).

Material Resistance

Ensure that the cleaning and disinfection agent do not contain the following:

- Aldehydes (this will fix blood-based soiling)
- Strong organic, mineral or oxidising acids (lowest permissible pH value 5.5)
- Strong alkalis (highest permissible pH 11, neutral / enzymatic or slightly alkaline cleaners are recommended)
- Organic solvents (e.g. ether, ketones, benzines), fluorinated alcohols
- Oxidising agents (e.g. hydrogen peroxide)
- Halogens (chlorine, iodine, bromine)
- Aromatic / halogenated hydrocarbons
- Formamide
- Trichloroethylene / perchloroethylene

Pre-treatment / Manual Cleaning

1. Immersed the instrument in cleaning solution within maximum 2 hours after use
2. Remove any soiling visible on the jaws using a clean soft brush (or a clean, soft and lint-free cloth) which is used solely for this purpose. Never use metal brushes or wire wool.
3. Open and close the jaws, immersed in cleaning solution at least 3 times
4. Immediately place the instrument winded with a diameter not less than 20 cm into an ultrasound bath filled with cleaning solution for 30 minutes (time and concentration may vary according to the solution manufacturer's specifications. Cleaning method has been validated with ultrasound bath using 0.5% neodisher® Mediclean forte (Dr. Weigert GmbH) at 35 kHz)
5. Rinse the instrument for at least 1 minute under running water (temperature < 35°C/85°F)

Mechanical Cleaning / Disinfection

Proceed as followed:

1. Place the instrument winded with a diameter not less than 20 cm into the cleaning and disinfection equipment that fulfils the following requirement:
 - Equipment is effective and suitable for the cleaning and disinfection of medical instruments (i.e: FDA approved, CE marked)
 - Only sterile or low-germ (max. 10 bacteria/ml) and low-endotoxin (max. 0.25 endotoxin units /ml) water is used (e.g. purified water/highly purified water)
 - The air used for drying is filtered (oil free, low-germ and low-particle)
 - The device is maintained and checked regularly
 - The cleaning and disinfection solutions used are compatible with the instrument (See "Material Resistance")
 - The concentration, temperatures and exposure times specified by the cleaning / disinfectant agent manufacturer is followed at all time, this includes the rinsing specifications
2. Start the cleaning and disinfection program:
 - The programme must be suitable for the instrument
 - Thermal disinfection programme should be used (90°C / 194°F for at least 5 minutes)
 - The program contains at least 3 rinsing cycles (depleting steps) after cleaning (or neutralization if applied) or conductance control recommended to effectively prevent detergent residues
3. Remove the instruments from the equipment when the program is finished
4. Check the functionality of the device by opening and closing the jaws several times
5. Check the device for signs of corrosion, damaged surfaces, cracks, soiling or discolouration
6. Remove any damaged device from use
7. Any soiled device must be cleaned and disinfected again

8. Package the device into a single-use sterilisation pack or a sterilisation container (suitable for steam sterilisation of at least 134°C / 270 °F) as soon as possible after removal from the equipment

Sterilisation

Only steam sterilisation must be used.

Flash sterilisation, hot air sterilisation, radiation sterilisation, formaldehyde, ethylene oxide sterilisation or plasma sterilisation are not permitted.

The following aspects shall be taken into account:

- Fractionated vacuum method (at least 3 vacuum steps and with adequate product drying not less than 20 minutes)
- Steam steriliser compliant with DIN EN 13060/DIN EN 285 or ANSI AAMI ST79
- Validated in accordance with DIN EN ISO 17665 (valid IQ/OQ and product-specific performance assessment (PQ))
- Maximum sterilisation temperature 134°C / 273°F plus tolerance as defined in DIN EN ISO 17665
- Sterilisation time (exposure time at sterilisation temperature) at least 4 min at 134°C/270°F).
- For prion inactivation (depending on the local recommendation) at least 18 minutes at 134°C/273°F

12. Storage

Do not put any objects on the instrument or its package!

Do not store the instruments near aggressive chemical products!

Do not expose the instruments to direct or indirect sunlight or other ultra-violet rays!

Keep in dry area

Complaints will not be taken under consideration if the instruments have been stored improperly.

13. Risk in case of reuse

The quantity of reprocessing cycles the device can withstand is closely related to the care given to the device during its intended use and the care taken during the cleaning and sterilization. We have validated 200 reprocessing cycles but the user can use for longer looking, in between each use, for signs of performance degradation such as kinks on the spiral, damaged to the jaws and stain. Shall any of these be present, it is recommended to segregate the device and use another one.

14. Serious Incident, Complaints and repairs

In case of a serious incident, please report to the manufacturer and to the competent authority of the European Member State where the user is established.

Shall the incident occur outside the Europe, please report to the manufacturer and follow the local applicable regulation.

In health protection of our employees only disinfected or sterilised instruments will be accepted for analysis or repairing. In any case disinfection or sterilisation dates and validity must be labelled outside the package. If this requirement has not been fulfilled the instruments will be returned to the sender without analysis or repairing.

15. Declarations

We, hereby, declare that the medical device neither contain nor is produced with tissue of animal origin.

We, hereby, declare that the medical device neither contain nor is produced with tissue or extracted material of human origin.

We, hereby, declare that the medical device do not contain medicinal substance.

We, hereby, declare that the medical device do not contain human blood derivative.

We, hereby, declare that the medical device do not contain phthalates or other products that could leak from the device that may be carcinogenic, mutagenic or toxic.

We, hereby, declare that the medical device is not composed of substances or combination of substances intended to be introduced into the human body.

16. Product disposal

After use, these products may be a potential biohazard. Handle and dispose of in accordance with accepted medical practice and applicable local, state and federal laws and regulations.

17. Legal foundation

The Law of the European Union applied.

18. Product support

In case of questions or difficulties concerning our instruments please contact your local distributor or G-FLEX directly during regular working hours.

Times available: Monday to Friday; 9am to 5pm (MET)

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